EN

PRESSIO® CABLES

For use with Pressio Monitoring Systems only

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Table of Contents

EN

Pressio® Cables	7
1. Preamble	7
2. Intended use	7
3. Description	7
4. Cleaning and disinfection procedure	7
4.1. Introduction	8
4.2. Prerequisites	8
4.3. Cleaning procedure	8
4.4. Disinfection procedure	8
4.5. Inspection	8
5. Using the cables	8
5.1. Catheter Extension Cable	8
5.2. Patient Monitor Cable	9
5.3. Cables connection	9
6. Maintenance of the cables	9
7. Environmental conditions, storage and shipping	9
7.1. Environmental conditions	9
7.2. Storage	9
7.3. Shipping	9
8. Processing the product after use	9
8.1. Product return	9
8.2. Product elimination	9
9. Warranty	9
10. Symbols	10
11. Electromagnetic compatibility	10
12. References	11

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CAUTION

Federal (USA) law restricts this device to sale by or on the order of a physician.

1. Preamble

These Instructions for Use detail all the information required to identify, connect, and clean the cables in use with a Pressio Interface PSO-IN00, the Pressio intracranial temperature module (PSO-MT00), the Pressio Monitor (PSO-3000 or PSO-4000), and the Pressio Catheters range ; hereinafter referred to as Monitor and Catheter.

Any request for information or modification relating to this guide should be sent to: Sophysa – 5, rue Guy Moquet – 91400 Orsay – France.

2. Intended use

Sophysa provides three types of cables to be used exclusively with its Monitoring systems.

– Catheter Extension Cable

The Catheter Extension Cable is intended to be used exclusively with the Pressio 2 ICP Monitor (PSO-4000). It makes it possible to electrically connect the Pressio Catheter to the Pressio 2 ICP Monitor. It conducts the analog signal from ICP and ICT sensors to the Pressio 2 ICP Monitor.

It allows the Pressio 2 ICP Monitor to communicate with the Pressio Catheter memory in the form of a digital signal.

– Pressure cable

The pressure cable is intended to be used with Pressio Monitoring systems (the Pressio ICP Monitor (PSO-3000), the Pressio 2 ICP Monitor (PSO-4000), the Pressio intracranial temperature module (PSO-MT00) and the Pressio interface (PSO-IN00)).

It makes it possible to connect monitoring systems to the patient monitor and transmit the ICP value in the form of an analog signal.

– Temperature cable

The temperature cable is intended to be used with Pressio Monitoring systems (the Pressio ICP Monitor (PSO-3000), the Pressio intracranial temperature module (PSO-MT00) and the Pressio 2 ICP Monitor (PSO-4000)). It makes it possible to connect monitoring systems to the patient monitor and transmit the ICT value in the form of an analog signal.

3. Description

The cables are indicated for use with the Monitor and Catheters only.

The Monitor must be used with the Sophysa-supplied cables listed below.

NOTE

The Catheter Extension Cable (PSO-EC30) must only be used with the catheters listed in the *Pressio® Catheter Kits* table (Section 12. References (p. 11)).

NOTE

Cables connections illustrated in this document are examples provided with the use of a Pressio 2 ICP Monitor (PSO-4000).

Figure 1. Pressio 2 system configuration


Table 1. Cables references

Item	Product code	Description
A	-	To mains power
B	-	To patient bedside monitor
C	-	To PC
D	-	To patient
1	Power cable	<i>Included in the Monitor package</i>
2	PSO-MCxx	Patient monitor/pressure cable Length: 2.9m <i>Supplied separately</i>
3	PSO-MCT-y	Patient monitor/temperature cable Length: 2.9m <i>Supplied separately</i>
4	Mini-USB	PC/Pressio 2 ICP Monitor Length: 1.8m <i>Included in the Monitor package</i>
5	PSO-EC30	Catheter Extension Cable Length: 2 m <i>Included in the Monitor package, also available separately</i>

4. Cleaning and disinfection procedure**CAUTION**

Do not use a cable if there is any soiling or visible residues on the external surface or inside the wires.

4.1. Introduction

The cables are delivered clean, but not disinfected.

Clean all the cables and then disinfect the Catheter Extension Cable before first use and between each patient, as described below.

WARNING

Do not clean the Monitor or the Catheter, nor disinfect the cables, when in use on the patient.

WARNING

Do not immerse, autoclave or soak the cables in a liquid. Their performance (including drift and electrical safety) could be affected.

NOTICE

Do not use solvents or cleaning agents which could damage the Monitor cables, such as:

- cleaning/disinfection agents (based on phenols),
- cleaning/disinfection by boiling,
- cleaning/disinfection by hot air/steam,
- acetone, ammonia, benzene, bleaching agent, chlorine, chlorine water, water above 60°, paint solvents, trichloroethylene.

For more information, contact Sophysa Customer Service at contact@sophysa.com, or contact your local distributor.

4.2. Prerequisites

1. Put on gloves and keep them throughout the procedure.
2. Take pre-soaked wipes with 70% isopropyl alcohol (IPA).

4.3. Cleaning procedure

The purpose of this procedure is to remove any soiling and visible residues on the cables.

NOTE

Do not put excessive pressure on the product labels.

1. Clean the cable surfaces for at least 1 minute using pre-soaked wipes with 70% isopropyl alcohol (IPA), to remove any visible residues. Thoroughly wipe all the cables, using a new wipe for each cable.
2. Inspect the components.
If residues remain, take a new pre-soaked wipe with 70% isopropyl alcohol (IPA) and wipe the surfaces again.

NOTE

Repeat this step until all visible residues are removed from all components.

3. Disinfect the Catheter Extension Cable as described in the following section.

4.4. Disinfection procedure

The purpose of this procedure is to remove any microorganisms present on the Catheter Extension Cable.

1. Thoroughly wipe the Catheter Extension Cable. It must remain visibly wet for at least 2 minutes.

If needed, use additional wipes to ensure continuous 2 minutes of wet contact time.

2. Let the Catheter Extension Cable completely air dry for 1 hour before using it again.

4.5. Inspection

After each cleaning or disinfection procedure, visually inspect the cables for any damage. See *Section 6. Maintenance of the cables* (p. 9).

5. Using the cables

NOTE

Pressio Cables should only be used by trained personnel in medical devices.

Check that neither part of the cable has been damaged during transport.

CAUTION

If any part of the cable is damaged, do not use it. Contact Sophysa or your local distributor.


WARNING

Check proper cables connection before using your monitoring system. Bad connection could lead to the wrong diagnosis and thus risk patient injury.

5.1. Catheter Extension Cable

The Catheter Extension Cable is supplied with every PSO-4000 Monitor and is also available separately.

Figure 2. Connecting system

[6] Plug [7] Fixation clips [8] Dongle Inlet

The plug (6) connects to the Monitor. A guide pin and color code (white) facilitate proper connection.

The fixation clips (7) are tied to the bed sheets, or to the patient's clothing. Proper use of these clips limits traction on the implanted Catheter, and reduces the risk of Catheter disconnection.

The dongle inlet (8) connects to the Catheter. A guide pin and color code (blue arrows) facilitate proper connection.

At the end of the cable, insert the Catheter until you reach the blue line. The line on the Catheter dongle indicates the point of complete insertion.

CAUTION

Do not allow any connectors to come into contact with liquid.

WARNING

The Catheter Extension Cable is not protected against the effects of a cardiac defibrillator.

5.2. Patient Monitor Cable**NOTE**


Follow cable compatibility detailed in the *Pressio® Accessories* table (*Section 12. References (p. 11)*).

The pressure and temperature cables are fitted with guide pins, to facilitate proper connections.

5.3. Cables connection

The three cable sockets are placed on the Monitor right panel, with a specific color code for each cable.

- The Catheter Extension Cable socket (white) is used to connect the Catheter Extension Cable plug, to allow Catheter/Monitor connection.
- The ICP socket (blue) is used to connect the Monitor pressure cable plug to the patient bedside monitor for the transmission of ICP values.
- The ICT socket (red) is used to connect the Monitor temperature cable plug to the patient bedside monitor for the transmission of ICT values.

To connect the cables, align the arrows located on the plug with the matching marker [9]. Then, push the cable plug in. It should enter easily without forcing.

To disconnect the cables, pull the cable plug away from the Monitor socket.

6. Maintenance of the cables

Sophysa recommends continuous vigilance to ensure that the Monitor, its cables and accessories are in good operating condition for every new patient.

Check the visual aspect of the cables before use. Make sure that:

- there is no foreign body in the cable plug,
- electrical connections of the plug are not twisted,
- there are no visible cracks on the cable,
- cable markings are still visible.

Contact Sophysa Customer Service at contact@sophysa.com to organize the shipment of your equipment.

7. Environmental conditions, storage and shipping

Do not place the Monitor cables in an MRI magnetic field.

WARNING

Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to the Monitor cables. Otherwise, degradation of the performance of this equipment could result.

7.1. Environmental conditions

The Monitor and its accessories are designed to withstand the following environmental conditions:

- Temperature: between +10 °C (50 °F) and +40 °C (104 °F).
- Relative humidity without condensation: between 15% and 95%.
- Altitude: between -500 m and +3,000 m (corresponding to a flight altitude up to 12,000 m in a pressurized cabin).

7.2. Storage

The Monitor, its accessories and cables are designed to withstand the following storage conditions:

- Temperature: between -20 °C (-4 °F) and +60 °C (140 °F).
- Relative humidity: between 15% and 95%.
- Altitude: between -500 m and +4,600 m (corresponding to a flight altitude up to 12,000 m in a pressurized cabin).

7.3. Shipping

Protect the Monitor, its accessories and cables from shocks and vibrations during shipping.

The Monitor, its accessories and cables are designed to withstand the following shipping conditions:

- Temperature: between -20 °C (-4 °F) and +60 °C (140 °F).
- Relative humidity: between 15% and 95%.
- Altitude: between -500 m and +4,600 m (corresponding to a flight altitude up to 12,000 m in a pressurized cabin).

8. Processing the product after use**8.1. Product return**

Do not use a cable that is damaged.

As part of its continual improvement program, Sophysa asks its customers to inform it and the legal authority of the country of any unexpected and serious problems that occur with the product.

Return the faulty cable to Sophysa for proper investigation. In order to properly assess the returned cable, it must be accompanied by an explanatory *Return to Manufacturer Authorization* form.

8.2. Product elimination

In order to eliminate a cable, clean it as instructed in *Section 4. Cleaning and disinfection procedure (p. 7)* and then send the cable back to Sophysa for proper elimination. Contact Sophysa customer service at contact@sophysa.com or your local distributor to organize the shipment of your equipment.

9. Warranty

The performances and safety of the cables are ensured only with the Monitor and the compatible components (such as

the Catheter and the patient bedside monitor) designed, tested and manufactured by Sophysa.

Sophysa warrants the performance and safety of this medical device under the normal conditions of the intended use of the device, adapted to its intended purpose and use, and in accordance with these Instructions for Use.

The medical device must be stored and transported in an environment and under conditions that also comply with the information in these Instructions for Use. These storage and transport conditions have been tested and validated by Sophysa. Thus, Sophysa does not grant any other express or implicit guarantee as for the good conservation and the safety of the product in other premises than its own which would not respect these conditions. Likewise, no express or implicit guarantee is granted by Sophysa as to the suitability of the product for the use which will be made of it, or its adaptation to a particular use, except within the indications and the intended purpose of the product, or when it has been transformed, modified or repaired except within the instructions of Sophysa.

Under no circumstances, Sophysa can be held responsible in case of damages, for any incident and/or complication, resulting from damage or prejudice arising directly or indirectly from the unsuitable use of the device and/or a use of the device which fail to conform or the non-respect of its conditions of maintenance, cleaning, maintenance, storage or transport.

The cables have an expected lifetime of 2 years in normal use (disconnecting and reconnecting between each monitoring). However, the results of the visual checks prevail (see Section 6. *Maintenance of the cables* (p. 9)).



These visual checks will indicate whether or not the cables can still be used.

Table 2. Date of first use

Date	
------	--

10. Symbols

	Catalog number
	Lot number
	Manufacturer
	Date of manufacture
	CE conformity marking
	Consult Instructions for Use

	Humidity limitation: 15% to 95%
	Temperature limits: -20°C (-4° F) to 60°C (140°F)
	The Catheter Extension Cable is not protected against the effects of a cardiac defibrillator.
	MR unsafe The cables are not suitable for use in an MRI environment. Do not use during an MRI examination.
	Medical Device
	Unique Device Identification
	Keep dry
	Keep away from sunlight
	By prescription only
	Discarding this type of product with other waste is prohibited.

11. Electromagnetic compatibility

In case of electrostatic discharge or electromagnetic disturbance, find below the behavior that the Monitor, the cables or the sensor might have:

- If a strong electrostatic discharge occurs on the Monitor, the cables or the sensor:
 - the mean ICP or ICT value can be temporarily affected,
 - the ICP curve can show temporary peaks of high amplitude.
- If an electromagnetic disturbance occurs on the Monitor, the cables or the sensor:
 - the mean ICP or ICT value can be affected,
 - the ICP curve can show abnormal ripples.

NOTE

The emission characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in a residential environment (for which CISPR 11 class B is normally required) this equipment might not offer adequate protection to radio-frequency communication services.

You might need to take mitigation measures, such as relocating or re-orienting the equipment.

12. References

Table 3. Pressio® Catheter Kits

PSO-PB	Pressio® ICP monitoring kit, parenchymal with bolt
PSO-PBT	Pressio® ICP and ICT monitoring kit, parenchymal with bolt
PSO-PT	Pressio® ICP monitoring kit, parenchymal tunneling
PSO-PTT	Pressio® ICP and ICT monitoring kit, parenchymal tunneling
PSO-VT	Pressio® ICP monitoring kit, ventricular tunneling with external CSF drainage function
PSO-VTT	Pressio® ICP and ICT monitoring kit, ventricular tunneling with external CSF drainage function

Table 4. Pressio® Monitoring System

PSO-4000	Pressio® 2 ICP Monitor Power Cable and Catheter Extension Cable (PSO-EC30) included
----------	----------------------------------------------------------------------------------------

Table 5. Pressio® Accessories

PSO-EC30	Catheter extension cable (included with PSO-4000, also available separately)
PSO-MCxx	Pressure cable <ul style="list-style-type: none"> – MC01: PHILIPS (AGILENT) - 12 pins <i>Also compatible with Philips (Intellivue, Efficia), Mindray</i> – MC02: SIEMENS (SIRECUST) - 10 pins <i>Also compatible with Dräger (Infinity)</i> – MC03: SPACELABS & MINDRAY - 6 pins – MC04: GE DATEX-Ohmeda - 10 pins – MC05: GE Solar (MARQUETTE) - 11 pins <i>Also compatible with GE Carescope</i> – MC06: HELLIGE - 10 pins – MC07: SIEMENS - 7 pins <i>Also compatible with Dräger (Infinity)</i> – MC08: NIHON KOHDEN - 5 pins – MC10: DATASCOPE - 6 pins
PSO-MCT-y	Temperature cable <ul style="list-style-type: none"> – MCT-A: PHILIPS (AGILENT) - 2 pins <i>Also compatible with Philips (Intellivue, Efficia)</i> – MCT-B: SIEMENS - 7 pins <i>Also compatible with Dräger (Infinity)</i> – MCT-C: SPACELABS - 10 pins – MCT-E: GE Solar (MARQUETTE), GE DATEX-Ohmeda - 11 pins – MCT-F: HELLIGE, DATEX-Ohmeda, NIHON KOHDEN, MINDRAY & DATASCOPE - JACK 6.35 mm

Table 6. Year of first CE marking

PSO-EC30	2016
PSO-MCxx	2005
PSO-MCT-y	2010

Pressio® Cables comply with the following standards:

- EN ISO 14971: 2019
- IEC 60601-1: 2012
- IEC 60601-1-2: 2014
- NF EN 62366-1:2008
- IEC 60601-1-6: 2010
- ISO 10993-1: 2018

Technical specifications and list of product references may be modified without notice.

Availability may vary according to country.

This page is intentionally left blank.

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Pressio® is a registered trademark of Sophysa.
©2020 Sophysa. All rights reserved.