PATIENT IDENTIFICATION CARD
CARTE IDENTIFICATION PATIENT

For Hydrocephalus
SOPHY® Adjustable Valve

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The adjustment spring of Sophy® valves is calibrated in order to make the magnetic rotor insensitive to the influence of everyday magnetic fields. Magnetic fields generated by airport security doors, microwave oven, mobile phone, high tension wire, TV, etc., should not affect the pressure setting.

The permanent magnets in domestic appliances, toys, audio headphones, loudspeakers, and electromagnetic fields created by electric motors in shavers, hair dryers, clippers, etc., while turned on, can modify the pressure setting. Thus, in order to prevent any risk of pressure alteration, they must not be brought close to the valve site. MRI exposure can alter the pressure setting of the Sophy® valve.

It is the responsibility of the physician to inform the patient or the patient's family to avoid all activity exposing the device to shocks (violent sports, etc.) or to a strong magnetic or electromagnetic field, likely to cause dysadjustment of the pressure setting. The valve pressure setting must systematically be confirmed before and after exposure to MRI, after strong magnetic or electromagnetic field exposure, and in case of a shock on the implantation site.

Unexplained change of the pressure setting has been reported. Abnormal clinical signs between two appointments must lead to urgent neurosurgical advice.

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SOPHYSA USA Inc. : 303 S Main Street, Crown Point IN 46307 - USA