Instructions for use

Sophy® Adjustable Pressure Valve

Sterile, disposable valve for Cerebrospinal Fluid drainage

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
Sterile, single use

Table of contents

1.	Indications	6
2.	Contraindications	6
3.	Description and Operating Principle of the Sophy [®] valve	6
3.1.	Adjustable pressure	6
4.	Principle of the Sophy [®] valve adjustment	7
4.1.	Locating Instrument	8
4.2.	Reading Instrument	8
4.3.	Setting Instrument	8
4.4.	Demonstration valve	8
5.	Configurations of the adjustable Sophy [®] valve	8
6.	Measurement Unit and Calibration of Operating Pressures	9
7.	Behavior during Magnetic Resonance Imaging (MRI)	9
8.	Sterilization – Decontamination of Sophy [®] Valves and Valve Kits	9
9.	Instructions	10
9.1.	Choice of valve model	10
9.2.	Adjustment of a Sophy [®] valve BEFORE IMPLANTATION: selection of the operating pressure	10
9.3.	Implantation Technique	11
9.4.	Reading and/or changing the operating pressure of a Sophy [®] valve AFTER IMPLANTATION	12
9.5.	Post-operative X-ray control: Identification of the valve model and pressure reading	14
9.6.	Patency Test (post-operative)	15
9.7.	Sampling the CSF and injection	15
9.8.	Potential Causes of Adjustment Difficulties and Suggested Measures	15
10.	Precautions for the Daily Life of the Patient	16
11.	Complications / Side effects	17
11.1.	Obstruction	17
11.2.	Infection	17
11.3.	Overdrainage	17
11.4.	Complications specific to lumboperitoneal shunts	17
11.5.	Other	18
12.	Storage conditions	18
13.	Processing of the products after use	18
13.1.	Destruction after use	18
13.2.	Return of products	18
14.	Warranty	18
15.	Symbols	18
16.	References	20

FIGURE 1 – ADJUSTABLE SOPHY® MINI VALVE (MODELS SM8, SM8A, AND SM8B). TOP AND SIDE VIEWS.




FIGURE 2A - COMPOSITION OF A COMPLETE ADJUSTABLE SOPHY® MINI KIT (MODEL SM8-2010).

- A** SM8 valve with pre-connected peritoneal catheter.
- B** Right angle adapter.
- C** Straight ventricular catheter.
- D** Introducing stylet.

FIGURE 2B - COMPOSITION OF A COMPLETE ADJUSTABLE SOPHY® MINI KIT (MODEL SM8-2040).

- A Lumbar catheter.
- B Asymmetrical metal connector.
- C Intermediate catheter.
- D SM8 adjustable Sophy® Mini valve.
- E Peritoneal catheter.
- F 14G Tuohy needle.

WARNING

Federal (USA) Law restricts this device to sale by or on the order of a physician.

Read the Instructions for Use carefully before adjusting the valve operating pressure.

1. Indications

The adjustable Sophy® valve is intended for the treatment of hydrocephalus by shunting the Cerebrospinal Fluid (CSF) to the abdominal cavity or the right atrium of the heart.

2. Contraindications

The contra-indications are the following :

- established or suspected infections along the length of the shunt (meningitis, ventriculitis, peritonitis, septicemia or bacteremia) or any infection present in any part of the body whatsoever,
- patients on anticoagulant therapy, or presenting with bleeding diathesis, or with hemorrhagic CSF, as the presence of blood in the system could lead to an obstruction in the system),
- ventriculo-atrial shunts in patients suffering from congenital cardiopathies or other malformations of the cardio-pulmonary system.


However, if the clinical benefit outweighs these contraindications, the implantation of a shunt in such cases is performed under the responsibility of an experienced neurosurgeon. The clinical condition of the patient must, therefore, be subject to increased surveillance.

WARNING

Do not use an external shunting device (drainage bag, etc.) in series with a valve as the two systems could interfere with each other and disrupt control of the drainage.

3. Description and Operating Principle of the Sophy® valve

FIGURE 3: THE SOPHY® VALVE

The Sophy® adjustable pressure valve (*Figures 1&3*) is a single use implantable device.

The Sophy® adjustable pressure valve allows single direction drainage of the Cerebrospinal Fluid (CSF).

The CSF arrives in the valve through the inlet connector [1], passes through the body of the valve [5] and leaves through the outlet connector [6].

The connectors are made of stainless steel and the body of the valve is in polysulfone.

On either side of the inlet connector [1] are suture holes [2] which make it possible to attach the valve to the sub-cutaneous tissues to prevent the valve from migrating.

On the upper surface of the valve an arrow [4] shows the direction of CSF flow through the valve. This helps to position the valve correctly during implantation.

On the lower surface of the valve there is an individual serial number.

The valve body contains a ball-in-cone mechanism which determines the operating pressure of the valve.

In conditions of normal use this mechanism provides an anti-reflux function and is not sensitive to temperature variations.

The valve body, which cannot be deformed, protects the mechanism from mechanical shocks.

It also prevents attempts at pumping or puncturing the valve and makes it insensitive to variations in percutaneous pressure.

The Sophy® adjustable valve is phthalate-free and does not contain natural or synthetic latex.

3.1. ADJUSTABLE PRESSURE

The resistance of the Sophy® valve can be adjusted without a further operation, in order to adapt it to the clinical

evolution of the patient.

The principle is based on the pressure variation exerted on a ruby ball [7] by a flat, semi-circular spring [8] at different points on its curvature.

The spring is attached to a rotor [9] which is able to pivot inside the valve body around its central axis [11].

A notched indexing system makes it possible to select different positions corresponding to different operating pressures. Two safety stops included in the valve body limit any possibility of rotation beyond the pre-determined positions.

The operating pressure of the Sophy® valve is determined by the angular position of its rotor.

3.1.1. 8 pressures available for SM8

Eight pressures are available for each SM8 Sophy® model. These correspond to the eight possible rotor positions.

For each SM8 Sophy® model, position No. 1 corresponds to the lowest pressure and position No. 8 to the highest pressure.

3.1.2. 3 pressures available for SM3

Three pressures are available for each SM3 Sophy® model. These correspond to the three possible rotor positions.

For each SM3 Sophy® model, position No. 1 corresponds to the lowest pressure and position No. 3 to the highest pressure.

3.1.3. Pressure ranges

Titanium radio-opaque studs [10] inserted on the right-hand side of the valve body indicate the 8 or 3 positions of the rotor/spring assembly (See §9.5 - Post-operative X-ray control).

The SM8 Sophy® valve exists in 4 different pressure ranges (see Table of pressure ranges).

To the left of the inlet connector there are radio-opaque points [3] which are used to identify the pressure range of the valve (Figures 3 & 9, and §9.5 - Post-operative X-ray control).

Table of Operating Pressures for SM8 valves

Note on the Locator	SM8-140	SM8	SM8-300	SM8-400	
References concerned	SM8-140	SM8, SM8-2010 SM8A, SM8A-2010 SM8B, SM8B-2010 SM8-2020, SM8-2021 SM8-2030, SM8-2031 SM8-2040	SM8-300	SM8-400	
Identification of the range	0 point	1 point	2 points	3 points	
Available pressures (mmH ₂ O)	1	10	30	50	80
	2	25	50	75	120
	3	40	70	100	150
	4	60	90	125	190
	5	80	110	150	230
	6	100	140	180	270
	7	120	170	220	330
	8	140	200	300	400

Table of Operating Pressures for the SM3 valve

Note on the Locator	SM3	
References concerned	SM3	
Available pressures (mmH ₂ O)	1	50
	2	110
	3	170

4. Principle of the Sophy® valve adjustment

For more information, refer to the instruction for use of the Adjustment Kit.

The operating pressure of a valve can be adjusted in order to adapt to the clinical evolution of the patient.

Changing the operating pressure of the Sophy® valve is carried out percutaneously using an adjustment kit made up of 3 parts:

- Locating Instrument (Locator),
- Reading Instrument (Compass),
- Setting Instrument (Magnet).

The components of the Sophy® adjustment kit are packaged non-sterile in a re-usable box with magnetic shielding.

The Sophy® adjustment kit is not made with natural or synthetic latex.

The Sophy® adjustment kit is specifically for setting adjustable Sophy® valve pressures.

4.1. LOCATING INSTRUMENT

The Locating Instrument (**Locator**) is used to locate the valve under the skin, an essential step in reading and adjusting the operating pressure.

The Locator houses the Reading Instrument and the Adjustment Instrument in turn.

It enables Sophy® valves to be adjusted whatever their pressure range (see Table of Operating Pressures).

Two locators are available:

- SAK-LI1 for the adjustment of SM8 valves
- SAK-LI2 for the adjustment of SM3 valves

WARNINGS

Do not use a Sophy® adjustment kit without having identified the valve model and making sure that the **Locator** and the pressure range shown on the **Locator** corresponds to this model.

Over-drainage or under-drainage can result from taking a reading and/or an adjustment with a pressure range on the **Locator** different from that of the implanted valve model.

To adjust a valve the **Locator** must display the pressure range for the valve model to be adjusted. The **Locator** must be perfectly positioned above the implant.

The red arrow shows the direction of the CSF flow, and makes it possible to position the **Locator** correctly in relation to the valve.

The central cut-out in the shape of an imprint of the valve is used to position the **Locator** as close as possible to the implant.

4.2. READING INSTRUMENT

The Reading Instrument (**Compass**) fits into the **Locator** and enables a reliable, accurate direct reading of the operating pressure value for the valve.

PRECAUTION

When reading the pressure setting, make sure that the **Magnet** or any other ferromagnetic object is not less than 0.5 m away from the compass so that this reading cannot be falsified by the influence of magnetic fields.

The **Compass** shows the position of the valve rotor. This corresponds to the operating pressure read.

The operating pressure reading is determined by the alignment of the **Compass** needle with one of the lines on the contour of the **Compass** and the corresponding pressure value in the reading area of the **Locator**.

4.3. SETTING INSTRUMENT

The Setting Instrument (**Magnet**) is made up of a specific assembly of magnets designed to enable the operating pressure of the Sophy® valve to be adjusted.

PRECAUTION

Changing the valve pressure must only be done by a neurosurgeon, or other qualified person.

The **Magnet** fits into the **Locator**.

The insertion of the **Magnet** into the **Locator** makes it possible:

- to adjust the operating pressure by rotating the **Magnet**.

A red marker pointing to the pressure values in the **Locator** reading area and the “clicks” produced when the **Magnet** rotates make it possible to determine the successive operating pressures set by the **Magnet**.

Clockwise rotation increases the operating pressure of the valve, while counter-clockwise rotation reduces it.

PRECAUTIONS

Do not store or handle the **Magnet** near any device likely to be altered by a powerful magnetic field.

Do not move the **Magnet** close to a powerful magnetic field (e.g. MRI). It could become a dangerous projectile or its performance could be affected.

Handle the **Magnet** with care near metal objects (e.g. hospital furniture). They could be violently attracted.

4.4. DEMONSTRATION VALVE

A demonstration valve is supplied in the adjustment kit to provide training and display of the changes in the position of the rotor during the adjustment procedure.

5. Configurations of the adjustable Sophy® valve

The adjustable Sophy® valve is available in 4 different pressure ranges, each of which offers 8 pressures (see §3 - Table of Operating Pressures).

The adjustable Sophy® valves in the “standard” pressure ranges (30-200 mmH₂O) are available in 3 models:

- without a reservoir (SM8),
- with an integral reservoir of the antechamber type (SM8A),
- with an integral reservoir of the burr hole type (SM8B).

These 3 models are available in the form of:

- valve only,
- complete kit, including a proximal catheter (straight or elbow) and a pre-connected peritoneal catheter (Figures 2a & 2b).

The adjustable Sophy® valves in the “special” pressure ranges (10-140, 50-300, and 80-400 mmH₂O) are available without a reservoir (SM8-140, SM8-300 and SM8-400).

These “special pressure” valves are only available in the valve only form and are not available as complete kits or with a reservoir.

Sophysa offers a complete range of ventricular and distal radio-opaque catheters which allow the CSF to flow to the valve and from the valve to the peritoneum or right atrium respectively, depending upon the type of shunt chosen by the neurosurgeon.

To be complete an adjustable Sophy® shunt system must consist of a ventricular catheter, an adjustable Sophy® valve and a distal catheter (atrial or peritoneal).

Each adjustable Sophy® valve is packaged with a Patient Identification Card (PIC) and traceability labels.

The surgeon is responsible for filling in this Patient Identification Card and giving it back to the patient.

6. Measurement Unit and Calibration of Operating Pressures

The pressures mentioned are in mmH₂O.

1 mmH₂O corresponds to 9.807 Pa or 0.074 mmHg.

The valves are calibrated on the basis of a flow rate of 10 ml/h.

Each valve is tested individually: the measurement concerns the upstream pressure of a 10 ml/h flow of water passing through the valve and the Sophysa proximal and distal catheters.

The calibration is performed disregarding the resistance of the catheters.


Thus the pressures given on the valve labels correspond to the resistance of the valve alone.

The catheters add their own resistance to the shunt.

The reservoirs are not considered to have any particular resistance.

The calibration of Sophy® valves is performed with a tolerance of -10/+15 mmH₂O on the measured pressures.

FIGURE 4: FLOW RATE-PRESSURE CURVES FOR THE SM8 MODEL

This curve is obtained by varying the applied pressure for each pressure setting and measuring the resulting flow rate. The values are given disregarding the resistance of the catheters.

The effects on the device operating pressure of changes in the position of the patient and of sub-cutaneous pressure are negligible.

7. Behavior during Magnetic Resonance Imaging (MRI)

A shunt consisting of a Sophy® valve (including connectors and possible reservoirs) and its catheters, is considered to be “MR Conditional” in accordance with the definition in the standard, ASTM F2503.

A patient fitted with a Sophy® valve can undergo an MRI examination, even immediately after the device has been implanted.

The results of in vitro tests have demonstrated that the Sophy® valve does not present any danger in the following examination conditions:

- MRI with a static magnetic field limited to 3-Tesla and

- with a spatial gradient magnetic field limited to 19 T/m;
- Whole body averaged SAR (Specific Absorption Rate) limited to 2 W/kg for 15 minutes exposure at 3 Tesla.
- Whole body averaged SAR (Specific Absorption Rate) limited to 2 W/kg for 15 minutes exposure at 1.5 Tesla.
- No gradient magnetic fields limitation.

NOTE

The patient must be informed that he/she is likely to feel slight discomfort, which is completely inoffensive, during an MRI examination.

PRECAUTIONS

For the Sophy® valve choose an implantation site away from an area of specific clinical interest, such as a tumor, which may subsequently require repeated MRI examinations. In fact, the Sophy® valve micro-magnets are a potential source of artifacts on MRI images. The size of these artifacts could be very large in size in relation to the size and shape of the Valve.

WARNING

The adjustment of the pressure must be checked routinely before and after MRI examinations. Please see additional recommendation under section §9.4.1-Reading of operating pressure after implantation.

The tests conducted in accordance with the standard, ASTM F2182, showed that the rise in temperature caused by exposure to 3-Tesla MRI (+ 2.2°C +/- 0.2°C) was negligible and had no physiological consequences for the patient.

Tests conducted in accordance with standards, ASTM F2213 and ASTM F2052, showed that the torque and displacement force induced by a magnetic field of 3 Tesla or less did not present any risks for the patient.

The tests conducted in accordance with the standard, ASTM F2119-07, showed that image artifacts may reach a surface of 44 cm² for Spin Echo sequences and 128 cm² for Gradient Echo sequences on a 3 Tesla MRI.

8. Sterilization – Decontamination of Sophy® Valves and Valve Kits

The Sophy® valves and valve kits are packed individually in double peel-off, sterile, pyrogen-free packaging, sterilized with ethylene oxide.

WARNING

Do not use the valves or valve kits if the sterile packaging is open or damaged, or if the expiry date has passed.

This product is intended for SINGLE USE ONLY. It is intended to be used once only for a single patient. Do not re-sterilize or re-use after unpacking and/or explantation.

Resterilization can damage the product, potentially leading to patient injury. Reuse of this device may change its mechanical or biological features and may cause device failure, allergic reactions or bacterial infections.

NOTE

Sophysa cannot be held responsible for the performance of any product that has been re-sterilized and/or re-used, or for any complications which might result from this.

9. Instructions

9.1. CHOICE OF VALVE MODEL

The choice of the Sophy[®] valve model is left to the initiative of the neurosurgeon and depends upon the clinical needs of the patient.

PRECAUTIONS

Do not use pre-connected valve kits (valves with pre-attached distal catheter) for ventriculo-atrial shunts.

There could be implantation difficulties relating to the shunt length adaptation at the atrium.

Use an SM8A (antechamber), an SM8B (burr hole reservoir) model, or one of the SM8 models combined with a ventricular catheter with reservoir if the use of the shunt system is desired to check the patency of the shunt, to sample the CSF or for injections.


9.2. ADJUSTMENT OF A SOPHY[®] VALVE BEFORE IMPLANTATION: SELECTION OF THE OPERATING PRESSURE

Refer to §4 – *Principle of the Sophy[®] valve adjustment*. For more information refer to the instructions for the adjustment kit.

Before opening the packaging which guarantees the sterility of the valve, it is necessary to adjust the rotor in the position that corresponds to the pressure selected for the implantation by the surgeon, depending on the patient's needs.

The double sterile packaging for each Sophy[®] valve includes a housing provided specially for the components of the Sophy[®] adjustment kit, which makes it possible to adjust the Sophy[®] valve through the packaging.

FIGURE 5: DETAILS OF THE DOUBLE STERILE PACKAGING FOR A SOPHY[®] VALVE (SM8)

PRECAUTION

Do not unpack or implant the valve without previously adjusting it to the initial operating pressure defined by the neurosurgeon.

WARNING

Do not use the Sophy[®] adjustment kit in the operative field. As the Sophy[®] adjustment kit cannot be sterilized, using it during the operation would cause a high risk of infection for the patient.

9.2.1. Identification of the valve model (Step A)

Identify the reference and the pressure range for the Sophy[®] valve to be adjusted: refer to the labeling on the packaging.

Confirm that the valve model and the **Locator** match.

9.2.2. Positioning the Compass and Reading the Pressure (Step B)

If necessary, refer to §4 for the description of the Reading Instrument.

Place the **Compass** in the **shell of the external packaging** in the place provided for it, centered above the valve.

On the **shell of the external packaging** read the operating pressure value shown by the **Compass** needle.

NOTE

Reading of the setting position can also be performed visually without using the Sophy[®] adjustment kit, as the rotor is visible through the valve and the double sterile packaging.

9.2.3. Positioning the Magnet and Adjusting the Pressure (Step C)

If necessary, refer to §4 for the description of the Setting

Instrument.

Remove the **Compass**, align the mark on the **Magnet** with the pressure read in the previous step, and slide the **Magnet** vertically into the shell of the external packaging without changing the orientation of the **Magnet**.

Turn the **Magnet** until it reaches the position chosen as the operating pressure for the valve to be implanted.

Remove the **Magnet** vertically.

NOTE

Do not attempt to turn the **Magnet** if it is not centered on the shell of the external packaging, or if it is not completely inserted into the shell of the external packaging. The rotation of the valve rotor could be compromised.

9.2.4. Verification of the pressure setting (Step D)

Verification is performed using the **Compass** in accordance with the procedure described in Step B, or visually on the valve.

9.2.5. Recording the operating pressure chosen for the implantation (Step E)

Note the pressure value read in the previous step on the Sophy® Patient Identification Card (NT120 or PC-SM8).

9.3. IMPLANTATION TECHNIQUE

Implantation of a Sophy® adjustable valve refers to current aseptic neurosurgical practices.

The implantation of a shunt including a Sophy® valve may be performed in several ways.

The surgeon will choose the technique depending upon his experience and the clinical status of the patient.

The final implantation of the device must satisfy the conditions for optimal drainage of the CSF.

The surgeon must select the implantation area taking into account the fact that the valve is a potential source of artifacts when an MRI examination is performed (see § 7 - *Behaviour during Magnetic Resonance Imaging (MRI)*).

PRECAUTION

Do not perform the implantation of a shunt without having a replacement shunt system available in case it is required.

WARNING

Do not carry out an additional test before implantation: each valve has been individually calibrated and checked. Any pre-operative pressure tests will increase the risk of infection.

9.3.1. Ventricular Catheter

- Introduce the catheter into the ventricle using the introducing stylet supplied for this purpose.
- If necessary, adjust the implantation depth of the ventricular catheter with the right angle adapter supplied. Position it in the axis of the burr hole.
- Purge the catheter of air with the CSF.
- If necessary, check that the reservoir is properly filled, and then clamp.

- Connect and delicately ligate the catheter to the inlet connector of the valve (or that of the reservoir for valve models with integrated reservoirs). Check that the arrow located on the upper surface of the valve is correctly oriented in the direction of the flow. The clamp can then be released.

WARNING

Ensure that the arrow on the upper surface of the valve is correctly oriented in the direction of the flow: assembly of the valve in the opposite direction would prevent any drainage.

PRECAUTION

Due to the fragility of the silicone, using metallic forceps for inserting catheters and for ligating them onto connectors is not recommended. This would create a risk of cutting or piercing the catheters.

9.3.2. Valve

The Sophy® models may be implanted either on the skull or in the sub-clavicular fossa (pectoral implantation). This latter implantation may be preferred in paediatrics to avoid any risk of skin erosion. Thoracic implantation may also be chosen in patients with a pathology of the posterior fossa, who are likely to be subject to MRIs frequently, in order to prevent artifacts on the images.

PRECAUTIONS

Sophy® models must be implanted under sub-cutaneous tissue less than 10 mm thick.

The selected implantation site must be sufficiently unencumbered to facilitate detection of the valve by the **Locator** when adjustments are made later (space of at least 4 cm all around the valve, away from the ear and the neck).

If these precautions are not observed, there may be difficulties with adjustment or it may even be impossible to make an adjustment, due to poor localization of the valve by the **Locator** or the **Magnet** moving the rotor insufficiently.

WARNING

In the case of valves with an integrated reservoir or catheter, do not attempt to detach the reservoir or catheter from the valve. Detachment of the reservoir or catheter may unscrew the connector closure screw and uncalibrate the valve.

1. Purge the valve of air. To prevent any risk of introducing an air bubble, it is recommended that the valve be left to fill directly with the patient's CSF. In the majority of cases, the valve fills immediately.

However, in patients with low intracranial pressure or if the valve is set to a Medium or High pressure, the valve cannot fill spontaneously. In this case :

- place a piece of catheter on the outlet connector and slowly aspirate the CSF using a syringe fitted with a Luer connector,
- or even:
- press the dome of the reservoir to allow the CSF

to fill the valve (for SM8A and SM8B models).

PRECAUTION

Do not fill or purge the valve with any liquid other than the patient's CSF or water for injection (WFI) before implantation to avoid any risk of deposits in the valve, which could lead to an obstruction in the shunt system or a blockage in the valve mechanism.

2. Check that the valve is correctly filled with CSF and there are no air bubbles inside the valve. If this is not the case, continue to purge. The presence of air bubbles could cause a significant change to the operating pressure initially chosen.
3. Check that the arrow on the upper surface of the valve is visible and correctly oriented in the direction of the CSF flow.

PRECAUTION

Do not implant the valve without suturing it to the underlying tissues by its two connectors or by the suture holes provided for this. If the shunt system migrates the drainage may stop and other complications ensue.

WARNING

Before suturing the valve to the underlying tissues, check that the arrow on the upper face of the valve is visible.

If the arrow is not visible this means that the valve has been implanted the wrong way up (upside – down). In this case, direct reading with the **Compass** will be reversed, e.g. reading of a high position (No.8) for a valve set to a low position (No.1). This could cause serious clinical consequences (over- or underdrainage). In this situation, contact Sophysa for adjustment instructions.

9.3.3. Peritoneal Catheter

1. Make a short peritoneal incision in the peri-umbilical region.
2. Tunnel the distal catheter.
3. Connect the proximal end of the catheter to the valve outlet connector and ligate it delicately.
4. Check the flow of CSF.
5. Adapt the length of the catheter.
6. Bury the distal end of the catheter in the peritoneal cavity.

PRECAUTION

Due to the fragility of the silicone, using metallic forceps for inserting catheters and for ligating them onto connectors is not recommended. This would create a risk of cutting or piercing the catheters.

9.3.4. Specific characteristics linked to the implantation of the SM8-2040 pre-connected lumboperitoneal kit

9.3.4.1. Lumbar Catheter

1. Make a small vertical skin incision between two spinous processes at the chosen level.
2. Insert the 14G Tuohy needle into the centre of the incision, push the needle into the lumbar

subarachnoid space using an appropriate technique.

3. Pass the open end of the catheter with multiple perforations through the Tuohy needle. The first graduation placed at 11 cm corresponds to the length of the needle. The following graduations indicate the length of catheter effectively inserted into the lumbar sub-arachnoid space.
4. Check the flow of CSF.
5. Remove the Tuohy needle while leaving the catheter in place.
6. Make an incision in the lower side part of the thorax on the axillary axis.
7. Tunnel the catheter from the lumbar incision to the thoracic incision.
8. Purge the catheter of air with the CSF and then clamp.

WARNING

Avoid implanting a lumboperitoneal shunt in children.

PRECAUTION

Do not use any Tuohy needle other than that supplied in each kit by Sophysa.

9.3.4.2. Valve with pre-connected peritoneal catheter

1. Make a small peritoneal incision in the peri-umbilical region.
2. Form a subcutaneous pocket for the valve from the thoracic incision.
3. Tunnel the peritoneal catheter from the thoracic incision to the peritoneal incision.
4. Adapt the length of the lumbar catheter. Connect it to the intermediate catheter connector. Adapt the length of the intermediate catheter and connect it to the valve.
5. Ligature delicately.
6. Remove the clamp.
7. Purge the valve of air. To prevent any risk of introducing an air bubble, it is recommended that the valve be left to fill directly with the patient's CSF. In the majority of cases, the valve fills immediately.

However, in patients with low intracranial pressure or if the valve is set to a Medium or High pressure, the valve cannot fill spontaneously. In this case place a piece of catheter on the outlet connector and slowly aspirate the CSF using a syringe fitted with a Luer connector.

PRECAUTION

Do not fill or purge the valve with any liquid other than the patient's CSF or water for injection (WFI) before implantation to avoid any risk of deposits in the valve, which could lead to an obstruction in the shunt system or a blockage in the valve mechanism.

9.4. READING AND/OR CHANGING THE OPERATING PRESSURE OF A SOPHY® VALVE AFTER IMPLANTATION

Refer to §4 – Principle of the Sophy® valve adjustment. For more information refer to the instructions for the

adjustment kit.


WARNING

The operation scar and/or post-operative edema can make positioning the **Locator** painful and inaccurate. This could then cause difficulties in adjustment or even make it impossible.

9.4.1. Positioning the patient (Step 0)

For the adjustment, the patient should be positioned so that the valve implantation site is easy to access and the valve is as horizontal as possible.

FIGURE 6 – RECOMMENDED POSITIONING OF THE PATIENT (VALVE HORIZONTAL)

9.4.2. Reading the pressure

9.4.2.1. Identification of the valve model (Step 1)

Identify the reference and the pressure range for the Sophy® valve to be adjusted: refer to the Patient Identification Card or to the patient's medical record (traceability label) and/or perform an X-ray control of the valve (see §9.5. Post-operative X-ray control).

Confirm that the valve model and the **Locator** match.

9.4.2.2. Display of the pressure range on the Locator (Step 2)

If necessary, refer to §4 for the description of the Locating Instrument (**Locator**)

Check that the pressure range displayed in the **Locator** window corresponds to the previously identified valve.

If this is not the case, display the pressure range corresponding to the implanted valve on the **Locator**.

- Hold the **Locator** with one hand.
- Turn the ring until a “click” is heard. This click guarantees that a complete pressure range is displayed.
- Repeat the operation until the range required is found amongst the 4 ranges available.

WARNING

Cases of over-drainage or under-drainage can result from a reading and/or an adjustment with a **Locator** that does not correspond to the type of valve implanted.

9.4.2.3. Positioning the Locator (Step 3)

Palpate the valve implantation site in order to determine both the location and orientation of the valve.

The valve inlet and outlet connectors at each end of the valve and the reservoir (if fitted) are the easiest components to locate. Find these first.

Place the **Locator** on the implantation site with its axis aligned with that of the valve connectors, and with the arrow of the **Locator** pointing in the direction of the CSF flow.

Centre the **Locator** above the valve as well as possible: locate the valve by palpating it through the central cut-out of the **Locator**.

Press the **Locator** firmly against the valve so as to centre and immobilize the valve.

9.4.2.4. Positioning the Compass and Reading the Pressure (Step 4)

If necessary, refer to §4 for the description of the Reading Instrument

Hold the **Locator** with one hand.

Place the **Compass** in the **Locator** using the guide pins: align the lines on the **Compass** with the pressure values engraved on the **Locator**.

To ensure the **Locator** is perfectly centred above the valve, centre the shaft of the **Compass** needle in the circular target by planar movements of the **Locator-Compass** assembly parallel to the valve.

If the **Compass** needle is aligned between two positions, turn the **Locator** slightly with the **Compass** so as to align the **Compass** needle on the nearest position (if the needle is exactly between two positions, choose one and align the needle on it; the real orientation of the **Locator** will be adjusted during Step 5).

On the **Locator** read the operating pressure shown by the **Compass** needle.

9.4.2.5. Confirmation of the reading

In the recommended implantation conditions, X-ray control is optional because the pressure can be read directly using the adjustment kit. However, X-ray confirmation is recommended in the following cases:

- if there is a discrepancy between the pressure read in Step 4b and the value read in Step 1b on the Patient Identification Card (PC-SM8) and/or on the X-ray taken to identify the valve model,
- if the valve has been implanted too deeply, under subcutaneous tissue more than 10 mm thick (see §9.3.2 - Implantation technique - Valve),
- if the user is not familiar with the use of the Sophy® adjustment kit
- after an MRI examination.

9.4.3. Setting a new pressure

After reading the pressure level, as described above, keep the **Locator** in place to perform the valve adjustment.

9.4.3.1. Orientation of the Locator – "Extreme position" method (Step 5)

If a difference is observed between the pressure value read in the previous step and the patient monitoring documents (Patient Identification Card, X-ray films, etc.), or if the **Compass** needle pointed between two positions, it is recommended that the angular positioning of the **Locator** in relation to the valve axis is refined. To do this, use the "extreme position" method described below:

- Without moving the **Locator** remove the **Compass**, align the mark on the **Magnet** on the pressure read to insert the **Magnet** vertically in the **Locator**.
- Turn the **Magnet** until it just passes the extreme position the furthest away (No. 1 or No. 8).
- Without moving the **Locator**, remove the **Magnet** vertically.
- Place the **Magnet** at least 0.5 m from the valve.
- Insert the **Compass** into the **Locator**.

If the **Locator** is correctly oriented and if the adjustment is correctly performed, the **Compass** will confirm that the extreme position has been reached.

If the needle does not point to this extreme pressure value:

- Turn the **Locator-Compass** assembly until the needle is perfectly aligned with the extreme pressure value selected.
- If necessary recentre the axis in the target.

9.4.3.2. Setting the pressure (Step 6)

If necessary, refer to §4 for the description of the Setting Instrument

Up to the last adjustment step (end of Step 7b), hold the **Locator** firmly in position with one hand.

Remove the **Compass** memorising the pressure it indicates. Insert the **Magnet** aligning the mark on the memorised position.

Turn the **Magnet** until it reaches the position selected as the new operating pressure for the valve.

NOTE

The design of the valve makes adjustment easier when it is performed clockwise. Consequently, when the **Locator** orientation step has been performed on Position No.8 a return to Position No.1 is recommended followed by setting the desired pressure.

Confirm the adjustment by a brief press on the **Magnet**. Remove the **Magnet** vertically.

9.4.3.3. Verification of the pressure setting (Step 7)

If necessary, refer to §4 for the description of the Reading Instrument

Verification is performed using the **Compass** in accordance with the procedure described in Step 4.

NOTE

Setting the selected operating pressure cannot always be performed at the first attempt. It is possible that the procedure will have to be repeated once or several times, in order to ensure that the desired pressure setting is obtained.

Paragraph 9.8 is devoted to the description of measures to be taken if adjustment is difficult.

PRECAUTION

Adjusting the pressure using a Sophy® adjustment kit requires experience. If the user is not familiar with the use of the kit, it is then recommended that the pressure adjustment made is checked with an X-ray.

9.4.3.4. Recording the new operating pressure (Step 8)

Record the pressure value read during verification (Step 7) on the Sophy® Patient Identification Card (NT120 or PC-SM8).

9.4.3.5. Monitoring after an adjustment

It is advisable to monitor the patient carefully for 24 hours following any change in the valve operating pressure setting.

9.5. POST-OPERATIVE X-RAY CONTROL: IDENTIFICATION OF THE VALVE MODEL AND PRESSURE READING

NOTE

During the radiological examination orientate the patient so that the ray source points perpendicularly onto the valve body.

In this way identifying the valve by its radio-opaque point is made easy.


The radio-opaque points to the left of the inlet connector are used to identify the pressure range of the SM8 Sophy® valve:

Number of radio-opaque points	zero	one	two	three
		•	••	•••
Maximum operating pressure (mmH ₂ O)	140	200	300	400

NOTE


The pressure range of the implanted valve is shown on the Patient Card, in the description of the valve model used.

FIGURE 7 – RADIOGRAPHIC IDENTIFICATION OF THE DIFFERENT MODELS OF THE SM8 SOPHY® VALVE.

It is easy to read the pressures on a radiological image (Figure 8).

FIGURE 8 – CONFIRMATION OF THE OPERATING PRESSURE (X-RAY FILMS OF THE SM8 VALVE MODEL)

On the Sophy® valve the main pressures of the valve are located by 1 radio-opaque point for Low pressure, 2 radio-opaque points for Medium pressure and 3 radio-opaque points for High pressure. Intermediate pressures are determined by the position of the rotor in relation to these radio-opaque markings for the main pressures.

Locate the valve inlet connector, wider due to the presence of a nut.

The position for the lowest operating pressure (Position No. 1) is the position located nearest to this connector.

If the right-hand valve rotor micro-magnet is facing the radio-opaque point nearest the inlet connector, the valve is set at the lowest pressure (Position No. 1).

Then, moving clockwise away from the inlet connector, each of the following positions corresponds to a higher pressure.

The radio-opaque marking furthest away from the valve inlet connector corresponds to the highest pressure (Position No. 8).

9.6. PATENCY TEST (POST-OPERATIVE)

There are two steps for the post-operative test on the patency of the shunt :

9.6.1. Testing the patency of the ventricular catheter

NOTE

This control is possible with SM8A (antechamber) and SM8B (burr hole reservoir) type models. For the SM8, SM3 type models, a ventricular catheter with reservoir must be used.

Pinch the catheter with a finger just after the valve outlet connector.

With another finger, press the reservoir to make the CSF flow back into the ventricular catheter. A reservoir that cannot be compressed easily or does not fill quickly may indicate there is an obstruction in the ventricular catheter.

9.6.2. Patency test downstream of the reservoir (valve and distal catheter)

NOTE

This control is impossible with a SM8B type model (burr hole reservoir) because there is no access to the ventricular catheter upstream of the reservoir.

Pinch the catheter with a finger just before the reservoir, then with another finger press the reservoir to push the CSF through the valve and distal catheter. A reservoir that cannot be compressed easily may indicate an obstruction either of the valve or the distal catheter.

PRECAUTION

Do not rely only on the characteristics of the patency test to diagnose an obstruction in the shunt system. Obstruction of a shunt system can occur in any of its components and should be diagnosed first of all by the clinical data and additional examinations.

9.7. SAMPLING THE CSF AND INJECTION

Access to the CSF is obtained by pricking the reservoir with a 24G (or smaller diameter) Huber needle.

The integral reservoir on the SM8A and SM8B models is designed for occasional use.

Its watertight performance is reduced after very frequent pricking into the dome.

- To inject in the proximal direction, compress the catheter just after the valve outlet connector.
- To inject in the distal direction, compress the catheter upstream of the reservoir.

NOTE

Elective injection in the distal direction is not possible with a model of the SM8B type (burr hole reservoir) because there is no access to the ventricular catheter upstream of the reservoir.

Ensure that the base of the reservoir is not crossed with the needle.

PRECAUTIONS

Do not inject into, or take samples from, the CSF without having tested the shunt patency. Significant overpressure could damage the shunt if it is known that there is an obstruction.

Do not inject too fast or inject too great a volume. The increase in pressure could damage the shunt.

Do not use a syringe with a volume of less than 10cc for injections or taking samples. Too great a pressure could damage the shunt.

9.8. POTENTIAL CAUSES OF ADJUSTMENT DIFFICULTIES AND SUGGESTED MEASURES

Poor positioning of the **Locator** in relation to the centre of the valve, or in relation to the **valve's** flow path, or in relation to the flat surface/plane of the **valve**, may lead to inaccurate reading of the operating pressure by the **Compass**, and/or to difficulty adjusting with the **Magnet**.

Before doing anything else check that the **Locator** is correctly positioned and oriented above the valve (see *Steps 0* (Positioning of the Patient), *3* (Positioning the Locator) and *5* (Orientation of the Locator) in §9.4 - *Setting the operating pressure after implantation*).

NOTE

Any off-centring of more than 2 mm or any inclination in relation to the valve could cause an incorrect pressure reading making its rotation impossible, and so making it impossible to change the pressure.

However, there may be other causes for difficulties in adjustment:

- Excessively thick sub-cutaneous tissue above the valve (10 mm) caused, for example, by post-operative edema, a hematoma or cicatricial tissue would make the localization of the valve inaccurate, thus causing complications in the pressure adjustment procedure.
- The presence of debris or deposits inside the valve could also make adjustment using the usual procedure difficult.

Special procedures may then be applied to facilitate the adjustment of the valve, possibly using fluoroscopic control.

9.8.1. Adjustment without the Locator

If the implantation is deep, the standard procedure described above may be adapted so that the use of the **Locator** is no longer needed when adjusting a new pressure. The **Magnet** is thus closer to the valve.

NOTE

This technique without the Locator concerns pressure Adjustment only. Reading and confirmation of the pressure level must be done in accordance with the standard procedure and/or by X-ray control.

1. Follow the usual procedure (Steps 1 to 4) to position the **Locator** using the **Compass**.
2. On the skin locate the axis corresponding to the current pressure of the valve (indicated by the **Compass** needle) and the axis corresponding to the desired pressure (known with the **Locator**).
3. Remove the **Locator-Compass** assembly and place the **Magnet** in the same place, directly on the implantation site, oriented in the direction corresponding to the current pressure located in the previous step.
4. Turn the **Magnet** until it reaches the orientation corresponding to the desired pressure.
5. Confirm the adjustment by a brief press on the **Magnet**.
6. Remove the **Magnet** vertically.
7. Check the setting using the **Compass** and **Locator** or with an X-ray.
8. If the pressure has not changed repeat the adjustment holding the **Magnet** in the bottom position (hold down continuously) during rotation. Then check the new setting using the **Compass** and **Locator** or with an X-ray.

9.8.2. Adjustment of valves implanted in

conditions other than those recommended

The implantation of a valve not performed under the recommended conditions (see §9 – *Implantation technique*) may lead to a **Compass** reading of a pressure value inconsistent with the patient record or the clinical status of the patient.

In this case, an X-ray removes any doubt. It is the absolute proof of correct adjustment and the correct direction of the implantation.

9.8.3. Specific case of valves implanted upside down.

Upside-down implantation (but with the flow direction respected) will be recognizable on X-ray :

With valve correctly implanted, the 5 radio-opaque dots must be located near the nose if the valve is implanted on the patient's right side, and opposite the nose if the valve is implanted on the patient's left side.

If the valve is implanted upside down, setting is performed using the following sequence:

1. Position the **Locator** in the opposite direction to the flow (red arrow towards the inlet connector).
2. Centre the **Locator-Compass** assembly using the **Compass** target (and not the shape of the Locator).
3. Carry out the reading and then the adjustment, in accordance with Steps 4 and 6, using the opposite indication on the **Locator** (e.g. the indication 65 mmH₂O corresponds to a value of 150 mmH₂O).
4. Verify the new setting with an X-ray.

10. Precautions for the Daily Life of the Patient

A Patient Identification Card (PC-SM8) is supplied with the Sophy® valve. It enables the neurosurgeon to consult and update information relating to the implanted device (reference, operating pressure, implantation site, etc.) systematically and to ensure that the illness is properly monitored.

PRECAUTION

The patient should be warned that it is important to carry his/her Patient Identification Card (PC-SM8) at all times. This card gives information on the medical situation of the patient to all medical personnel.

The doctor is responsible for informing the patient or his/her family that any activity exposing the patient to a powerful magnetic field may cause the device settings to be disrupted.

Magnetic fields generated by walk-through detectors in airports, microwave ovens, cordless telephones, high tension cables, and TV etc. must not affect the operating pressure.

Permanent household magnets, in some toys, audio headphones, loudspeakers as well as magnetic fields created by electric motors in razors, hairdryers, clippers, etc. may upset the valve adjustment when they are in operation. It is therefore recommended that they are not allowed too close to the valve implantation site to prevent any risk of changing the valve pressure.

The doctor is responsible for informing the patient or his/her family that the person fitted with a shunt must avoid any activity that may subject this shunt to direct shocks (violent sports, etc.) as these are likely to damage it.

WARNING

The pressure settings must be checked routinely before and after exposure to a powerful magnetic field or if there is an impact at the implantation site.

PRECAUTION

The patient must be warned that vibrations due to the CSF flow may possibly be felt because of the implantation of the valve on the skull.

11. Complications / Side effects

Complications which may result from the implantation of a CSF shunt system include the inherent risks in the use of drugs, any surgical intervention and the insertion of a foreign body.

PRECAUTION

Patients treated with a shunt system must be closely monitored post-operatively in order to detect any signs of complications early.

The doctor is responsible for educating the patient or his/her family about CSF shunt systems, in particular describing the complications linked to implanted shunt systems as well as giving explanations about possible alternative therapies.

The main complications of shunts are obstruction, infection and over-drainage. These complications require the rapid intervention of a doctor.

11.1. OBSTRUCTION

Obstruction is the most frequent complication in shunt systems. It can occur at any point in a shunt.

The ventricular catheter can be obstructed by a blood clot, cerebral tissue or even tumoral cells.

The end of the ventricular catheter can also become embedded in the choroid plexus or in the ventricular wall, either directly or following a collapse of the walls, a consequence of over-drainage.

The cardiac catheter can be colonized by a thrombus while the appearance of a clot around the catheter could cause an embolism in the pulmonary circulation.

The peritoneal catheter may become obstructed by the peritoneum or by intestinal loops.

Loss of patency in a shunt may also be the result of an obstruction by fragments of cerebral tissue or by biological deposits (protein deposits, etc.).

Obstruction of the shunt will quickly result in the reappearance of the signs and symptoms of intracranial hypertension.

These signs and symptoms vary from patient to patient and over time.

In infants and young children, the symptoms may be an abnormal increase in the size of the skull, a bulge in the

fontanelles, dilation of the scalp veins, vomiting, irritability with a lack of attention, downward deviation of the eyes, and sometimes convulsions.

In older children and adults, intracranial hypertension due to hydrocephalus may be the cause of headaches, vomiting, blurred vision, diplopia, drowsiness, slowing of movements, gait disorders or psychomotor slowing which could lead to total invalidity.

If an obstruction is confirmed and a patency test does not make it possible to reduce the obstruction, revision surgery or removal of the device must be envisaged.

11.2. INFECTION

Chronic malfunction of the shunt could cause a leak and a discharge of CSF along its length increasing the risk of infection.

Local or systemic infection is another possible complication of CSF shunt systems. It is generally secondary to the colonization of the shunt by cutaneous germs. Nevertheless, as for all foreign bodies, any local or systemic infection can colonize the shunt. Erythema, edema and skin erosions along the length of the shunt may be an indication of an infection of the shunt system.

Prolonged, unexplained fever may also be the result of a shunt system infection.

Septicemia, favored by an alteration in general status, can start from a shunt infection.

If there is infection, removal of the system is indicated in conjunction with the start of a specific treatment by a general or intrathecal route.

11.3. OVERDRAINAGE

Overdrainage can result in a collapse of the ventricles (slit ventricle syndrome) and the appearance of a subdural hematoma.

In children, depression of the fontanelles, overlapping of the scalp bones, even a craniostenosis or a change from communicating hydrocephalus to obstructive hydrocephalus by stenosis of the Aqueduct of Sylvius could occur.

Adults can present with a variety of symptoms such as vomiting, auditory or visual disorders, drowsiness or even headaches in the upright position but which improve in the supine position.

Depending on clinical observations and medical imaging, the doctor can reduce the symptoms of over-drainage and correct the ventricle size by changing the operating pressure of the Sophy[®] valve.

However, immediate drainage of a subdural hematoma may be indicated.

11.4. COMPLICATIONS SPECIFIC TO LUMBOPERITONEAL SHUNTS

The lumbar catheter may migrate in the sub-arachnoid space. The peritoneal catheter may migrate in the peritoneal cavity under the action of intestinal peristaltic waves. Perforation or occlusion of abdominal viscera by the peritoneal catheter could occur.

Growth of the body may progressively cause the catheters

to come out of their insertion sites.

These malfunctions require the shunt to be put back in place immediately.

Cases of tonsillar hernias linked to a lumboperitoneal shunt have been reported, particularly in children. Lumboperitoneal shunts are not advised in children.

11.5. OTHER

Failure of a shunt system may also be linked to disconnection of its various components.

The ventricular catheter may migrate inside the ventricle. The peritoneal catheter may migrate into the peritoneal cavity under the action of the peristaltic waves of the intestine, while an atrial catheter may migrate into the right-hand cavities of the heart following the blood flow.

Perforation or occlusion of abdominal viscera by the peritoneal catheter could occur.

Growth of the body may progressively cause the catheters to exit their insertion sites.

These malfunctions require the shunt to be repositioned immediately.

Cases of cutaneous necrosis over the implantation site are possible.

Over time fibrous adhesions may fix the ventricular catheter in the choroid plexus or the cerebral tissue. If removal is being considered, gentle rotation of the catheter about its axis may make it possible to free it. The catheter should never be withdrawn forcibly. If it cannot be taken out without forcing, it is preferable to leave it in place rather than risk an intra-ventricular hemorrhage.

Cases of allergy to silicone have been described.

Cases of epilepsy after implantation of a ventricular shunt have been described.

The ruby ball in the valve can potentially take up an off-center position on its housing due to the presence of a cluster of cells or protein deposits. Among others, such situations can cause :

- loss of regulatory function in the valve potentially increasing the risk of overdrainage.
- an impaired anti-reflux function.

The mobility of the rotor could be impeded by an aggregation of cells or a protein deposit. This could make it impossible to adjust the valve with the Magnet.

Blood clots, cerebral cells or tumoral cells contained in the CSF could lodge in the valve mechanism, which would have the potential to cause changes in the operating characteristics of the valve.

12. Storage conditions

Keep in the original packaging.

Keep in a cool, dry place, away from sun light and heat.

13. Processing of the products after use

13.1. DESTRUCTION AFTER USE

An unpacked, used or explanted Sophy® valve must be destroyed in accordance with the procedures in force in the

medical establishment.

13.2. RETURN OF PRODUCTS

If an explanted valve needs to be returned to Sophysa for analysis, it must be returned immersed in water, indicating if necessary whether cleaning has been performed.

Never use a saline solution likely to form deposits in the valve body which could block the rotor.

In order to assess the returned product properly, it must be accompanied by an explanatory Return to Manufacturer Authorisation form.

14. Warranty

The performance of the Sophy® adjustment kit is only ensured with the Sophy® range of valves, designed, tested and manufactured by Sophysa.

Sophysa guarantees that this medical device is free of any material and manufacturing defects. Apart from this guarantee, Sophysa does not provide any other guarantee, express or implicit, including commercialization or adaptation for a particular use. Sophysa cannot be held responsible for any incident, complication, damage or prejudice occurring directly or indirectly from the use of this device. Sophysa does not authorize anyone whomsoever to take responsibility on its behalf for its products.

The performance of adjustable Sophy® valves is only guaranteed with the range of silicone catheters and accessories designed, tested and manufactured by Sophysa. However, it is possible to use other brands of catheters provided that their internal diameter is identical to that of the catheters recommended by Sophysa.

15. Symbols

REF

Catalog reference

Manufacturer

Date of manufacture

Caution, see the Instructions for Use

STERILE EO

Sterilized using Ethylene Oxide

Do not re-use

Do not re-sterilize

Use by

Batch code

Serial number

CE Conformity Marking

MR Conditionnal

16. References

Sophy® ADJUSTABLE PRESSURE VALVES FOR CSF SHUNTS

SM3*	Sophy® Mini adjustable valve - 50-170 mmH₂O Adjustable pressure valve, 3 pressures: 50 (Low), 110 (Medium), 170 (High) mmH ₂ O.
SM8	Sophy® Mini adjustable valve - 30-200 mmH₂O Adjustable pressure valve, 8 pressures: 30 (Low), 50, 70, 90, 110 (Medium), 140, 170, 200 (High) mmH ₂ O.
SM8-140*	Sophy® Mini adjustable valve - 10-140 mmH₂O Adjustable pressure valve, 8 pressures: 10 (Low), 25, 40, 60, 80 (Medium), 110, 120, 140 (High) mmH ₂ O.
SM8-300*	Sophy® Mini adjustable valve - 50-300 mmH₂O Adjustable pressure valve, 8 pressures: 50 (Low), 75, 100, 125, 150 (Medium), 180, 220, 300 (High) mmH ₂ O.
SM8-400*	Sophy® Mini adjustable valve - 80-400 mmH₂O Adjustable pressure valve, 8 pressures: 80 (Low), 120, 150, 190, 230 (Medium), 270, 330, 400 (High) mmH ₂ O.
SM8A	Sophy® Mini adjustable valve SM8/ Antechamber SPV valve, 8 pressures: 30 (Low), 50, 70, 90, 110 (Medium), 140, 170, 200 (High) mmH ₂ O, with integrated antechamber.
SM8B	Sophy® Mini adjustable valve SM8/Burr hole reservoir SPV valve, 8 pressures: 30 (Low), 50, 70, 90, 110 (Medium), 140, 170, 200 (High) mmH ₂ O, with integrated burr hole reservoir.

SHUNT KITS FOR CRANIAL IMPLANTATION

SM8-2000	Complete Sophy® Mini adjustable kit; 30-200 mmH₂O SM8 valve, 8 pressures: 30 (Low), 50, 70, 90, 110 (Medium), 140, 170, 200 (High) mmH ₂ O with pre-connected peritoneal Catheter (B905S).
SM8-2010	Complete Sophy® Mini adjustable kit; 30-200 mmH₂O SM8 valve, 8 pressures: 30 (Low), 50, 70, 90, 110 (Medium), 140, 170, 200 (High) mmH ₂ O with pre-connected peritoneal Catheter (B905S). Straight ventricular catheter (BO19-10)
SM8-2020	Complete Sophy® Mini adjustable kit; 30-200 mmH₂O SM8 valve, 8 pressures: 30 (Low), 50, 70, 90, 110 (Medium), 140, 170, 200 (High) mmH ₂ O with pre-connected peritoneal Catheter (B905S). Elbow ventricular catheter with reservoir (PRO6)
SM8-2021	Complete Sophy® Mini adjustable kit; 30-200 mmH₂O SM8 valve, 8 pressures: 30 (Low), 50, 70, 90, 110 (Medium), 140, 170, 200 (High) mmH ₂ O with pre-connected peritoneal Catheter (B905S). Elbow ventricular catheter with reservoir (PRO7)
SM8-400-2010*	Complete Sophy® Mini adjustable kit; 80-400 mmH₂O SM8 valve, 8 pressures: 80 (Low), 120, 150, 190, 230 (Medium), 270, 330, 400 (High) mmH ₂ O with pre-connected peritoneal Catheter (B905S). Straight ventricular catheter (BO19-10)
SM8A-2010	Complete Sophy® Mini adjustable kit; 30-200 mmH₂O SM8A valve, 8 pressures: 30 (Low), 50, 70, 90, 110 (Medium), 140, 170, 200 (High) mmH ₂ O with pre-connected peritoneal Catheter (B905S). Straight ventricular catheter (BO19-10)
SM8B-2010	Complete Sophy® Mini adjustable kit ; 30-200 mmH₂O/Burr hole reservoir SM8B valve, 8 pressures: 30 (Low), 50, 70, 90, 110 (Medium), 140, 170, 200 (High) mmH ₂ O with pre-connected peritoneal Catheter (B905S). Straight ventricular catheter (BO19-10)

SHUNT KIT FOR LUMBOPERITONEAL IMPLANTATION

SM8-2040	Complete Sophy® Mini adjustable kit; 30-200 mmH₂O SM8 valve, 8 pressures: 30 (Low), 50, 70, 90, 110 (Medium), 140, 170, 200 (High) mmH ₂ O, with complete set of catheters for the lumboperitoneal shunt (LPS50). Pre-connected peritoneal catheter. Pré-connected intermediate catheter
-----------------	--

SHUNT KITS FOR PECTORAL IMPLANTATION

SM8-2030	Complete Sophy® Mini adjustable kit; 30-200 mmH₂O SM8 valve, 8 pressures: 30 (Low), 50, 70, 90, 110 (Medium), 140, 170, 200 (High) mmH ₂ O with pre-connected peritoneal Catheter (B905S). Elbow ventricular catheter with reservoir (PLO6)
SM8-2031	Complete Sophy® Mini adjustable kit; 30-200 mmH₂O SM8 valve, 8 pressures: 30 (Low), 50, 70, 90, 110 (Medium), 140, 170, 200 (High) mmH ₂ O with pre-connected peritoneal Catheter (B905S). Elbow ventricular catheter with reservoir (PLO7)
SM8-2032	Complete Sophy® Mini adjustable kit; 30-200 mmH₂O SM8 valve, 8 pressures: 30 (Low), 50, 70, 90, 110 (Medium), 140, 170, 200 (High) mmH ₂ O with pre-connected peritoneal Catheter (B905S). Elbow ventricular catheter with reservoir (PLO6)

SM8A-2033

Complete Sophy[®] Mini adjustable kit; 30-200 mmH₂O/Antechamber

SM8A valve, 8 pressures: 30 (Low), 50, 70, 90, 110 (Medium), 140, 170, 200 (High) mmH₂O with pre-connected peritoneal Catheter (B905S). Elbow ventricular catheter with reservoir (PLO5)

DEVICE FOR READING AND ADJUSTING THE OPERATING PRESSURE OF THE Sophy[®] ADJUSTABLE VALVE

SAK

Sophy[®] Adjustable valve adjustment kit


3 elements: SAK-LI1 (SM8), SAK-SI and SAK-RI
with the option of the element: SAK-LI2 (SM3 and SU8)

** Not available for the US market*

Technical specifications and List of product references may be modified without notice. For all product dimensions indicated in this document, tolerances are: $\pm 5\%$.

Availability may vary according to country.

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