Instructions for use

SOPH-A-PORT[®] IMPLANTABLE ACCESS PORT

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
1. Description

Soph-A-Port is a totally implantable device designed to provide repeated access to the vascular system or selected body site to avoid the trauma associated with multiple punctures or the inconvenience of an externalized catheter.

Access to Soph-A-Port is performed by percutaneous needle insertion, using a non-coring needle.

Soph-A-Port consists of a polysulfone port body and a self-sealing silicone septum. Sutures holes are integrated into the port base to facilitate fixation to the underlying fascia. The port catheter is made of radiopaque silicone or polyurethane (See reference table).

A patented self-locking system with a titanium mobile connector and a screw-lock, insures a safe and efficient connection of the catheter to the port (Schema 1).

A serial number is engraved on each port.

SOPHYSA manufactures various types of Soph-A-Port.

They differ in regard to port body size, catheter size and catheter material.

Some Soph-A-Port are available with preattached or unattached catheter configurations. All Soph-A-Port are provided with non-coring needle 22 G.

WARNING

Soph-A-Port is not designed to deliver imaging solutions in high pressure conditions. Should high pressure delivery of such solutions be necessary, on a patient with a Soph-A-Port access port, an alternative pathway must be used.

2. Indications

Soph-A-Port is indicated when patient therapy requires repeated access to the vascular system or selected body site for:

- Medical Injection
- Infusion therapy
- Blood sampling
- Transfusion
- Delivery of imaging solutions
- Parenteral nutrition.

3. Contraindications

Soph-A-Port is not designed, sold or proposed for use in any other indications than above mentioned.

The following contraindications should be observed:

Soph-A-Port should not be prescribed for patients in which the planned drug therapy requires the use of substances known to be incompatible with the materials of construction.

It is recommended to check the compatibility of the drugs with the materials of construction prior to injection.

Soph-A-Port should not be implanted in patient whose body size is too small to support the size of the unit, or in patient suspected to have an allergic reaction to the components material.

Implantation of Soph-A-Port should be avoided in the presence of known or suspected local or general infection, unless access to the vascular system with Soph-A-Port is the solution by which the infection or sepsis has to be treated.

4. SOPH-A-PORT instructions of implantation

Soph-A-Port implantation and catheter placement may be accomplished through a variety of surgical techniques and the physician must choose the appropriate method according to his own experience and patient's need.

As for any foreign body implantation, a strict aseptic technique should be observed.

The port catheter should be positioned at the selected site of therapy and secured by accepted surgical technique to prevent catheter migration.

Position should be confirmed by appropriate X-Ray procedures. Sufficient slack should be left between the catheter insertion point and the port body to preclude strain on the catheter.

The port body should be secured to the underlying fascia in a subcutaneous pocket created by blunt or sharp dissection. Integrated suture holes should be used for fixation.

Care should be taken so that the incision does not cross the septum of the port after closure.

Consideration in choosing the port implant position is advised to avoid locations which interfere with patient

mobility or clothing result in potential pressure point and with associated therapy.

Implantation site should be chosen taking care of the pathology.

For ease in accessing, the port body should be placed over a bony structure (thoracic area). The port should not be placed too deeply as septum identification and needle placement will be difficult.

Prior to closure of the wounds, system patency should be verified.

4.1. IMPLANTATION OF SOPH-A-PORT VENOUS

Equipment required:

- 1 subclavian introducer set (if percutaneous procedure)
- 1 Soph-A-Port
- 1 20 Gauge blunt tip needle (if unattached model)
- 1 22 Gauge Huber non-coring needle or infusion set
- 1 10 ml syringe filled with heparinized normal saline (50 to 500 U/ml).

4.1.1. Patient installation

Patient should be placed in Trendelenburg position, head turned in the opposite direction to implanted side.

Aseptically prepare, drape, and anesthetize operative site.

4.1.2. Port and Catheter Preparation

In case of unattached model, port and catheter are flushed separately:

1. Connect 20 G blunt tip needle to syringe containing heparinized saline. Attach to distal tip of catheter and dispel air. Clamp catheter proximal tip with atraumatic or vascular clamp.
2. Attach the Huber non-coring needle to syringe of heparinized saline. Penetrate port body septum. Gently flush air from port body while holding port such that outlet connector is straight up.

In case of preattached model, port and catheter are flushed simultaneously directly through the body septum with the syringe of heparinized saline attached to the non-coring needle. Clamp catheter tip.

Intravenous placement of the catheter can be performed by percutaneous puncture or by standard surgical procedure.

In case of preattached model, the length of the catheter must be adjusted precisely prior to insertion.

4.1.3. Intravenous catheter placement

4.1.3.1. Percutaneous procedure

1. In case of preattached model, a dissection pocket and a subcutaneous tunnel to entry site is recommended prior to catheter insertion.
2. Choose appropriate size introducer. Puncture the Subclavian vein using 18 G needle with syringe. Gently aspirate as needle is inserted. Remove syringe leaving 18 G needle in place.
3. Insert flexible metallic guide wire through needle. Use appropriate X-ray technique to verify correct position of guide wire. Remove needle.

4. Vessel dilator with sheath are advanced over exposed flexible metallic wire as far as possible. Withdraw vessel dilator and flexible guide wire, leaving sheath in place.

CAUTION: Prevent air embolism by placing thumb over exposed orifice of sheath.

5. Insert catheter into sheath. Firmly advance distal tip of the catheter through sheath until properly positioned. Peel away sheath while withdrawing it from vessel. Peel both tabs of sheath. Care should be taken not to withdraw catheter as sheath is removed.

Catheter position in superior Vena Cava should be confirmed by X-ray. Secure catheter in place.

CAUTION:

During percutaneous insertion into the subclavian vein, place the catheter away from the angle formed by the clavicle and the first rib, in order to avoid crushing and rupture of the catheter ("Pinch-off").

4.1.3.2. Standard Surgical Procedure

1. Incision is made in supra clavicular fossa to expose cephalic vein or a small transverse incision in neck to expose external jugular vein. Isolate vessel.
2. Perform a venotomy. Introduce and advance catheter previously heparinized and clamped to desired location. Confirm catheter placement by appropriate X-ray technique. Suture of the catheter to the vein is not recommended.
3. Close the venotomy.

CAUTION:

When using external jugular vein, carefully position the catheter over clavicle to avoid kinking or compression.

4.1.4. Dissection of Port Pocket

Port body should be positioned on anterior chest wall, approximately 1 cm under the surface of the skin.

A pocket is formed either from a small thoracic incision or by extension from the cervical incision, by blunt dissection down to underlying fascia.

The pocket is appropriately prepared with regard to the size of the port so that looseness is prevented.


Perfect hemostasis must be obtained prior to port fixation.

A subcutaneous tunnel is performed between the pocket and the venous entry site.

4.1.5. Connection Port-Catheter


1. Trim proximal end of catheter to desired length. Slide the screw-lock on to proximal end of catheter. Leave 1 or 2 cm of catheter protruding (Schema 2).

SCHEMA 2


- Advance catheter over the connector up to the abutment (Schema 3).

SCHEMA 3


- Screw the lock to port body (Schema 4).

SCHEMA 4

CAUTION: The connector is mobile and kept on its seat by the screw-lock; traction exerted on the catheter, before screwing the lock or when disconnection is needed, will pull out the connector from its seat (Schema 5). In that case replace the connector correctly on its seat. Never remove the silicone joint from the connector.

SCHEMA 5

After screwing the lock, exert traction on the catheter to check the correct catheter connection. If the lock is not totally screwed, disconnection and/or leakage can occur.

As the size of the connector and the screw-lock are specific to each catheter, Sophysa does not warranty the performance of the catheter connection to Soph-A-Port when using different catheter than provided with each port by Sophysa.

There is high risk of disconnection and/or leakage when using different catheter than provided with the port by Sophysa.

- Insert the port body in the pocket. It should be secured to underlying fascia with non-absorbable sutures by using integrated sutures holes. Avoid percing catheter with the suture needle. Care should

be exercised so that the septum of port is placed at distance of the incision after closure, and so that there is no kinking of the catheter.

- Prior to wound closure, verify patency of the system. Puncture the septum by using a 22 G Huber non-coring needles and a 10 ml syringe filled with 5 ml of heparinized saline. Aspirate to confirm ability to withdraw blood and flush port with the 5 ml heparinized saline (50 to 500 U/ml). Maintain positive pressure on syringe plunger to avoid reflux of blood into catheter tip while withdrawing needle.
- Close incision by appropriate surgical technique. Dress wound.

4.2. IMPLANTATION OF SOPH-A-PORT ARTERIAL

Prior to implantation, prepare port and catheter as described in Soph-A-Port intravenous implantation procedure.

The gastroduodenal artery is the most common location for intra-arterial catheter placement, with the port body placed at bottom of the chest wall.

Other sites may be used depending on the organ targeted for drug delivery.

The catheter insertion site is determined by direct visualization due to variation in patient anatomy.

4.2.0.1. Standard Surgical Procedure

- Clamp gastroduodenal artery upstream and ligature the artery downstream of the chosen entry site.
- Trim catheter to necessary length ensuring that at least two suture rings are left on catheter.
- Insert catheter into artery through small arteriotomy. Catheter is advanced up to, but not into, lumen of common hepatic artery.
- Verify correct position of the catheter and ligate both sides of fixation rings.

CAUTION:

Use care when placing suture around artery to avoid cutting or occluding the catheter.

- Remove the clamp.

4.2.0.2. Percutaneous Procedure

Percutaneous puncture of the femoral artery by Sedlinger method is possible. Advancing of the catheter tip up to the gastroduodenal artery will be followed under X-ray control.

Polyurethane catheter without suture rings are used.

For port body placement and catheter connection, follow guidelines as described under Soph-A-Port intravenous implantation procedure.

4.3. IMPLANTATION OF SOPH-A-PORT SPINAL

4.3.1. Patient installation

Patient is placed in lateral decubitus position.

Aseptically prepare, drape and anesthetize operative site.

4.3.2. Catheter and port placement

1. Make a small vertical incision between two spines at the chosen level.
2. Insert the Tuohy needle into the center of incision, push the needle into the peridural space using appropriate technique.
3. Place through the needle the catheter with its guide up to the desired level under X-ray control.
4. Pull out the guide of the catheter.
5. Remove the Tuohy needle while keeping the catheter in place.

Port body is usually positioned in the lower thoracic antero-lateral region.

The spinal catheter is tunneled over the thoracic wall up to the pocket prepared for the port.

For port body placement and catheter connection, follow guidelines as described under Soph-A-Port intravenous implantation procedure.

1. Prior to wound closure verify the patency of the system; puncture the septum by using a 22 G Huber non-coring needle attached to a 10 ml syringe filled with physiological saline and flush the system.
2. Close incisions by appropriate surgical technique. Dress wound.

4.4. IMPLANTATION OF SOPH-A-PORT PERITONEAL

4.4.1. Patient installation

Patient is placed in decubitus position. Aseptically prepare, drape and anesthetize operative site.

4.4.2. Catheter and port placement

1. Make a small midline vertical incision close to umbilicus. Grasp peritoneum with forceps; elevate and carefully incise under direct vision.
2. Insert index finger into peritoneal incision to check for adhesions.
3. Insert catheter into selected area. Place purse-string suture around catheter to prevent leakage though not too tight so as to occlude catheter.
4. Confirm catheter positioning by X-ray.

CAUTION:

To avoid injection into subcutaneous tissue, ensure that all catheter holes remain within peritoneal cavity.

The port is usually placed in the lower thoracic antero-lateral region.

The peritoneal catheter is tunneled from the umbilicus to the pocket prepared for the port.

Follow the guidelines as described under Soph-A-Port intravenous implantation procedure for port body placement and catheter connection.

1. Attach 22 G Huber non-coring needle to 50 ml syringe filled with physiological saline. Flush system with saline solution and immediately aspirate to confirm that there is neither obstruction nor leakage.
2. Close incisions by appropriate surgical technique. Dress wound.

5. Sterilization

Soph-A-Port is individually packaged in a double sterile packaging. It is sterilized with ethylene oxide.

Soph-A-Port must not be implanted if the packaging is opened or damaged or beyond the expiry date.

It is for single use only. It must not be resterilized and reused after explantation.

6. Precautions

To minimize the possibility of contamination, very careful aseptic technique should be observed during inspection and placement of the port.

When suturing around the catheter, avoid excessive suture tightness to prevent occlusion of the catheter.

Avoid piercing catheter with suture needle. Do not use toothed forceps for clamping.

For direct subclavian access, lateral placement of the catheter which avoids the costal clavicular angle is recommended to minimize the possibility of the catheter shear.

Soph-A-Port may be used immediately after verification of catheter tip placement.

7. Complications

The complications which can result from implantation of an access port include the risk inherent to the use of drugs, any surgical procedure and insertion of a foreign body.

The potential risks normally associated with any access port are applicable to Soph-A-Port.

The patient or the patient's family must be informed about these possible complications.

Complications related to the insertion technique of the system; they include, but are not limited to:

- Surgical Complications (Haemorrhage, Haematoma..)
- Air Embolism
- Pneumothorax, Haemothorax
- Arterial puncture
- Arterio-venous fistula
- Infection.

Complications related to the system; they include, but are not limited to:

- Device Rejection
- Device Rotation or Extrusion
- Catheter Occlusion
- Rupture of the catheter
- Catheter Disconnection or Migration
- Cardiac Arrhythmia
- Local Infection, Peritonitis, endocarditis, Septicemia
- Vessel Thrombosis, Thrombo-embolism
- Extravasation
- Air Embolism.

8. Use and maintenance of SOPH-A-PORT

8.1. GENERAL GUIDELINES

Access to Soph-A-Port should be performed using aseptic technique, by qualified personnel on the order of a physician.

Do not use a syringe size smaller than 10 ml which can induce overpressure in the system.

Only Huber non-coring needles 22 G or smaller diameter, straight or 90° angled should be used with Soph-A-Port.

Using bigger diameter than 22 G can induce decreasing of the septum performance.

As air should never penetrate in the system, stopcock should be used to minimize risk of air embolism.

The non-coring needle should be inserted perpendicular to the implanted port. The needle should be advanced through the septum until it contacts the base of the port.

Once correctly positioned in the septum, do not move the needle.

Verify compatibility of the drug with heparinized saline prior to injection.

Before injection in Soph-A-Port the physician should observe contraindications, cautions, warnings and instructions as notified on the drug labelling.

8.1.0.1. Bolus Injection - Infusion

Locate Soph-A-Port by palpating the port body.

Observing aseptic technique, using surgical gloves and sterile accessories, prepare injection site.

Attach syringe with normal saline with stopcock and non-coring needle. Flush to remove air.

Aseptically immobilize the port between two fingers.

Insert needle into the center of Soph-A-Port septum at perpendicular angle. Push the needle until it contacts the base of port body.

Note: In case of vascular access, verify patency of the system by aspiration of 3-5 ml of blood.

Inject 20 ml of normal saline to flush port catheter. Put stopcock in closed position, remove empty syringe and attach drug syringe. Put stopcock in open position and inject drug slowly.

Note: In case of multiple drug injections, rinse with 5 ml saline between each injection.

For continuous infusion, it is advisable to use a right angle anti-coring needle with extension tubing or a winged infusion set.

Connect infusion pump to extension tubing. Start infusion pump, then open tubing clamp.

Always apply occlusive dressing (preferably a transparent one) over needle and port site.

Note: Once the infusion finished, clamp tubing, then switch off the pump.

CAUTION:

Examine injection site closely.

Discontinue injection if it does not flow correctly. If patient feels an abnormal sensation or pain at injection site, it may indicate that drug has extravasated. In that case discontinue infusion immediately and notify physician.

8.1.1. At the end of drug injection

Clamp tubing and/or put stopcock in closed position and disconnect drug syringe.

Re-attach syringe filled with 10 ml saline; unclamp tubing and/or put stopcock in open position and rinse site and catheter.

Attach syringe containing 5 ml of heparinized saline (50 to 500 U/ml).

Flush the catheter by injection at a flow rate of 5 ml/mn maximum, creating an "heparin-lock".

Slowly remove the needle while injecting for keeping a positive pressure in the site to prevent blood reflux.

Apply adhesive bandage.

8.1.1.1. Blood Sampling

Insert non-coring needle with stopcock into prepared port site and flush catheter with 5 ml of normal saline.

Note: Blood sampling may be done by syringe aspiration or by use of a vacutainer.

Withdraw and discard the first 5-10 ml blood. Perform required blood sampling. Immediately flush catheter with 20 ml of normal saline.

Re-establish an "heparin lock".

8.1.1.2. Intrapitoneal Infusion Procedure

Due to flow rate limitations, Soph-A-Port is not recommended for continuous ambulatory peritoneal dialysis.

Prepare the solution bottle and its tubing set.

Locate Soph-A-Port by palpating the port body.

Observing aseptic technique, using surgical gloves, antiseptic solutions and sterile accessories, prepare injection site.

Attach 10 ml syringe filled with normal saline to right angle non-coring needle and extension tubing or winged infusion set. Flush tubing to remove air. Clamp tubing.

Aseptically immobilize the port between two fingers.

Insert needle into the center of Soph-A-Port septum at a perpendicular angle until it contacts the base of body port.

Inject contents of syringe filled with normal saline to flush port catheter. Close clamp. Remove syringe.

Connect tubing set to needle. Open clamps to let the fluid flow.

CAUTION:

Examine injection site closely.

Discontinue injection if it does not flow correctly.

If patient feels an abnormal sensation or pain at delivery site, it may indicate that drug has extra-vasated. In that case discontinue delivery immediately and notify physician.

Apply occlusive dressing (preferably a transparent one) over needle and port site.

8.1.2. At the end of fluid instillation

Clamp tubing set of solution bottle and extension tubing of the needle.

Remove dressing.

Disconnect tubing set from needle.

Attach syringe filled with 10 ml saline to the needle and rinse site and catheter.

Remove needle from port septum, and apply adhesive bandage.

8.2. MAINTENANCE

Vascular access

A "heparin lock" should be administered after each injection, infusion and blood sampling.

In addition a "heparin lock" should also be administered:

- Once a month for venous access
- Once a week for arterial access
- Before each session of chemotherapy.

8.3. TROUBLESHOOTING

Injection out of the chamber with extravasation inducing pain and local swollen:

1. Stop the injection.
2. Re-aspirate locally the product injected while removing the needle.
3. Dress the site.
4. Notify immediately to the physician.

Able to infuse through Soph-A-Port but unable to aspirate:

The catheter tip may be abutting the vessel wall.

Aspiration causes the vessel wall to be drawn into the catheter thus blocking blood withdrawal. An infusion, however, pushes the wall away from the tip and restores patency.

Slowly inject heparinized saline, check there is no extravasation; the site can then be used for treatment.

However some attempts may be helpful to restore the ability to aspirate from Soph-A-Port such as:

- placing the patient in Trendelenburg position
- turning the head in the opposite direction from the implanted site
- placing the arm over head.

Unable neither to aspirate nor to inject through Soph-A-Port:

Verify correct positioning of the needle in the chamber; check the needle is long enough and check the needle tip can contact the base of the chamber.

Incompatible drugs may have been infused creating a precipitate which can obstruct the system.

Undertake a delicate mobilisation of the port through the skin to release a possible kinking at the connexion port-catheter.

Change head and/or arm position of the patient to check kinking of the catheter in the costo-clavicular angle.

In case of obstruction of the port or catheter, notify to the physician.

Obstruction of the system may also be caused by a thrombosis encapsulating the catheter tip. A contrast study done through the port will confirm the thrombosis. A fibrinolytic agent such as urokinase may be used on the order of a physician to restore patency.

The procedure should be performed as outlined by the drug manufacturer's labelling.

CAUTION:

Do not force fluids through Soph-A-Port to clear an obstruction. A high pressure situation may cause irreversible catheter damage. Do not use a syringe smaller than 10 ml.

If there is any doubt regarding patency of the catheter, the system must be removed.

9. Patient identification card

A Patient Identification Card and a Patient Information Booklet is provided with each Soph-A-Port. The card completed by the physician should be given to the patient who should be advised to carry it all times.

This card contains pertinent information about Soph-A-Port.

10. Return to Sophysa

Ports and attached catheters which are explanted due to suspected malfunction should be returned to Sophysa sunk in sterile water, for analysis.

The clinician must contact Sophysa before shipment.

11. Warranty

Sophysa hereby warrants solely to original implanting purchaser of a Soph-A-Port that it will provide, at no charge, a replacement for the system that fails for any reasons attributable to defects in materials and workmanship, for a period of 90 days from the date of implant provided the following conditions are satisfied:

- The use of the system must be in accordance with Instructions Manual which is included with the product
- The explanted system must be returned to Sophysa within 30 days of explant, and becomes the property of Sophysa.

Occlusion of the catheter resulting in the explantation of the system is not considered failure of the device, and therefore is not covered under terms of this Limited Warranty.

This limited warranty applies only to the use of the system for the initial patient and the sole liability of Sophysa is limited to the replacement of the system under the conditions set forth herein.

No representation or warranty is made that failure of the system will not occur, or that the body will not react adversely to the implantation of the system or that the medication will be safe and effective for the treatment of

the user's condition; Sophysa disclaim the responsibility for any medical complications, including death, resulting from the use of the system.

This warranty is expressly in lieu of all other warranties. There shall be no liability on the part of Sophysa for any general, special, incidental, or consequential damages arising out of the sale, manufacture or use of the product sold hereunder.

Sophysa makes no warranties expressed or implied, including, but not limited to, any warranty of merchantability or fitness of the products for any purpose or reason, with respect to the Soph-A-Port implantable drug delivery system sold under this warranty; except as in this statement otherwise expressly provided.

12. Symbols

REF

Catalog reference

Manufacturer

Date of manufacture

Caution, see the Instructions for Use

STERILE EO

Sterilized using Ethylene Oxide

Do not re-use

Do not re-sterilize

Use by

LOT

Batch code

SN

Serial number

CE

CE Conformity Marking

13. References

VENOUS

- SP 20110** **Soph-A-Port Standard venous unattached**
Port. Base: 33x27 mm. Height: 13 mm. Septum diameter: 13 mm. Septum thickness: 6.5 mm. Residual volume: 0.45 ml. Weight: 8.2 g.
Catheter. ID: 1.1 mm. OD: 2.2 mm. Internal volume: 0.113 ml / 100 mm. Length: 700 mm. Material: Silicone XRO. Graduations at 5, 10, 15, 20 cm.
- SP 20111** **Soph-A-Port Standard venous preattached**
Port. Base: 33x27 mm. Height: 13 mm. Septum diameter: 13 mm. Septum thickness: 6.5 mm. Residual volume: 0.45 ml. Weight: 8.2 g.
Catheter. ID: 1.1 mm. OD: 2.2 mm. Internal volume: 0.113 ml / 100 mm. Length: 700 mm. Material: Silicone XRO. Graduations at 5, 10, 15, 20 cm.
- SP 10110** **Soph-A-Port Mini venous unattached**
Port. Base: 26x21 mm. Height: 10 mm. Septum diameter: 10 mm. Septum thickness: 4.8 mm. Residual volume: 0.24 ml. Weight: 3.8 g.
Catheter. ID: 1.1 mm. OD: 2.2 mm. Internal volume: 0.113 ml / 100 mm. Length: 700 mm. Material: Silicone XRO. Graduations at 5, 10, 15, 20 cm.
- SP 10111** **Soph-A-Port Mini venous preattached**
Port. Base: 26x21 mm. Height: 10 mm. Septum diameter: 10 mm. Septum thickness: 4.8 mm. Residual volume: 0.24 ml. Weight: 3.8 g.
Catheter. ID: 1.1 mm. OD: 2.2 mm. Internal volume: 0.113 ml / 100 mm. Length: 700 mm. Material: Silicone XRO. Graduations at 5, 10, 15, 20 cm.
- SP 20120** **Soph-A-Port Standard venous high flow unattached**
Port. Base: 33x27 mm. Height: 13 mm. Septum diameter: 13 mm. Septum thickness: 6.5 mm. Residual volume: 0.45 ml. Weight: 8.2 g.
Catheter. ID: 1.5 mm. OD: 3.0 mm. Internal volume: 0.201 ml / 100 mm. Length: 700 mm. Material: Silicone XRO. Graduations at 5, 10, 15, 20 cm.
- SP 20121** **Soph-A-Port Standard venous high flow preattached**
Port. Base: 33x27 mm. Height: 13 mm. Septum diameter: 13 mm. Septum thickness: 6.5 mm. Residual volume: 0.45 ml. Weight: 8.2 g.
Catheter. ID: 1.5 mm. OD: 3.0 mm. Internal volume: 0.201 ml / 100 mm. Length: 700 mm. Material: Silicone XRO. Graduations at 5, 10, 15, 20 cm.
- SP 10120** **Soph-A-Port Mini venous high flow unattached**
Port. Base: 26x21 mm. Height: 10 mm. Septum diameter: 10 mm. Septum thickness: 4.8 mm. Residual volume: 0.24 ml. Weight: 3.8 g.
Catheter. ID: 1.5 mm. OD: 3.0 mm. Internal volume: 0.201 ml / 100 mm. Length: 700 mm. Material: Silicone XRO. Graduations at 5, 10, 15, 20 cm.
- SP 10121** **Soph-A-Port Mini venous high flow preattached**
Port. Base: 26x21 mm. Height: 10 mm. Septum diameter: 10 mm. Septum thickness: 4.8 mm. Residual volume: 0.24 ml. Weight: 3.8 g.
Catheter. ID: 1.5 mm. OD: 3.0 mm. Internal volume: 0.201 ml / 100 mm. Length: 700 mm. Material: Silicone XRO. Graduations at 5, 10, 15, 20 cm.
- SP 20130** **Soph-A-Port Standard venous small diameter unattached**
Port. Base: 33x27 mm. Height: 13 mm. Septum diameter: 13 mm. Septum thickness: 6.5 mm. Residual volume: 0.45 ml. Weight: 8.2 g.
Catheter. ID: 1.05 mm. OD: 1.65 mm. Internal volume: 0.095 ml / 100 mm. Length: 700 mm. Material: Polyurethane XRO. Graduations every 5 cm from 5 to 50 cm.
- SP 20131** **Soph-A-Port Standard venous small diameter preattached**
Port. Base: 33x27 mm. Height: 13 mm. Septum diameter: 13 mm. Septum thickness: 6.5 mm. Residual volume: 0.45 ml. Weight: 8.2 g.
Catheter. ID: 1.05 mm. OD: 1.65 mm. Internal volume: 0.095 ml / 100 mm. Length: 700 mm. Material: Polyurethane XRO. Graduations every 5 cm from 5 to 50 cm.
- SP 10130** **Soph-A-Port Mini venous small diameter unattached**
Port. Base: 26x21 mm. Height: 10 mm. Septum diameter: 10 mm. Septum thickness: 4.8 mm. Residual volume: 0.24 ml. Weight: 3.8 g.

- Catheter.** ID: 1.05 mm. OD: 1.65 mm. Internal volume: 0.095 ml / 100 mm. Length: 700 mm. Material: Polyurethane XRO. Graduations every 5 cm from 5 to 50 cm.
- SP 10131** **Soph-A-Port Mini venous small diameter preattached**
Port. Base: 26x21 mm. Height: 10 mm. Septum diameter: 10 mm. Septum thickness: 4.8 mm. Residual volume: 0.24 ml. Weight: 3.8 g.
Catheter. ID: 1.05 mm. OD: 1.65 mm. Internal volume: 0.095 ml / 100 mm. Length: 700 mm. Material: Polyurethane XRO. Graduations every 5 cm from 5 to 50 cm.

SPINAL

- SP 20210** **Soph-A-Port Standard spinal unattached**
Port. Base: 33x27 mm. Height: 13 mm. Septum diameter: 13 mm. Septum thickness: 6.5 mm. Residual volume: 0.45 ml. Weight: 8.2 g.
Catheter. ID: 0.76 mm. OD: 1.65 mm. Internal volume: 0.058 ml / 100 mm. Length: 800 mm. Material: Silicone XRO. Graduations at 11, 16, 21, 26 cm.
- SP 10210** **Soph-A-Port Mini spinal unattached**
Port. Base: 26x21 mm. Height: 10 mm. Septum diameter: 10 mm. Septum thickness: 4.8 mm. Residual volume: 0.24 ml. Weight: 3.8 g.
Catheter. ID: 0.76 mm. OD: 1.65 mm. Internal volume: 0.058 ml / 100 mm. Length: 800 mm. Material: Silicone XRO. Graduations at 11, 16, 21, 26 cm.

ARTERIAL

- SP 20310** **Soph-A-Port Standard arterial unattached**
Port. Base: 33x27 mm. Height: 13 mm. Septum diameter: 13 mm. Septum thickness: 6.5 mm. Residual volume: 0.45 ml. Weight: 8.2 g.
Catheter. ID: 0.6 mm. OD: 2.2 mm. Internal volume: 0.038 ml / 100 mm. Length: 500 mm. Material: Silicone XRO. Retention rings at 5, 6, 8 cm.
- SP 10310** **Soph-A-Port Mini arterial unattached**
Port. Base: 26x21 mm. Height: 10 mm. Septum diameter: 10 mm. Septum thickness: 4.8 mm. Residual volume: 0.24 ml. Weight: 3.8 g.
Catheter. ID: 0.6 mm. OD: 2.22 mm. Internal volume: 0.038 ml / 100 mm. Length: 500 mm. Material: Silicone XRO. Retention rings at 5, 6, 8 cm.
- SP 20330** **Soph-A-Port Standard arterial small diameter unattached**
Port. Base: 33x27 mm. Height: 13 mm. Septum diameter: 13 mm. Septum thickness: 6.5 mm. Residual volume: 0.45 ml. Weight: 8.2 g.
Catheter. ID: 1.05 mm. OD: 1.65 mm. Internal volume: 0.095 ml / 100 mm. Length: 900 mm. Material: Polyurethane XRO. Graduations every 5 cm from 5 to 70 cm.
- SP 10330** **Soph-A-Port Mini arterial small diameter unattached**
Port. Base: 26x21 mm. Height: 10 mm. Septum diameter: 10 mm. Septum thickness: 4.8 mm. Residual volume: 0.24 ml. Weight: 3.8 g.
Catheter. ID: 1.05 mm. OD: 1.65 mm. Internal volume: 0.095 ml / 100 mm. Length: 900 mm. Material: Polyurethane XRO. Graduations every 5 cm from 5 to 70 cm.

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