Instructions for use

SOPH-A-PORT® IMPLANTABLE ACCESS PORT

Sterile, for single use only

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Page intentionally left blank


Instructions for use

SOPH-A-PORT[®] IMPLANTABLE ACCESS PORT
Sterile, single use

Table of content

| | | |
|-------|--|----|
| 1. | Description | 3 |
| 2. | Indications | 3 |
| 3. | Contraindications | 3 |
| 4. | SOPH-A-PORT implantation instructions | 3 |
| 4.1. | Identification of the implantation site | 4 |
| 4.2. | Venous Implantation of SOPH-A-PORT | 4 |
| 4.3. | Arterial implantation of SOPH-A-PORT | 6 |
| 4.4. | Spinal implantation of SOPH-A-PORT | 6 |
| 4.5. | Peritoneal implantation of SOPH-A-PORT | 6 |
| 5. | Sterilisation | 7 |
| 6. | Precautions | 7 |
| 7. | Complications | 7 |
| 8. | Use and maintenance of SOPH-A-PORT | 7 |
| 8.1. | General guidelines | 7 |
| 8.2. | Maintenance | 9 |
| 8.3. | Troubleshooting | 9 |
| 9. | Replacement of the SOPH-A-PORT port and catheter | 10 |
| 9.1. | Replacement of the port | 10 |
| 9.2. | Replacement of the catheter | 10 |
| 10. | Patient identification card | 10 |
| 11. | Treatment of the product after use | 10 |
| 12. | Return to Sophysa | 11 |
| 13. | Warranty | 11 |
| 14. | Symbols | 11 |
| 15. | References | 12 |
| 15.1. | VENOUS | 12 |
| 15.2. | SPINAL | 12 |
| 15.3. | ARTERIAL | 12 |
| 15.4. | PERITONEAL | 12 |

SCHEMA 1

1- Body; 2- Septum; 3- Silicone joint; 4- Abutment; 5- Thread; 6- Catheter; a- Port; b- Mobile connector; c- Screw-lock

1. Description

SOPH-A-PORT® is a totally implantable device designed to allow repeated access to the vascular system or selected area of the body by preventing the trauma associated with multiple punctures or the inconvenience of an external catheter.

Injection into the SOPH-A-PORT is performed by percutaneous puncture using a special needle.

SOPH-A-PORT consists of a polysulfone injection port and a self-sealing silicone septum. Suture holes are integrated into the port base to facilitate fixation to the underlying fascia. The port catheter is in radiopaque silicone.

A patented self-locking system with a titanium mobile connector and a screw-lock, ensures safe, efficient connection of the catheter to the port (Diagram 1).

A serial number is engraved on each port.

SOPHYSA manufactures various SOPH-A-PORT models.

These models have different port and catheter sizes (see table of product references).

Some SOPH-A-PORT are available with preattached or unattached catheter configurations. All SOPH-A-PORT are supplied with a Huber Point (non-coring) 22 G needle.

WARNING

SOPH-A-PORT is not designed for high pressure injection of imaging solutions. Should high pressure injection of such solutions be necessary, in a patient with a SOPH-A-PORT, another access should be used.

2. Indications

SOPH-A-PORT® is indicated when treatment of the patient requires repeated access to the vascular system or a selected area of the body for:

- Injection of a drug
- Infusion
- Blood sampling
- Transfusion
- Injection of imaging solutions
- Parenteral nutrition.

3. Contraindications

SOPH-A-PORT is not designed, sold or offered for use in any indications other than those mentioned above.

The following contraindications should be observed:

SOPH-A-PORT should not be prescribed for patients for whom the planned medical treatment requires the use of substances known to be incompatible with the manufacturing materials.

It is recommended that the compatibility of the drugs with the manufacturing materials is checked prior to injection.

SOPH-A-PORT should not be implanted in a patient whose body is too small to support the dimensions, or in a patient suspected of being allergic to the material of the components.

Implantation of a SOPH-A-PORT should be avoided in the presence of a known or suspected local or general infection, unless access to the vascular system with SOPH-A-PORT is the solution for treating the infection.

4. SOPH-A-PORT implantation instructions

WARNING

The doctors responsible for the implantation must have suitable experience of port and catheter implantation procedures, as well as access procedures, and must read all the instructions for use and labels carefully before using the device.

Implantation of the SOPH-A-PORT and catheter placement may be performed using various surgical techniques; the physician must choose the appropriate method according to his experience and patient's needs.

As for any implantation of a foreign body, the technique must be strictly aseptic.

WARNING

The site and implantation technique must be carefully chosen to avoid any site of anatomical interest or to prevent the propagation of tumours.

The port catheter should be positioned at the selected treatment site and fixed using an appropriate surgical technique to prevent catheter migration.

Its position should be checked by an appropriate X-ray.

Sufficient slack should be left between the catheter insertion point and the port body to prevent any strain on the catheter.

The port body should be secured to the underlying fascia in a subcutaneous pocket created by blunt or sharp dissection. Fixing should be carried out using the suture holes incorporated into the port in order to prevent any distortion and traction on the catheter.

Care should be taken so that the incision does not cross the septum of the port after closure.

When selecting the port implantation position, it is recommended that locations which interfere with patient mobility or clothing by creating a potential pressure point, and which would also interfere with the treatment, be avoided.

The implantation site should be chosen taking the pathology into account.

To facilitate access, the port body should be placed over a bony surface (rib cage). The port should not be placed too deeply as this would make septum identification and needle placement difficult.

The port should be placed about 1 cm under the skin so that it can be palpated for subsequent accesses.

PRECAUTION

Do not expose the SOPH-A-PORT catheter or body to any substance containing povidone iodine.

Before closing the wound, check that the system is completely patent.

NOTES

Consider routine administration of peri-operative antibiotic therapy, and not via the catheter.

The antibiotic therapy regimen established by the hospital for CSF shunts may be used as a guide.

WARNING

As the size of the mobile connector and the screw-lock are specific to the catheter, only use the catheter supplied by Sophysa for connection to the SOPH-A-PORT, in order to avoid any dislodgement or leakage.

4.1. IDENTIFICATION OF THE IMPLANTATION SITE

WARNING

Before proceeding with the implantation of the device, the patient must be assessed for the presence of other existing devices, from which the new device must be clearly differentiated for example by its location, external marking or other means of identification.

To determine a suitable implantation site, the following factors should be taken into account:

- Normal and abnormal variations, particularly the anatomical markers and abnormalities (particularly congenital abnormalities)
- Patient's level of activity
- Patient mobility

- Maturity of the patient, especially the maturity or immaturity of the organ systems including the immune system
- Growth rate, weight, height, body mass and body surface of the patient
- Development of the patient, particularly hormonal influences and the stage of puberty (e.g. the development of the mammary buds at the pre-adolescent and adolescent stages may affect the location for the device)
- Expected behaviour of the patient
- Home environment and level of assistance provided by the family or the care personnel
- Pressure points due to mobility, immobility or the patient's clothes being too tight (in places where the device is likely to produce an irritation or skin erosion)
- Location of other implanted devices

4.2. VENOUS IMPLANTATION OF SOPH-A-PORT

Equipment required:

- 1 subclavian introducer system (if percutaneous procedure)
- 1 SOPH-A-PORT
- 1 x 20 Gauge blunt tip needle
- 1 x 22 Gauge Huber Point (non-coring) needle or infusion set
- 1 x 10 ml syringe filled with heparinized saline (50 to 500 U/ml).

4.2.1. Installation of the patient

Place the patient in the Trendelenburg position, head turned in the opposite direction to the implanted side.

Disinfect, place the surgical drapes and anaesthetise the area concerned.

4.2.2. Port and catheter preparation

Rinse the port and catheter separately.

1. Fit the 20 G blunt tip needle to the syringe containing heparinized saline. Place it at the distal end of catheter and inject to purge the air. Clamp the proximal end of the catheter with an atraumatic or vascular clamp.
2. Attach the Huber Point (non-coring) needle to the syringe of heparinized saline. Inject into the port septum. Carefully evacuate the air from the port while holding this so that the outlet connector is vertical.

Intravenous placement of the catheter can be performed by percutaneous puncture or by standard surgical procedure.

4.2.3. Intravenous catheter placement

4.2.4. Percutaneous procedure

1. Choose an appropriate size introducer. Locate the subclavian vein using an imaging system (such as ultrasound) and inject with an 18 G needle with syringe. Carefully aspirate as the needle is inserted. Remove the syringe leaving the 18 G needle in place.
2. Insert the flexible metal guide into the needle. Check

that this guide is correctly positioned using an X-ray. Remove the needle.

3. Advance a vessel dilator with its sheath as far as possible on the exposed metal guide. Withdraw the vessel dilator and the wire, leaving the sheath in place.

CAUTION

To prevent an air embolism place the thumb over the exposed orifice of the sheath.

4. Insert the catheter into the sheath. Firmly advance the distal end of the catheter into the sheath until it is correctly positioned. Peel the sheath while withdrawing it from the vessel. Peel both tabs of the sheath. Care should be taken not to withdraw the catheter as the sheath is removed. Check the catheter position in the superior vena cava by X-ray. Fix the catheter in place.

CAUTION

During the percutaneous introduction into the subclavian vein, place the catheter away from the angle formed by the clavicle and the first rib, in order to avoid crushing and rupturing the catheter ("Shearing").

4.2.4.1. Standard surgical procedure

1. Make an incision in the supraclavicular fossa to expose a cephalic vein or a small transverse incision in the neck to expose the external jugular vein. Isolate the vein.
2. Perform a phlebotomy. Introduce and advance the previously heparinized catheter to the desired location and clamp it. Check its placement by an X-ray. Suturing the catheter to the vein is not recommended.
3. Close the phlebotomy.

CAUTION

When operating on the external jugular vein, carefully position the catheter above the clavicle to avoid kinking or compressing it.

4.2.5. Dissection of a pocket for the port

The body of the port should be positioned on the anterior chest wall, approximately 1 cm under the surface of the skin.

Form a pocket either by a small thoracic incision or by extension of the cervical incision, by blunt dissection down to underlying fascia.

Prepare the pocket in a suitable manner, depending on the size of the port, to prevent any mobility.


Haemostasis must be correctly established before fixing the port.

Make a subcutaneous tunnel between the pocket and the introduction site in the vein.

4.2.6. Port-catheter connection

1. Trim the proximal end of catheter to the desired length. Slide the screw-lock along the catheter so that 1 or 2 cm of catheter are left protruding from the proximal end (Diagram 2).


DIAGRAM 2

1- Thread; 2- Catheter

2. Advance the catheter on the connector up to the abutment (Diagram 3).


DIAGRAM 3

1- Abutment; 2- Catheter


3. Screw the lock on the port body (Diagram 4).

DIAGRAM 4

CAUTION: The connector is mobile and kept on its seat by the screw-lock system; any traction exerted on the catheter before screwing the locking system or when disconnection is needed will pull the connector out of its joint (Diagram 5). In this case replace the connector correctly on its seat. Never remove the silicone seal from the connector.

DIAGRAM 5

After screwing the lock, exert traction on the catheter to check it is properly connected. If the screw-lock is not completely tightened, disconnection and/or leakage may occur.

As the size of the connector and the sleeve are specific to each type of catheter, Sophysa does not guarantee the quality of the catheter connection to the SOPH-A-PORT if a catheter other than that supplied with each port by Sophysa is used.

There is high risk of disconnection and/or leakage if a catheter other than that supplied with the port is used.

4. Insert the port into the pocket. It should be fixed to the underlying fascia with non-absorbable sutures

- using the integrated sutures holes. Avoid piercing the catheter with the suture needle. Make sure that the septum is placed away from the incision after closure, and that there is no kinking of the catheter.
- Prior to wound closure, verify the patency of the system. Puncture the septum using a 22 G Huber Point (non-coring needle) and a 10 ml syringe filled with 5 ml of heparinized saline. Aspirate to check it is possible to sample blood and flush the port with the 5 ml heparinized saline (50 to 500 U/ml). Maintain positive pressure on the syringe plunger to avoid backflow of blood into end of the catheter while withdrawing the needle.
 - Close the incision by an appropriate surgical technique. Dress the wound.

4.3. ARTERIAL IMPLANTATION OF SOPH-A-PORT

Prior to implantation, prepare the port and catheter as described for SOPH-A-PORT venous implantation.

The gastroduodenal artery is the most common location for an intra-arterial catheter, with the port body being placed at the bottom of the chest wall.

Other sites may be used depending on the organ to which the drug should be administered.

The catheter insertion site should be determined by direct visual examination, given the variations in patient anatomy.

4.3.1. Standard surgical procedure

- Clamp the gastroduodenal artery upstream and ligature it downstream of the chosen introduction point.
- Cut the catheter to the necessary length ensuring that at least two fixing rings are left above this.
- Insert the catheter into the artery using a small arteriotomy. Advance it as far as, but not into, the lumen of the common hepatic artery.
- Verify the correct position of the catheter and ligate on either side of the fixing rings.

CAUTION

Place the suture carefully around the artery to avoid cutting or obstructing the catheter.

- Remove the clamp.

4.4. SPINAL IMPLANTATION OF SOPH-A-PORT

4.4.1. Installation of the patient

The patient is placed in the lateral decubitus position.

Disinfect, place the surgical drapes and anaesthetise the area concerned.

4.4.2. Catheter and port placement

- Make a small vertical incision between the two spinous processes at the chosen level.
- Insert the Tuohy needle into the centre of the incision and push it into the peridural space using an appropriate technique.
- Place the catheter through the needle using its metal guide up to the desired level under X-ray control.

- Pull out the catheter guide.
- Remove the Tuohy needle while holding the catheter in place.

PRECAUTIONS

If the catheter must be withdrawn during positioning, withdraw the Tuohy needle and the catheter as a single unit. Then carefully move the Tuohy needle back from the catheter and start the procedure again.

If there is any doubt about the integrity of the catheter use a new kit.

The port body is usually positioned in the lower antero-lateral thoracic region.

The spinal catheter is tunneled above the chest wall up to the pocket prepared for the port.

Place the port and make the catheter connection following the instructions given for the venous implantation of SOPH-A-PORT.

- Prior to wound closure verify the patency of the system; puncture the septum by using a 22 G Huber Point (non-coring) needle fitted on a 10 ml syringe filled with saline and flush the system.
- Close the incisions by an appropriate surgical technique. Dress the wound.

NOTE

When connecting the catheter to the port, do not introduce air bubbles into the catheter because this risks causing post-spinal anaesthesia headaches.

4.5. PERITONEAL IMPLANTATION OF SOPH-A-PORT

4.5.1. Installation of the patient

The patient is placed in the decubitus position. Disinfect, place the surgical drapes and anaesthetise the area concerned.

4.5.2. Catheter and port placement

- Make a small midline vertical incision close to the umbilicus. Pick up the peritoneum with the forceps; lift it and carefully make an incision under direct vision.
- Insert the index finger into the peritoneal incision to check there are no adhesions.
- Insert the catheter into the chosen area. Place the purse-string suture around the catheter to prevent leakage, but not pulling too tight so as not to obstruct the catheter.
- Confirm catheter positioning by X-ray.

CAUTION

To avoid injection into the subcutaneous tissue, ensure that all the catheter holes remain within the peritoneal cavity.

Plan a sufficient length of catheter (12 cm or more) to prevent the risk of catheter retraction.

The port is usually placed in the lower antero-lateral thoracic region.

The peritoneal catheter is tunneled from the umbilicus to the pocket prepared for the port.

Place the port and make the catheter connection following the instructions given for the venous implantation of SOPH-A-PORT.

1. Attach a 22 G Huber point (non-coring) needle to a 50 ml syringe filled with saline. Flush the system with this solution and immediately aspirate to check that there is no obstruction or leakage.
2. Close the incisions by an appropriate surgical technique. Dress the wound.

PRECAUTION

Prevent any pressure from the catheter on the intestines or organs in the pelvis and do not inject or aspirate from the port violently.

5. Sterilisation

Each SOPH-A-PORT is supplied in a sterile, apyrogenic single unit double package. It is sterilized with ethylene oxide.

The SOPH-A-PORT must not be implanted if the packaging has been opened or damaged or the expiry date has been exceeded.

It is for single use only. It must not be resterilised or reused after explantation.

6. Precautions

To minimize the possibility of contamination, the aseptic technique should be very meticulous during the control and placement of the port.

When suturing around the catheter, when applicable, avoid making the suture too tight so as not to obstruct the catheter.

Avoid piercing the catheter with the suture needle. Do not use crocodile clamps for clamping.

For direct subclavian access lateral placement of the catheter is recommended, thus avoiding the costoclavicular angle and the possibility of the catheter being sheared.

SOPH-A-PORT may be used immediately after verifying the positioning of the catheter end.

7. Complications

The possible complications from implantation of an access port are those from the inherent risk of the use of drugs, any surgical procedure and the introduction of a foreign body.

The potential risks normally associated with any access port are applicable to SOPH-A-PORT.

The doctor responsible for the implantation should inform the patient or the family, the care personnel or the patient's healthcare legal representative of the possible complications, some of which may cause serious lesions or even death.

Complications related to the system insertion technique they include, but are not limited to:

- Surgical complications (haemorrhage, haematoma, etc.)

- Air embolism
- Pneumothorax, haemothorax
- Arterial puncture
- Arteriovenous fistula
- Infection
- Seromas
- Meningitis
- Post-spinal anaesthesia headaches
- Thrombosis
- Pinching of the catheter
- Arterial dissection
- Extrahepatic perfusion
- Perforation of the intestine by the catheter

Complications related to the system; they include, but are not limited to:

- Local inflammation
- Device rejection
- Rotation or extrusion of the device
- Catheter obstruction
- Rupture of the catheter
- Disconnection or migration of the catheter
- Cardiac arrhythmia
- Local Infection, peritonitis, endocarditis, septicemia
- Embolisation and related symptoms:
 - Thrombosis, thrombo-embolism
 - Cerebral infarction
 - Extrahepatic perfusion
 - Liver failure
- Extravasation
- Air embolism
- Seromas
- Retraction of the catheter outside the peritoneal cavity
- Perforation of the intestine by the catheter
- Pain and discomfort during the treatment phase
- Complications related to the reactions of the vessels:
 - Angiospasm
 - Vasculitis
 - Aneurysm, false aneurysm
 - Venous stenosis

NOTES

If there is local inflammation, treat with antibiotic therapy and anti-inflammatory treatment.

If there is cardiac arrhythmia explant the port.

Depending on the type of complication the surgeon may have to consider explantation and implantation of a new port.

WARNING

If there is an arterial dissection stop the treatment and explant the port.

8. Use and maintenance of SOPH-A-PORT

8.1. GENERAL GUIDELINES

The use of the SOPH-A-PORT should be carried out using

an aseptic technique, by qualified personnel on the instructions of a physician.

Do not use a syringe of less than 10 ml likely to cause overpressure in the system.

Only 22 G or smaller diameter Huber Point (non-coring) needles, either straight or with a 90° angle should be used with SOPH-A-PORT.

Using a diameter larger than 22 G can cause a decrease in septum performance.

As air should never penetrate the port, a stopcock should be used to minimize risk of air embolism.

Adapt the length of the needle to the corpulence of the patient.

The needle should be introduced perpendicular to the implanted port. The needle should be pushed through the septum until it comes into contact with the base of the port.

Once it is correctly positioned in the septum, do not move the needle.

We recommend that the puncture points are varied in order to preserve the integrity of the skin and the impermeability of the septum.

Verify the compatibility of the drug with heparinized saline prior to injection.

Before injection into the SOPH-A-PORT the doctor should take into account the contraindications, precautions, warnings and instructions shown on the drug labelling.

WARNINGS

The device may not be compatible with organic solvents with very low polarity such as benzene, carbon tetrachloride, etc.

The device may not be compatible with docetaxel.

The compatibility of the device with a highly lipid solution has not been tested.

The device should not be used to administer drugs which are not compatible with silicon and sodium oligomers.

8.1.1. Bolus Injection - Infusion

Locate the SOPH-A-PORT by palpating to find the port body.

Prepare the injection area using an aseptic technique, wearing surgical gloves and using sterile accessories.

Fit the Huber Point (non-coring) needle and a stopcock onto a syringe containing saline. Purge the air.

Aseptically immobilize the port between two fingers.

Insert the needle into the centre of the SOPH-A-PORT septum perpendicularly. Push the needle until it comes into contact with the base of the port body.

Note: In the case of vascular access, aspirate 3 to 5 ml of blood to ensure the system is free.

Inject 20 ml of normal saline to flush the port catheter. Put the stopcock in the closed position, remove the empty syringe and fit the syringe containing the drug. Put the stopcock in the open position and inject the drug slowly.

Note: In the case of multiple drug injections, rinse with 5 ml saline between each injection.

For a continuous infusion, it is advisable to use a right

angled needle with extension tubing or a winged infusion set.

Connect the infusion pump to the extension tubing. Start the infusion pump, then open the tubing clamp.

Always apply an occlusive dressing (preferably transparent) on the needle and port.

Note: At the end of the infusion, clamp the tubing and then switch off the pump.

CAUTION

Meticulously examine the injection site.

Discontinue the injection if it the flow is not correct. If patient feels an abnormal sensation or pain at the injection site, it may indicate extravasation of the drug. In this case discontinue the infusion immediately and notify a doctor.

8.1.2. At the end of the drug injection

Clamp the tubing and/or put the stopcock in the closed position and then remove the drug syringe.

Fit a syringe filled with 10 ml saline; unclamp the tubing and/or put the stopcock in the open position and rinse the port and catheter.

Put the syringe containing 5 ml of heparinized saline (50 to 500 U/ml) in place.

Flush the catheter by injecting at a maximum flow rate of 5 ml/min, to create a "heparin-lock".

Slowly remove the needle while injecting in order to keep a positive pressure in the port to prevent a backflow of blood.

Apply an adhesive dressing.

8.1.2.1. Blood sampling

Insert a Huber Point (non-coring) needle with stopcock into the port and flush the catheter with 5 ml of normal saline.

Note: Blood sampling may be carried out by aspiration using a syringe or with a Vacutainer.

Take and discard the first 5-10 ml blood. Take the desired blood sample. Immediately flush the catheter with 20 ml of normal saline.

Re-establish a "heparin lock".

8.1.2.2. Intraperitoneal infusion method

As its flow rate is limited SOPH-A-PORT is not recommended for continuous ambulatory peritoneal dialysis.

Prepare the solution bottle and its tubing.

Locate the SOPH-A-PORT by palpating to find the port body.

Prepare the injection area using an aseptic technique, wearing surgical gloves, and using antiseptic solutions and sterile accessories.

Attach a 10 ml syringe filled with normal saline to a right-angled needle with extension tubing or a winged infusion set. Flush the tubing to remove the air. Clamp the tubing.

Aseptically immobilize the port between two fingers.

Insert the needle into the centre of the SOPH-A-PORT septum perpendicularly until it comes into contact with the base of the port.

Inject the saline contained in the syringe to flush the port catheter. Close the clamp. Remove the syringe.

Connect the infusion tubing to the needle. Open the clamp to allow the fluid to flow.

CAUTION

Examine the injection site carefully.

Discontinue the injection if the flow is not correct.

If patient feels an abnormal sensation or pain at the injection site, it may indicate extravasation of the drug. In this case discontinue the infusion immediately and notify a doctor.

Apply an occlusive dressing (preferably transparent) to the needle and port.

8.1.3. At the end of the fluid instillation

Clamp the solution bottle tubing and the needle extension tubing.

Remove the dressing.

Disconnect the tubing from the needle.

Fit the needle on the syringe filled with 10 ml saline and flush the port and catheter.

Remove the needle from the port septum, and apply an adhesive dressing.

8.2. MAINTENANCE

Vascular access

A "heparin lock" should be established after each injection, infusion and taking of blood samples, as well as:

- Once a month for venous access
- Once a week for arterial access
- Before each session of chemotherapy.

8.3. TROUBLESHOOTING

Injection outside the cavity with extravasation causing pain and local oedema:

1. Stop the injection.
2. Locally re-aspirate the product injected while removing the needle.
3. Apply a dressing.
4. Notify a doctor immediately.
5. Assess the need to stop the infusion, remove the device and put a local medical treatment in place.

Possible to infuse through the SOPH-A-PORT but impossible to aspirate:

The end of the catheter may be in contact with the vascular wall.

Aspiration causes the vessel wall to be drawn into the catheter thus blocking blood sampling. Conversely, the infusion pushes the wall away from the catheter end and clears the system.

Slowly inject the heparinized serum, check there is no extravasation; the site can then be used for treatment.

However some procedures may be useful to allow aspiration by the SOPH-A-PORT again:

- placing the patient in the Trendelenburg position

- head turned in the opposite direction to the implantation area
- arm above the head.

Impossible to either aspirate or inject via SOPH-A-PORT:

Check the needle is correctly positioned in the port, that the needle is long enough and that its tip can come into contact with the base of the port.

Incompatible drugs may have been injected and formed a precipitate which obstructs the system.

Mobilise the port gently through the skin to remove possible kinking at the port-catheter connection.

Change the head and/or arm position of the patient to check whether there is kinking of the catheter in the costoclavicular angle.

Take a new needle and repeat the injection procedure with standard, preservative-free saline.

If the port or catheter is obstructed, notify a doctor.

Obstruction of the system may also be caused by a thrombosis encapsulating the catheter end. Check for a thrombosis by an X-ray with contrast agent. A fibrinolytic agent such as urokinase may be used on the instruction of a physician to restore patency. Depending on the results, the surgeon may evaluate the need for explantation and replacement of the port.

Proceed as indicated in the instructions supplied by the drug manufacturer.

CAUTION

Do not try to force the fluid that must flow through the SOPH-A-PORT to clear an obstruction. A high pressure may cause irreversible catheter damage. The syringe should not be smaller than 10 ml.

If there is any doubt regarding the patency of the catheter, the system must be removed.

Warnings and precautions relating to the patient's daily life:

Inform the patient or the family, care personnel or the healthcare legal representative of the patient of the following warnings:

- Activities involving sudden, excessive or repeated bending, twisting, bouncing or stretching may damage or dislodge the components and require a further surgery or replacement. These activities may also bend or obstruct the catheter causing a loss of treatment, which may lead to a dosage error or an insufficient drug dose with clinically significant consequences.
- Patients must avoid rubbing the port or catheter excessively through the skin. A manipulation may cause the port to turn over so that the septum is no longer accessible. A manipulation may also cause a skin erosion, damage to the components, disconnection, bending or dislodging of the catheter, and may cause a leakage of the drug or CSF into the surrounding tissues and damage to these tissues. A manipulation may also cause a loss of or change to the treatment, which may lead to a dosage error or an insufficient drug dose with clinically significant

consequences, and require further surgery or replacement.

- The port and the catheter are subject to wear. Over time, a defect in the port or catheter may occur and require further surgery or replacement. A component defect may cause leakage of the drug or CSF into the surrounding tissues and damage to these tissues or a loss of treatment. This may lead to a dosage error, an insufficient drug dose or a drug overdose with clinically significant consequences.
- Before undergoing examinations or medical procedures, the patient or the family/the patient's legal representative must inform all the health personnel that the patient is fitted with an implantable port. If the health personnel is not informed, this may cause delays in procedures, lesions in the patient, component damage requiring further surgery or replacement, or lead to a dosage error, insufficient drug dose or drug overdose with clinically significant consequences.

9. Replacement of the SOPH-A-PORT port and catheter

9.1. REPLACEMENT OF THE PORT

Several causes may lead to a port defect, in particular:

- A septum damaged by too many punctures
- Punctures performed with a perforating needle or too big a Huber needle
- Disconnection of the port and catheter.

If the port does not work but the catheter works, the port may be replaced whilst keeping the catheter in place. Observe these recommendations if the port is defective but the catheter remains entirely functional.

WARNINGS

The SOPH-A-PORT should only be replaced by qualified doctors.

Use an aseptic technique.

During explantation there may be fibrous tissue on the port implantation site. The surgeon should adopt an appropriate technique to free the port.

1. Once the anaesthetic has been administered, place the patient as described in Section 4 "SOPH-A-PORT implantation instructions".
2. Prepare the port implantation site with antiseptic solution. Allow the antiseptic solution to dry.
3. Drape the area leaving the port implantation site exposed.
4. Position the patient's hips, legs and arms so that pressure points are prevented.
5. Make an incision and carefully dissect the subcutaneous tissues to expose the existing port and catheter.
6. Clamp the catheter close to the port with rubber-tipped forceps.
7. Cut the catheter between the forceps and the outlet pin to separate the catheter from the port.
8. Remove the sutures from the base of the port.
9. Remove the port from the pocket with the cut end of

the catheter.

Put the port with the fragment of catheter aside with a view to returning them to the manufacturer.

10. Remove the rubber-tipped forceps. Observe a fluid reflux to confirm the catheter is patent.
11. Put the rubber-tipped forceps back in place until the catheter is connected to a new port.
12. Continue the implantation of the new port following the steps in Section 4 "SOPH-A-PORT implantation instructions".
13. Return the port, the mobile connector with pin and the fragment of catheter to the manufacturer. Follow the instructions in Section 12 "Return to Sophysa".

9.2. REPLACEMENT OF THE CATHETER

NOTE

If the catheter does not work the entire system must be replaced.

Several causes may lead to a catheter defect, in particular:

- Loss of patency
- A leak
- Dislodgement of the catheter.

WARNINGS

The SOPH-A-PORT should only be replaced by qualified doctors.

Use an aseptic technique.

1. Once the anaesthetic has been administered, place the patient as described in Section 4 "SOPH-A-PORT implantation instructions".
2. Consult the instructions in Section 9.1 "Replacement of the port" above, to access the port and remove it.
3. Withdraw the catheter from the subcutaneous region exercising slight traction and using surgical dissection techniques as required. Put the catheter aside with a view to returning it to the manufacturer.
4. If at this stage using the port location for reimplantation is not planned, confirm haemostasis and close the incision by placing sutures.
5. Return the catheter and port to the manufacturer. Follow the instructions in Section 12 "Return to Sophysa".
6. If implanting a new SOPH-A-PORT is planned, follow the instructions in Section 4 "SOPH-A-PORT implantation instructions".

10. Patient identification card

A Patient Identification Card is provided with each SOPH-A-PORT. After filling in the card the physician should give it to the patient who should be advised him/her to carry it all times.

This card contains useful information about SOPH-A-PORT.

11. Treatment of the product after use

It is possible that the device has to be removed for the

following reasons:

- The device should be replaced
- The patient no longer needs the treatment
- There is an infection around the components of the implanted device
- A skin erosion has occurred on the device or the wound has deteriorated.

When a device is explanted or when the components are discarded, observe the following instructions:

WARNINGS

The SOPH-A-PORT port and catheter are **single use only**. Do not re-sterilise and do not re-implant explanted devices.

Never re-use an implantable device or an implantable accessory after exposure to tissues or biological fluids.

Return of the explanted components of the SOPH-A-PORT device

Please return all the explanted components as they are, handling or modifying them as little as possible.

If a SOPH-A-PORT has to be returned to Sophysa for evaluation:

- immerse the components in water (without saline)
- Indicate whether the components have been cleaned.

To enable the device to be evaluated do not autoclave or expose it to ultrasound cleaning appliances.

So that Sophysa can evaluate the returned product properly, it should be accompanied by a **Return to Manufacturer Authorization** form.

Destruction after use: Dispose of any component which is not returned in accordance with local environmental regulations.

12. Return to Sophysa

Ports and attached catheters which are explanted due to a suspected malfunction should be returned to Sophysa for analysis immersed in sterile water.

The doctor should contact Sophysa before returning them.

13. Warranty

Sophysa only guarantees to the initial purchaser of a SOPH-A-PORT for implantation the free replacement of this system in the case of failure for any reason attributable to defects in materials and manufacture. This guarantee is valid for 90 days from the date of implantation, subject to the following conditions:

- The use of the system must be in accordance with Instruction Manual which is included with the product
- The explanted system must be returned to Sophysa within 30 days after explantation, and becomes the property of Sophysa.

Obstruction of the catheter resulting in the explantation of the system is not considered as a failure of the device, and therefore is not covered by the clauses of this limited warranty.

This limited warranty applies only to the use of the system for the initial patient and the sole liability of Sophysa is

limited to the replacement of the system under the conditions stipulated here.

Sophysa does not declare or guarantee that there will be no failure of the system, that the body will not react adversely to the implantation of the system or that the medication will be safe and effective for the treatment of a patient's condition; Sophysa disclaims all responsibility for any medical complications, including death, resulting from the use of the system.

This warranty expressly replaces all other warranties. Sophysa also declines all responsibility for any general, special, secondary or immaterial damages arising from the sale, manufacture or use of the product sold under this warranty.

Sophysa gives no express or implicit guarantees, including, but not limited to, any warranty of merchantability or fitness of the products for any purpose or reason, for the SOPH-A-PORT implantable drug delivery system sold under this warranty, apart from other provisions expressly stipulated here.

14. Symbols

REF

Catalogue reference

Manufacturer

Date of manufacture

Caution, see the Instructions for Use

STERILE EO

Sterilized using Ethylene Oxide

Do not re-use

Do not re-sterilize

Use before

LOT

Batch code

SN

Serial number

CE Conformity Marking

15. References

15.1. VENOUS

SP-20110

SOPH-A-PORT® Standard venous unattached

Port. Base: 33x27 mm. Height: 13 mm. Ø Septum: 13 mm. Septum thickness: 6.5 mm. Residual volume: 0.45 ml. Weight: 8.2 g.

Catheter. ID: 1.1 mm. external Ø: 2.2 mm. Internal volume: 0.113 ml /100 mm. Length: 700 mm. Material: Silicone XRO. Graduations at 5, 10, 15, 20 cm.

SP-10110

SOPH-A-PORT® Mini venous unattached

Port. Base: 26x21 mm. Height: 10 mm. Ø Septum: 10 mm. Septum thickness: 4.8 mm. Residual volume: 0.24 ml. Weight: 3.8 g.

Catheter. ID: 1.1 mm. external Ø: 2.2 mm. Internal volume: 0.113 ml /100 mm. Length: 700 mm. Material: Silicone XRO. Graduations at 5, 10, 15, 20 cm.

15.2. SPINAL

SP-20210

SOPH-A-PORT® Standard spinal unattached

Port. Base: 33x27 mm. Height: 13 mm. Ø Septum: 13 mm. Septum thickness: 6.5 mm. Residual volume: 0.45 ml. Weight: 8.2 g.

Catheter. ID: 0.76 mm. external Ø: 1.65 mm. Internal volume: 0.058 ml /100 mm. Length: 800 mm. Material: Silicone XRO. Graduations at 11, 16, 21, 26 cm.

SP-10210

SOPH-A-PORT® Mini spinal unattached

Port. Base: 26x21 mm. Height: 10 mm. Ø Septum: 10 mm. Septum thickness: 4.8 mm. Residual volume: 0.24 ml. Weight: 3.8 g.

Catheter. ID: 0.76 mm. external Ø: 1.65 mm. Internal volume: 0.058 ml /100 mm. Length: 800 mm. Material: Silicone XRO. Graduations at 11, 16, 21, 26 cm.

15.3. ARTERIAL

SP-20310

SOPH-A-PORT® Standard arterial unattached

Port. Base: 33x27 mm. Height: 13 mm. Ø Septum: 13 mm. Septum thickness: 6.5 mm. Residual volume: 0.45 ml. Weight: 8.2 g.

Catheter. ID: 0.6 mm. external Ø: 2.2 mm. Internal volume: 0.038 ml /100 mm. Length: 500 mm. Material: Silicone XRO. Fixing rings at 5, 6, 8 cm.

SP-10310

SOPH-A-PORT® Mini arterial unattached

Port. Base: 26x21 mm. Height: 10 mm. Ø Septum: 10 mm. Septum thickness: 4.8 mm. Residual volume: 0.24 ml. Weight: 3.8 g.

Catheter. ID: 0.6 mm. external Ø: 2.2 mm. Internal volume: 0.038 ml /100 mm. Length: 500 mm. Material: Silicone XRO. Fixing rings at 5, 6, 8 cm.

15.4. PERITONEAL

SP-20410

SOPH-A-PORT® Standard peritoneal unattached

Port. Base: 33x27 mm. Height: 13 mm. Ø Septum: 13 mm. Septum thickness: 6.5 mm. Residual volume: 0.45 ml. Weight: 8.2 g.

Catheter. ID: 2.5 mm. external Ø: 4.0 mm. Internal volume: 0.573 ml /100 mm. Length: 500 mm. Material: Silicone XRO. Multi-perforation open end.

SP-20411

SOPH-A-PORT® Standard peritoneal preconnected

Port. Base: 33x27 mm. Height: 13 mm. Ø Septum: 13 mm. Septum thickness: 6.5 mm. Residual volume: 0.45 ml. Weight: 8.2 g.

Catheter. ID: 2.5 mm. external Ø: 4.0 mm. Internal volume: 0.573 ml /100 mm. Length: 500 mm. Material: Silicone XRO. Multi-perforation open end.

SP-10410

SOPH-A-PORT® Mini peritoneal unattached

Port. Base: 26x21 mm. Height: 10 mm. Ø Septum: 10 mm. Septum thickness: 4.8 mm. Residual volume: 0.24 ml. Weight: 3.8 g.

Catheter. ID: 2.5 mm. external Ø: 4.0 mm. Internal volume: 0.573 ml /100 mm. Length: 500 mm. Material: Silicone XRO. Multi-perforation open end.

SP-10411

SOPH-A-PORT® Mini peritoneal preconnected

Port. Base: 26x21 mm. Height: 10 mm. Ø Septum: 10 mm. Septum thickness: 4.8 mm. Residual volume: 0.24 ml. Weight: 3.8 g.

Catheter. ID: 2.5 mm. external Ø: 4.0 mm. Internal volume: 0.573 ml /100 mm. Length: 500 mm. Material: Silicone XRO. Multi-perforation open end.

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SOPH-A-POR^T® is a Sophysa SA registered trademark. - ©2019 Sophysa SA - All rights reserved.