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ENGLISH

FIGURE 1 - DE-108 EXTERNAL LUMBAR CATHETER.

A :EXTERNAL LUMBAR CATHETER. B : FIXATION TAB. C : TUOHY NEEDLE. D : COMPRESSION HUB (5FR).


FIGURE 2 : DE-109 EXTERNAL LUMBAR CATHETER WITH INTRODUCER GUIDE.

A : EXTERNAL LUMBAR CATHETER. B : FIXATION TAB. C : TUOHY NEEDLE. D : COMPRESSION HUB (5FR). E : INTRODUCER GUIDE.


Table of contents

1. Indications	3
2. Contra-indications	3
3. Device description	3
4. Sterilization - Decontamination	4
5. Patient education	4
6. Instructions	4
6.1. Drainage catheter placement	4
6.2. To connect catheter to drainage system	5
7. Complications	6
7.1. Infections	6
7.2. Overdrainage	6
7.3. Obstruction	6
8. Guarantee	6
9. Processing of the products after use	6
9.1. Destruction after use	6
9.2. Return of products	6
10. Symbols	6
11. References	7

WARNINGS :

Read all instructions, precautions and warnings prior to use.

Federal Law (USA) restricts this device to sale by or on the order of a physician.

1. Indications

These devices are intended for use as a means of temporary external diversion of CSF.

The main indication for drainage and monitoring of CSF is the treatment of raised intracranial pressure when the ventricles or sub-arachnoid spaces are accessible.

External drainage can also be used for the control of hydrocephalus in cases where an internal drainage system is infected.

External drainage can also be used either prior to implanting an internal shunt, or in an attempt to avoid the need for it.

External ventricular drainage is indicated in acute or sub-acute hydrocephalus (secondary to intra ventricular or massive meningeal haemorrhages, severe meningitis, etc.)

2. Contra-indications

This device is not intended for any use other than those indicated.

Intracranial pressure monitoring with a lumbar catheter is contraindicated in patients receiving anticoagulants or who are known to have a bleeding diathesis.

A patient undergoing external drainage and monitoring must be kept under continuous, close supervision especially to prevent complications (for more information see complications section). The use of a ventricular or lumbar drainage catheter is contraindicated where trained personnel are not available to supervise monitoring and drainage on a 24-hour-a-day basis.

The use of a lumbar catheter for drainage and monitoring of cerebrospinal fluid is not recommended :

- for patients with non-communicating hydrocephalus;
- for patients where lumbar puncture is contraindicated;
- in the presence of large intracranial mass lesions, tumors, hematomas, or cysts;
- in the presence of infections in the surrounding area which includes the skin, subcutaneous tissue, bone and the epidural space;
- for patients which have demonstrated blockage of cerebrospinal fluid to the subarachnoid space due to trauma, hematoma, fracture or tumor.

The use of a lumbar catheter under these conditions for external drainage and monitoring is at the discretion of the physician.

3. Device description

Sophysa external drainage lumbar catheters are available as individual components or included in complete drainage

system kits. Sophysa lumbar catheters are intended to be used with Sophysa drainage systems.

DE-108 External Lumbar Catheter (Figure1) and **DE-109 External Lumbar Catheter with Introducer Guide** (Figure 2) made of silicon tubing are radiopaque, 90 cm long and have an outer diameter of 1.6 mm and an inner diameter of 0.76 mm. They are multiperforated over a distance of 23 mm from the proximal tip. Each has depth markings at 11, 16, 21 and 26 cm from the proximal tip. DE-108 External Lumbar Catheter has an opened-end tip. DE-109 External Lumbar Catheter has a closed-end tip.

Each lumbar catheter is supplied with a fixation tab (for catheter fixation to the cutaneous tissue), a Tuohy needle (14G), 9 cm length (for introduction into the intrathecal space) and a compression hub 5Fr (for connecting the catheter to the tubing set). DE-109 Lumbar catheter is supplied with an introducer guide.

4. Sterilization - Decontamination

Products are packed individually in double peel-off, sterile, pyrogen-free packaging, sterilized with ethylene oxide.

WARNINGS :

Do not use the products if sterile packaging is open or damaged, or if the expiration date has passed.

This product is intended for SINGLE USE ONLY. It is intended to be used once only for a single patient. Do not re-sterilize or re-use after unpacking and/or explantation.

Resterilization can damage product, potentially leading to patient injury. Re-use of this device may change its mechanical or biological features and may cause device failure, allergic reactions or bacterial infections.

NOTE:

Sophysa cannot be held responsible for the performance of any product that has been re-sterilized, nor for any complications which might result from this.

5. Patient education

It is physician's responsibility to educate the patient and /or their representative(s) regarding external drainage and monitoring. This should include a description of associated complications, and explanation of potential alternative products and treatments.

6. Instructions

WARNINGS :

In order to minimize the possibilities of infection, meningitis or ventriculitis, aseptic technique is necessary in all phases of the use of this device.

It is imperative that this device not be handled with bare fingers or come into contact with lint. Silicone elastomer is highly electrostatic and therefore susceptible to contamination by airborne or surface particles, thus minimize exposure of silicone to air. The presence of these contaminants could cause adverse tissue reaction. Rubber-shod clamps or washed, gloved hands are the best means of handling implantable, silicone devices.

6.1. DRAINAGE CATHETER PLACEMENT

WARNING:

In order to ensure against ventricular collapse and the possible consequence of tentorial herniation, care should be taken so as little CSF as possible is lost during the insertion of the catheter when the lumbar subarachnoid space is first punctured during the insertion of the catheter.

Several techniques can be used to introduce the catheter into the sub-arachnoid spaces. The surgeon should choose the one most suited to the needs of the patient and to his experience.

What follows is given as general guidance for optimal use of the product.

6.1.1. Insertion site preparation

After preparing and draping the patient's lumbar area in a sterile fashion, local anesthetic is injected subcutaneously and intramuscularly in the midline, around the L3-4 or L4-5 area, then a small incision is made where the Tuohy needle will be inserted.

6.1.2. Lumbar catheter insertion

NOTE:

Wet the catheter with a sterile solution to facilitate placement.

DE-108 LUMBAR CATHETER INSERTION:

1. Insert the Tuohy needle in the interspinous space in the midline with the bevel parallel to the spine to make puncture easier (Cf. Figure 3 A).
2. Once ensured that the tip of the Tuohy needle is placed correctly, orient the bevel cephalad to allow passage of the catheter (Cf. Figure 3 B).
3. Remove the stylet from the Tuohy needle. A sampling of CSF can be done at this level.
4. Introduce the catheter into the Tuohy needle lumen with rubber-covered clamps. The catheter must go at least 10 cm beyond the tip of the Tuohy needle : there is only one black mark left at the Tuohy needle entry.
5. Carefully remove the Tuohy needle over the catheter.

DE-109 LUMBAR CATHETER INSERTION WITH INTRODUCER GUIDE:


1. Insert the Tuohy needle in the interspinous space in the midline with the bevel parallel to the spine to make puncture easier (Cf. Figure 3 A).
2. Once ensured that the tip of the Tuohy needle is placed correctly, the bevel is oriented cephalad to allow passage of the catheter (Cf. Figure 3 B).
3. Insert the guidewire into the lumen of the lumbar catheter, up to where the tip of the guidewire gets in contact with the closed tip of the lumbar catheter.
4. Remove the stylet from the Tuohy needle. A sampling of CSF can be done at this level.
5. Introduce the catheter into the Tuohy needle lumen. The catheter must go at least 10 cm beyond the tip of the Tuohy needle : there is only one black mark left at the Tuohy needle entry.
6. Carefully remove the Tuohy needle over the catheter.

- Remove the guidewire.

PRECAUTION :

If the catheter must be removed, withdraw the catheter and Tuohy needle simultaneously to avoid damaging or cutting the catheter.

FIGURE 3 - TUOHY NEEDLE INSERTION.

6.1.3. Lumbar catheter tunnelization


WARNING:

To minimize the possibilities of infection tunneling of the catheter should be approximately 3 to 5 cm.

Before tunnelization place again the stylet in the Tuohy needle.

- Place again the stylet in the Tuohy needle and introduce the needle from where the catheter should emerge towards its present location (Cf. Figure 4).

FIGURE 4 - TUNNELIZATION WITH THE TUOHY NEEDLE.

- Remove the stylet.
- Carefully introduce the lumbar catheter into the tip of the Tuohy needle with rubber-covered clamps and then carefully pull it out of the Tuohy needle.
- Carefully remove The Tuohy needle.
- Confirm the CSF is flowing and then adjust the lumbar catheter length by cutting it close to its exit point.
- Clamp the catheter.

6.1.4. Compression hub connection

- Insert and fix the compression hub into the distal tip of the catheter (Figures 5 and 6).
- The catheter is fixed to the cutaneous tissue using the fixation tab.

FIGURE 5 - COMPRESSION HUB. THE COMPRESSION HUB IS AN ALTERNATIVE TO THE CLASSICAL LUER-LOCK CONNECTOR. (A) CATHETER GUIDE. (B) COMPRESSION FITTING. (C) DISTAL CAP.



FIGURE 6 - TO USE A COMPRESSION HUB. (1) REMOVE THE COMPRESSION FITTING. (2) PASS THE LUMBAR CATHETER DISTAL END THROUGH THE COMPRESSION FITTING. (3) INSERT THE COMPRESSION HUB CATHETER GUIDE INTO THE LUMBAR CATHETER DISTAL END AND PUSH IT AS FAR AS POSSIBLE. (4) SCREW THE COMPRESSION FITTING.

6.2. TO CONNECT CATHETER TO DRAINAGE SYSTEM

WARNINGS:

Read all instructions, precautions and warnings mentioned in drainage system instructions for use, prior to connect catheter to drainage system.

Ensure that the drainage system patient line is equipped with a luer connector.

In order to ensure against ventricular collapse and the possible consequence of tentorial herniation, always perform a drainage set up against a positive pressure head on the order of 20 cm H2O or 15 mm Hg.

PRECAUTION:

Ensure flow of CSF from the distal end of the shunt system prior to closure.

After the catheter has been properly placed and the compression hub fixed into the distal tip of the catheter :

- Ensure the lumbar catheter is clamped at the exit site to allow only a minimal amount of CSF to escape during connector insertion and connection to system,
- Remove the end plug from the patient line,
- Attach the compression hub to the patient line,
- After ensuring there is no bubbles in the patient line set the patient line stopcock in the drain only position,

- Remove the clamp from the catheter.

WARNINGS:

Prior to use the drainage system :

- Read all instructions, precautions and warnings mentioned in its instructions for use.
- Ensure there is no leak on the tubing and at the connection between the drainage catheter and the drainage system.

7. Complications

There are three main complications associated with CSF external drainage system.

7.1. INFECTIONS

One major complication associated with ICP monitoring with a ventricular or lumbar catheter is the risk of infections, particularly meningitis and ventriculitis. The incidence of these infections can be reduced by care in inserting the catheter and stabilizing it:

- by passing it through a subgaleal (ventricular catheter) or subcutaneous (lumbar catheter) tunnel before it emerges;
- by fixing it to the scalp (ventricular catheter) or cutaneous tissue (lumbar catheter) with the fixation tab.

Wound infections may occur but usually subside when the catheter is removed.

Limiting the duration of monitoring from a single site to less than five days will reduce the infection rate. If monitoring must continue past five days, consider inserting a new catheter at a fresh site and changing the entire system.

In case of infection, the drainage system should be removed and specific treatment by general or intrathecal routes should be introduced.

7.2. OVERDRAINAGE

A second major complication associated with ventricular or lumbar drainage of CSF is overdrainage, which can lead to intracranial hemorrhage and permanent neurological deficit. Overdrainage can occur due to improper system testing or setup (resulting in system leakage or inappropriate system pressure head heights).

7.3. OBSTRUCTION

A third major complication associated with ventricular or lumbar drainage of CSF is obstruction. Poor recording of ICP will result if the catheter, patient line or other components of the monitoring system become clogged with blood clots, brain tissue fragments, or fibrinous debris.

Frequent punctures of the brain to insert the ventricular catheter can predispose to intracerebral hemorrhage and edema causing a further rise in ICP.

In patients with small ventricles, the ventricle walls may collapse around the tip of the catheter resulting in obstruction and predisposing to tentorial herniation. It is therefore extremely important to avoid excessive release of CSF before the catheter is attached to the patient line.

8. Guarantee

Sophysa guarantees that this medical device is free from defects in material or manufacturing. Apart from this guarantee, Sophysa does not grant any other guarantee, express or implicit, including marketing or adaptation for a specific use. Sophysa cannot be held responsible for any incident, complication, damage or prejudice resulting directly or indirectly from the use of this device. Sophysa does not authorize anyone to take responsibility on its behalf for its products.

The performance of lumbar catheter is only guaranteed with the range of devices and accessories designed, tested and manufactured by Sophysa. However, it is possible to use other brands of drainage system provided that they are suitable with Sophysa lumbar catheter.

9. Processing of the products after use

9.1. DESTRUCTION AFTER USE

An unpacked, used or explanted lumbar catheter must be destroyed in accordance with the procedures in force in the medical establishment.

9.2. RETURN OF PRODUCTS

If an explanted lumbar catheter needs to be returned to Sophysa for analysis, it must be returned immersed in water (non-saline), indicating if necessary whether cleaning has been performed.

In order to properly assess the returned product, it must be accompanied by an explanatory Return to Manufacturer Authorization form.

10. Symbols

Manufacturer

REF

Catalog number

SN

Serial number

LOT

Batch code

Use by

Caution, consult accompanying documents

STERILE EO

Sterilized using ethylene oxide

Do not resterilize

Do not reuse

CE conformity marking

11. References

Lumbar catheter for external drainage of CSF

DE-108

External lumbar catheter

DE-109

External lumbar catheter with introducer guide

External drainage systems

DE-210

External neurological drainage system

DE-410

External drainage and CSF monitoring system

Complete kit for external drainage

DE-308

Complete kit for external lumbar drainage (DE-210+DE-108)

DE-408

Complete kit for external lumbar drainage and CSF monitoring (DE-410+DE-108)

Product Description Subject To Change

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