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ENGLISH

FIGURE 1 - DE-106 EXTERNAL VENTRICULAR CATHETER.

A : EXTERNAL VENTRICULAR CATHETER. B : FIXATION TAB. C : STYLET. D : TROCAR. E : LUER-LOCK CONNECTOR.


FIGURE 2 : DE-105 VENTRICULAR CATHETER (LARGE LUMEN).

A : EXTERNAL VENTRICULAR CATHETER (LARGE LUMEN). B : FIXATION TAB. C : STYLET. D : TROCAR. E : LUER-LOCK CONNECTOR.


Table of contents

1. Indications	3
2. Contra-indications	3
3. Device description	3
4. Sterilization - Decontamination	4
5. Patient education	4
6. Instructions	4
6.1. Drainage catheter placement	4
6.2. To connect catheter to drainage system	5
7. Complications	5
7.1. Infections	5
7.2. Overdrainage	5
7.3. Obstruction	5
8. Guarantee	6
9. Processing of the products after use	6
9.1. Destruction after use	6
9.2. Return of products	6
10. Symbols	6
11. References	6

WARNINGS :

Read all instructions, precautions and warnings prior to use.

Federal Law (USA) restricts this device to sale by or on the order of a physician.

1. Indications

These devices are intended for use as a means of temporary external diversion of CSF.

DE-105 External Ventricular Catheter (large lumen) is particularly recommended in the cases of hemorrhagic CSF or subdural hematomas.

The main indication for drainage and monitoring of CSF is the treatment of raised intracranial pressure when the ventricles or sub-arachnoid spaces are accessible.

External drainage can also be used for the control of hydrocephalus in cases where an internal drainage system is infected.

External drainage can also be used either prior to implanting an internal shunt, or in an attempt to avoid the need for it.

External ventricular drainage is indicated in acute or sub-acute hydrocephalus (secondary to intra ventricular or massive meningeal haemorrhages, severe meningitis, etc.)

2. Contra-indications

This device is not intended for any use other than those indicated.

Intracranial pressure monitoring with a ventricular catheter is contraindicated in patients receiving anticoagulants or who are known to have a bleeding diathesis.

A patient undergoing external drainage and monitoring must be kept under continuous, close supervision especially to prevent complications (for more information see complications section). The use of a ventricular or lumbar drainage catheter is contraindicated where trained personnel are not available to supervise monitoring and drainage on a 24-hour-a-day basis.

The use of a ventricular catheter is contraindicated if infection is present.

3. Device description

Sophysa external drainage ventricular catheters are available as individual components or included in complete drainage system kits. Sophysa ventricular catheters are intended to be used with Sophysa drainage systems.

DE-106 External Ventricular Catheter (Figure1) is a transparent silicon tubing with a radiopaque stripe, 29 cm long and has an outer diameter of 3.0 mm and an inner diameter of 1.5 mm. This ventricular catheter is multiperforated over a distance of 23 mm from the proximal tip and has depth markings at 5, 10 and 15 cm from the proximal tip.

DE-105 External Ventricular Catheter (large lumen) (Figure 2) is a barium impregnated radiopaque silicon tubing, 35 cm long and has an outer diameter of 4.0 mm and an inner diameter of 2.3 mm. This ventricular catheter (large lumen) is multiperforated with 2.0 mm holes over a distance of 25 mm from the proximal tip and has depth markings at 5, 7.5, 10, 15 and 20 cm from the proximal tip.

Each ventricular catheter is supplied with a female luer-lock connector (for connecting the catheter to the tubing set), an introducer stylet (for introduction), a trocar (for tunneling), and a fixation tab (for catheter fixation to the scalp). DE-105 External Ventricular Catheter fixation tab is pre-attached to the catheter at 25 cm from the proximal tip.

4. Sterilization - Decontamination

Products are packed individually in double peel-off, sterile, pyrogen-free packaging, sterilized with ethylene oxide.

WARNINGS :

Do not use the products if sterile packaging is open or damaged, or if the expiration date has passed.

This product is intended for SINGLE USE ONLY. It is intended to be used once only for a single patient. Do not re-sterilize or re-use after unpacking and/or explantation.

Resterilization can damage product, potentially leading to patient injury. Re-use of this device may change its mechanical or biological features and may cause device failure, allergic reactions or bacterial infections.

NOTE:

Sophysa cannot be held responsible for the performance of any product that has been re-sterilized, nor for any complications which might result from this.

5. Patient education

It is physician's responsibility to educate the patient and /or their representative(s) regarding external drainage and monitoring. This should include a description of associated complications, and explanation of potential alternative products and treatments.

6. Instructions

WARNINGS :

In order to minimize the possibilities of infection, meningitis or ventriculitis, aseptic technique is necessary in all phases of the use of this device.

It is imperative that this device not be handled with bare fingers or come into contact with lint. Silicone elastomer is highly electrostatic and therefore susceptible to contamination by airborne or surface particles, thus minimize exposure of silicone to air. The presence of these contaminants could cause adverse tissue reaction. Rubber-shod clamps or washed, gloved hands are the best means of handling implantable, silicone devices.

6.1. DRAINAGE CATHETER PLACEMENT

WARNING:

In order to ensure against ventricular collapse and the possible consequence of tentorial herniation, care should be taken so as little CSF as possible is lost during the insertion of the catheter when the ventricle is first punctured during the insertion of the catheter.

Several techniques can be used to introduce the catheter into the ventricular spaces. The surgeon should choose the one most suited to the needs of the patient and to his experience.

What follows is given as general guidance for optimal use of the product

6.1.1. Insertion site preparation

Once the implantation site has been chosen, the area is shaved and prepared aseptically.

An incision of about one centimeter is made down to the bone. Make sure the bony plate is well exposed and perform hemostasis on the wound edges. Perforation is performed through the internal and external skull plates.

6.1.2. Ventricular catheter insertion

1. Puncture the dura (with a lumbar puncture needle for instance).
2. Implant the ventricular catheter in the direction of the ventricle thanks to its pre-inserted stylet. To get to the ventricle, catheter is guided in the direction of inner canthus of the opposite eye (Cf. Figure 3). After the catheter has been advanced about 5-6 cm, adequate flow of CSF is typically obtained
3. Once the catheter is properly positioned in the ventricle, remove the stylet while holding the catheter in place. A sampling of CSF can be done at this level.

FIGURE 3 - VENTRICULAR CATHETER INSERTION.

6.1.3. Ventricular catheter tunnelization


WARNING:

To minimize the possibilities of infection tunneling of the catheter should be approximately 3 to 5 cm.

1. Insert the trocar in the catheter (Cf. Figure 4 A).
2. Inserted the trocar between the scalp and the skull from the burr hole towards the emergence site (Cf. Figure 4 B).
3. Pull the catheter out of the tunnel, taking good care of maintaining it in place in the ventricle,
4. Carefully remove the trocar

5. Adjust catheter length by cutting it close to its exit point.
6. Clamp the catheter.

FIGURE 4 - VENTRICULAR CATHETER TUNNELIZATION.


6.1.4. Luer connector fixation

Luer-lock connecting fitting use:

Insert Luer connector into the distal tip of the catheter and suture it to the catheter (Cf Figure 5).

FIGURE 5 - LUER-LOCK CONNECTING FITTING USE.

PRECAUTION:

A double suture tie should be used to secure the ventricular catheter fixation to the connecting fitting. Check to ensure that the connection is tight prior to use. Do not tie ligatures too tightly or use non metallic ligatures since silicone elastomer is easily cut.

The catheter is fixed to the scalp using the fixation tab.

6.2. TO CONNECT CATHETER TO DRAINAGE SYSTEM

WARNINGS:

Read all instructions, precautions and warnings mentioned in drainage system instructions for use, prior to connect catheter to drainage system.

Ensure that the drainage system patient line is equipped with a luer connector.

In order to ensure against ventricular collapse and the possible consequence of tentorial herniation, always perform a drainage set up against a positive pressure head on the order of 20 cm H2O or 15 mm Hg.

PRECAUTION:

Ensure flow of CSF from the distal end of the shunt system prior to closure.

After the catheter has been properly placed and the Luer connector fixed into the distal tip of the catheter :

- Ensure the ventricular catheter is clamped at the exit site to allow only a minimal amount of CSF to escape during connector insertion and connection to system,
- Remove the end plug from the patient line,
- Attach the Luer connector to the patient line,
- After ensuring there is no bubbles in the patient line set the patient line stopcock in the drain only position and remove the clamp from the catheter.

WARNINGS:

Prior to use the drainage system :

- Read all instructions, precautions and warnings mentioned in its instructions for use.
- Ensure there is no leak on the tubing and at the connection between the drainage catheter and the drainage system.

7. Complications

There are three main complications associated with CSF external drainage system.

7.1. INFECTIONS

One major complication associated with ICP monitoring with a ventricular or lumbar catheter is the risk of infections, particularly meningitis and ventriculitis. The incidence of these infections can be reduced by care in inserting the catheter and stabilizing it:

- by passing it through a subgaleal (ventricular catheter) or subcutaneous (lumbar catheter) tunnel before it emerges;
- by fixing it to the scalp (ventricular catheter) or cutaneous tissue (lumbar catheter) with the fixation tab.

Wound infections may occur but usually subside when the catheter is removed.

Limiting the duration of monitoring from a single site to less than five days will reduce the infection rate. If monitoring must continue past five days, consider inserting a new catheter at a fresh site and changing the entire system.

In case of infection, the drainage system should be removed and specific treatment by general or intrathecal routes should be introduced.

7.2. OVERDRAINAGE

A second major complication associated with ventricular or lumbar drainage of CSF is overdrainage, which can lead to intracranial hemorrhage and permanent neurological deficit. Overdrainage can occur due to improper system testing or setup (resulting in system leakage or inappropriate system pressure head heights).

7.3. OBSTRUCTION

A third major complication associated with ventricular or lumbar drainage of CSF is obstruction. Poor recording of ICP will result if the catheter, patient line or other components of the monitoring system become clogged with blood clots, brain tissue fragments, or fibrinous debris.

Frequent punctures of the brain to insert the ventricular catheter can predispose to intracerebral hemorrhage and edema causing a further rise in ICP.

In patients with small ventricles, the ventricle walls may collapse around the tip of the catheter resulting in obstruction and predisposing to tentorial herniation. It is therefore extremely important to avoid excessive release of CSF before the catheter is attached to the patient line.

8. Guarantee

Sophysa guarantees that this medical device is free from defects in material or manufacturing. Apart from this guarantee, Sophysa does not grant any other guarantee, express or implicit, including marketing or adaptation for a specific use. Sophysa cannot be held responsible for any incident, complication, damage or prejudice resulting directly or indirectly from the use of this device. Sophysa does not authorize anyone to take responsibility on its behalf for its products.

The performance of ventricular catheter is only guaranteed with the range of devices and accessories designed, tested and manufactured by Sophysa. However, it is possible to use other brands of drainage system provided that they are suitable with Sophysa ventricular catheter.

9. Processing of the products after use

9.1. DESTRUCTION AFTER USE

An unpacked, used or explanted ventricular catheter must be destroyed in accordance with the procedures in force in the medical establishment.

9.2. RETURN OF PRODUCTS

If an explanted ventricular catheter needs to be returned to Sophysa for analysis, it must be returned immersed in water (non-saline), indicating if necessary whether cleaning has been performed.

In order to properly assess the returned product, it must be accompanied by an explanatory Return to Manufacturer Authorization form.

10. Symbols

Manufacturer

REF

Catalog number

SN

Serial number

LOT

Batch code

Use by

Caution, consult accompanying documents

STERILE EO

Sterilized using ethylene oxide

Do not resterilize

Do not reuse

CE conformity marking

11. References

Ventricular catheter for external drainage of CSF

DE-106

External ventricular catheter

DE-105

External ventricular catheter (large lumen)

External drainage systems

DE-210

External neurological drainage system

DE-410

External drainage and CSF monitoring system

Complete kit for external drainage

DE-306

Complete kit for external ventricular drainage (DE-210+DE-106)

DE-305

Complete kit for external ventricular drainage (large lumen ventricular catheter) (DE-210+DE-105)

DE-406

Complete kit for external ventricular drainage and CSF monitoring (DE-410+DE-106)

Product Description Subject To Change

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