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Table of contents

| | | |
|------|--|---|
| 1. | Indications | 2 |
| 2. | Contra-indications | 2 |
| 3. | Device description | 4 |
| 4. | Sterilization - Decontamination | 4 |
| 5. | Patient education | 4 |
| 6. | Instructions | 4 |
| 6.1. | External drainage | 4 |
| 6.2. | Injections and sampling | 6 |
| 6.3. | ICP monitoring with an external transducer | 6 |
| 7. | Complications | 7 |
| 7.1. | Infections | 7 |
| 7.2. | Overdrainage | 7 |
| 7.3. | Obstruction | 7 |
| 8. | Guarantee | 7 |
| 9. | Processing of the products after use | 7 |
| 9.1. | Destruction after use | 7 |
| 9.2. | Return of products | 7 |
| 10. | Symbols | 7 |
| 11. | References | 8 |

1. Indications

WARNINGS :

Read all instructions, precautions and warnings prior to use.

Federal Law (USA) restricts this device to sale by or on the order of a physician.

The main indication for drainage of CSF is the treatment of raised intracranial pressure when the ventricles or sub-arachnoid spaces are accessible.

External drainage can also be used for the control of hydrocephalus in cases where an internal drainage system is infected.

External drainage can also be used either prior to implanting an internal shunt, or in an attempt to avoid the need for it.

External ventricular drainage is indicated in acute or sub-acute hydrocephalus (secondary to intra ventricular or massive meningeal haemorrhages, severe meningitis, etc.)

2. Contra-indications

This device is not intended for any use other than those indicated.

The use of an external drainage system of CSF with lumbar catheter is not recommended :

- for patients with non-communicating hydrocephalus;
- for patients where lumbar puncture is contraindicated;
- in the presence of large intracranial mass lesions, tumors, hematomas, or cysts;
- in the presence of infections in the surrounding area which includes the skin, subcutaneous tissue, bone and the epidural space;
- for patients which have demonstrated blockage of cerebrospinal fluid to the subarachnoid space due to trauma, hematoma, fracture or tumor.

The use of a lumbar catheter under these conditions for external drainage is at the discretion of the physician.

The use of an external drainage system of CSF with ventricular catheter is contraindicated if infection of the scalp is present.


Use of this external drainage system is contraindicated in patients receiving anticoagulant therapy or in patients receiving anticoagulants or who are known to have a bleeding diathesis.

A patient undergoing external drainage and monitoring must be kept under continuous, close supervision especially to prevent complications (for more information see complications section).

The use of a ventricular or lumbar drainage catheter is contraindicated where trained personnel are not available to supervise monitoring and drainage on a 24-hour-a-day basis.

FIGURE 1 - DE-210 EXTERNAL DRAINAGE SYSTEM.

- (1) CLAMP. (2) TUBING (LENGTH 30CM). (3) DRIP CHAMBER. (4) MALE LUER CONNECTOR + PROTECTIVE CAP. (5) FEMALE LUER CONNECTOR + PROTECTIVE CAP. (6) CLAMP. (7) ANTI-REFLUX VALVE. (8) Y-SHAPED INJECTION PORT. (9) GRADUATED SCALE IN MMH₂O. (10) LUER-LOCK CONNECTION. (11) ICP ACCESS + PROTECTIVE CAP. (12) 3-WAY STOPCOCK. (13) MALE LUER CONNECTOR + PROTECTIVE CAP. (14) EMPTYING STOPCOCK WITH REMOVABLE CAP. (15) ADJUSTABLE CORD. (16) DRAINAGE BAG.

3. Device description

Sophysa external drainage system is available as single item or included in complete kits for external CSF drainage. Sophysa external drainage system is intended to be used with Sophysa external drainage catheters.

DE-210 External Drainage system (Figure 1) is a completely hermetic device, without vent, that allows aseptic shunting of CSF, from the ventricular cavities or the lumbar sub-arachnoid spaces to a collecting bag. In addition it allows monitoring of intracranial pressure (ICP), and injections of medication.

The DE-210 External Drainage system is divided in 2 parts (Figure 1).

A 600 ml graduated drainage bag with:

- an internal drip chamber (3), 30 cm tubing length (2) fitted with a clamp (1) and ending in male Luer connector (4),
- an access port with stopcock removable cap (14),
- 360 mmH₂O graduated scale (9).

A drainage tubing, 160 cm length, identified by a green stripe embedded in tubing wall, with female Luer-Lock connector (5), clamp (6), anti-reflux valve (7), Y-shaped injection port (9), 3-way stopcock with a Luer-Lock connector (11), ending in a male Luer-Lock connector (13) for connection to the drainage catheter.

4. Sterilization - Decontamination

Products are packed individually in double peel-off, sterile, pyrogen-free packaging, sterilized with ethylene oxide.

WARNINGS :

Do not use the products if sterile packaging is open or damaged, or if the expiration date has passed.

This product is intended for **SINGLE USE ONLY**. It is intended to be used once only for a single patient. Do not re-sterilize or re-use after unpacking.

Resterilization can damage product, potentially leading to patient injury. Re-use of this device may change its mechanical or biological features and may cause device failure, allergic reactions or bacterial infections.

NOTE:

Sophysa cannot be held responsible for the performance of any product that has been re-sterilized, nor for any complications which might result from this.

5. Patient education

It is physician's responsibility to educate the patient and /or their representative(s) regarding external drainage. This should include a description of associated complications, and explanation of potential alternative products and treatments.

6. Instructions

WARNINGS :

Drainage bags are delivered with the stopcock at the emptying site in the open position because of the method of sterilisation (ETO). In consequence, prior to use the drainage system, it is necessary to make sure that the drainage stopcock is set to the closed position.

6.1. EXTERNAL DRAINAGE

6.1.1. Drainage catheter placement

Implanting a catheter for external drainage is a surgical procedure.

Several techniques can be used to introduce the catheter into the ventricles or the sub-arachnoid spaces. The surgeon should choose the one most suited to the needs of the patient and to his experience.

WARNINGS :

Prior drainage catheter placement, read all instructions, precautions and warnings mentioned in drainage catheter instructions for use.

Ensure that the drainage catheter is equipped with a Luer connector.

6.1.2. Priming the system

The tubing has a dead volume of approximately 15ml, making it necessary to prime the system prior to use. The procedure should take into account the current guidelines on asepsis used in the establishment.

1. Remove the protective caps placed on both ends of the tubing, and on the entry to the drainage bag.
2. Attach the female Luer connector on the distal end of the tubing to the male Luer connector at the entry to the drainage bag.
3. Fill an appropriate size of syringe with sterile saline solution, flush to remove air, and connect it to the 3-way stopcock on the tubing.
4. External drainage systems are delivered with the 3-way stopcock in the position of drainage and monitoring (Figure 2). Set the position of the 3-way stopcock to connect the syringe with the extension tubing intended for connection to the drainage catheter (Figure 3.1). Inject the saline solution until it arrives at the male Luer connector on the tubing.

FIGURE 2 - 3-WAY STOPCOCK WHEN PACKAGED.




FIGURE 3.1 - 3-WAY STOPCOCK IN PRIMING POSITION (1ST STEP).


5. Close the clamp at the distal end of the tubing.
6. Position the 3-way stopcock so that the syringe is connected to the drainage bag (Figure 3.2).

FIGURE 3.2 - 3-WAY STOPCOCK IN PRIMING POSITION (2ND STEP).

7. Re-open the clamp on the tubing, then inject the saline solution until this flows into the drip chamber.
8. Close again the clamp at the distal end of the tubing.
9. Ensure that there is absolutely no air in the tubing, then put the 3-way stopcock in the drainage position (Figure 4).

FIGURE 4 - 3-WAY STOPCOCK IN DRAINAGE POSITION.

6.1.3. Connection of the catheter to the drainage system

After the catheter has been properly placed and the connection fitting fixed into the distal tip of the catheter :

- Ensure the drainage catheter is clamped at the exit site to allow only a minimal amount of CSF to escape during connection to system,
- Attach the catheter connection fitting to the male Luer connector of the patient line,

In order to ensure the system is correctly set up :

1. Check that there is no bubble in the tubing,
2. Ensure the tubing stopcock is positioned in the drain only position (Figure 4),
3. Remove the clamp from the tubing, then from the catheter,
4. Ensure flow of CSF from the drip chamber of the shunt system,
5. Ensure there is no leak on the tubing and at the connection between the drainage catheter and the drainage system.

Once the system has been correctly set up, if the patient has to be moved position the 3-way stopcock to off (Figure 3.2) and close catheter and tubing clamps.

When the patient reaches the desired location, the drainage system can be positioned.

6.1.4. Positioning the drainage bag


WARNINGS :

Before adjusting the height of the drainage bag relative to the ventricles, allow the CSF to flow along the length of the tubing so as to remove any bubble of air.

Intracranial pressure in mmH₂O is controlled by the height of the drip point on the drainage bag relative to the level of the frontal horns of the lateral ventricles of the patient, who is usually in the half upright position (Figure 5). This is valid for both ventricular and lumbar drainage.

Accordingly, to obtain a desired level of intracranial pressure, first suspend the drainage bag by its cord from a perfusion stand, then adjust its height so the desired level of pressure, directly read on the graduated scale, corresponds to the level of the frontal horns of the patient. For example, to maintain an intracranial pressure of 100 mmH₂O (1 mmH₂O = 9.807 Pa), raise the bag to bring the 100 mmH₂O mark on the scale to the assumed height of the frontal horns (the maximum difference between the real drainage pressure in the system and the pressure given by the bag height is +/- 25 mmH₂O).

FIGURE 5 - POSITION OF THE DRAINAGE BAG RELATIVE TO THE PATIENT).

WARNINGS :

Correct positioning of the drainage bag is essential. The intracranial pressure is controlled only by the height of the drainage bag relative to the frontal horns of the lateral ventricles.

PRECAUTIONS :

If a shunt equipped with a pressure regulation valve is temporarily externalised, the position of the drip chamber should take account of the valve operating pressure; if the drop counter chamber is raised any level above the system zero, this will result in a counter-pressure which will add to the resistance of the valve (for example, with a valve set to 90 mmH₂O, in order to maintain and ICP of 110 mmH₂O, the drip chamber must be positioned against the graduation mark 20 on the mmH₂O scale).

NOTE :

The approximate volumes are indicated on the drainage bag with a higher precision for the beginning of the filling (the maximum difference between the volume indicated and the real volume is +/- 5%).

6.1.5. Supervision during external drainage

Once the external drainage system is set up and the drainage catheter connected, following warnings must be taken into account.

WARNINGS :

A patient undergoing an external drainage an/or ICP monitoring must be kept under closed supervision with a qualified staff.

Make regular checks to ensure no kinking or compression of the tubing exists and that the clamps and the 3-way stopcock are open. Any obstruction of the tubing may stop drainage and have related pathological consequences.

Ensure there is no leak on the tubing and at the connection between the drainage catheter and the drainage system.

Regularly check that neither the patient, nor the bag are lowered or raised accidentally. Any change in the height of the drip chamber relative to the patient may lead to either overdrainage or under-drainage and their related pathological consequences.

Any change in the height of the bag should be made only on medical instructions, and only by qualified staff.

Do not allow alcohol or other solvents to come into contact with the tubing or collecting bag. Certain components will not tolerate exposure to such solvents. Malfunction or erasing of the markings of collecting bag could result if such solvents contact the outer surface of the device. Use only sterile saline solution to clean the device

In order to ensure against ventricular collapse and the possible consequence of tentorial herniation, always perform a drainage maneuver against a positive pressure head on the order of 20 cmH₂O or 15 mmHg.

If a patient who is undergoing external drainage needs to be moved, drainage system should be kept upright and correctly aligned with the patient. If it is not possible to move the drainage system in an upright manner, the following steps must be performed:

- Position the 3-way stopcock to off (Figure 3.2),
- Close catheter and tubing clamps,
- Move patient and system as required,
- Realign and readjust system and stopcock and clamp to initiate drainage when patient reaches new location.

6.1.6. To empty the drainage bag

1. Close the proximal clamp on the tubing.
2. Remove the cap from the emptying site.
3. Bring up a recipient for the liquid.
4. Open the emptying stopcock and empty the bag by allowing the liquid to flow out into the recipient.
5. Close the emptying stopcock and replace the cap.
6. Re-open the clamp to re-establish drainage.

6.1.7. To change the drainage bag

1. Close the two clamps located on the tubing on both sides of the Luer-Lock.
2. Disconnect the Luer-lock and remove the full bag.
3. Connect a new drainage bag to the tubing.
4. Correctly position the new bag at the desired height.
5. Open the clamps.

WARNINGS :

If one of the clamps is closed, no drainage takes place. It is essential to check that both clamps are open following a change of drainage bag.

6.2. INJECTIONS AND SAMPLING**WARNINGS :**

Ensure that the products to be injected are compatible with silicone and PVC prior to perform any injections through the port included with the tubing.

To perform injections or sampling check that the 3-way stopcock is in the drainage position.

1. Close the proximal clamp on the tubing.
2. Carefully disinfect the latex puncture site.
3. Perform the injection or sampling.
4. After injection or sampling, re-open the clamp to re-establish drainage.

NOTE :


Any injection should take into account the dead volume between the injection port and the proximal end of the catheter (approximately 4 ml).

6.3. ICP MONITORING WITH AN EXTERNAL TRANSDUCER

ICP may be monitored precisely by a pressure sensor linked to a display and connected to the free access of the 3-way stopcock.

Position the 3-way stopcock so that the transducer is in contact with the CSF, and the drainage is closed (Figure 6).

FIGURE 6 - POSITION OF THE 3-WAY STOPCOCK FOR MONITORING ICP WITH AN EXTERNAL TRANSDUCER.

In this way, monitoring can be carried out, with intermittent momentary interruption of the drainage. Do not forget to re-establish drainage as soon as the pressure has been measured.

When the monitoring is terminated, and the transducer has been removed, place the 3-way stopcock in the drainage

position, and be sure to replace the protective cap on the free access point of the 3-way stopcock.

7. Complications

There are three main complications associated with CSF external drainage system.

7.1. INFECTIONS

One major complication associated with ICP monitoring with a ventricular or lumbar catheter is the risk of infections, particularly meningitis and ventriculitis. The incidence of these infections can be reduced by care in inserting the catheter and stabilizing it:

- by passing it through a subgaleal (ventricular catheter) or subcutaneous (lumbar catheter) tunnel before it emerges;
- by fixing it to the scalp (ventricular catheter) or cutaneous tissue (lumbar catheter) with the fixation tab.

The incidence of these infections can also be reduced by reducing the number of times a drainage system is opened.

Wound infections may occur but usually subside when the catheter is removed.

Limiting the duration of monitoring from a single site to less than five days will reduce the infection rate. If monitoring must continue past five days, consider inserting a new catheter at a fresh site and changing the entire system.

In case of infection, the drainage system should be removed and specific treatment by general or intrathecal routes should be introduced.

7.2. OVERDRAINAGE

A second major complication associated with ventricular or lumbar drainage of CSF is overdrainage, which can lead to intracranial hemorrhage and permanent neurological deficit. Overdrainage can occur due to improper system testing or setup (resulting in system leakage or inappropriate system pressure head heights).

7.3. OBSTRUCTION

A third major complication associated with ventricular or lumbar drainage of CSF is obstruction. Poor recording of ICP will result if the catheter, patient line or other components of the monitoring system become clogged with blood clots, brain tissue fragments, or fibrinous debris.

Frequent punctures of the brain to insert the ventricular catheter can predispose to intracerebral hemorrhage and edema causing a further rise in ICP.

In patients with small ventricles, the ventricle walls may collapse around the tip of the catheter resulting in obstruction and predisposing to tentorial herniation. It is therefore extremely important to avoid excessive release of CSF before the catheter is attached to the patient line.

8. Guarantee

Sophysa guarantees that this medical device is free from defects in material or manufacturing. Apart from this guarantee, Sophysa does not grant any other guarantee, express or implicit, including marketing or adaptation for a specific use. Sophysa cannot be held responsible for

any incident, complication, damage or prejudice resulting directly or indirectly from the use of this device. Sophysa does not authorize anyone to take responsibility on its behalf for its products.

The performance of external drainage system is only guaranteed with the range of devices and accessories designed, tested and manufactured by Sophysa. However, it is possible to use other brands of external drainage catheters provided that they are suitable with Sophysa external drainage system.

9. Processing of the products after use

9.1. DESTRUCTION AFTER USE

An unpacked or used external drainage system must be destroyed in accordance with the procedures in force in the medical establishment.

9.2. RETURN OF PRODUCTS

If the used external drainage system needs to be returned to Sophysa for analysis, indicate if necessary whether cleaning has been performed.

In order to properly assess the returned product, it must be accompanied by an explanatory Return to Manufacturer Authorization form.

10. Symbols

Manufacturer

REF

Catalog number

SN

Serial number

LOT

Batch code

Use by

Caution, consult accompanying documents

STERILE EO

Sterilized using ethylene oxide

Do not resterilize

Do not reuse

CE conformity marking

11. References

Catheters for external drainage of CSF

DE-106

External ventricular catheter

DE-105

External ventricular catheter (large lumen)

DE-108

External lumbar catheter

DE-109

External lumbar catheter with introducer guide

External drainage systems

DE-210

External neurological drainage system

DE-220

Replacement drainage bag (x5)

Complete kit for external drainage

DE-305

Complete kit for external ventricular drainage (large lumen)
(DE-210+DE-105)

DE-306

Complete kit for external ventricular drainage
(DE-210+DE-106)

DE-308

Complete kit for external lumbar drainage (DE-210+DE-108)

Product Description Subject To Change

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