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## 1. Indications

### WARNINGS :

Read all instructions, precautions and warnings prior to use.

Federal Law (USA) restricts this device to sale by or on the order of a physician.

The main indication for drainage and monitoring of CSF is the treatment of raised intracranial pressure when the ventricles or sub-arachnoid spaces are accessible.

External drainage can also be used for the control of hydrocephalus in cases where an internal drainage system is infected.

External drainage can also be used either prior to implanting an internal shunt, or in an attempt to avoid the need for it.

External ventricular drainage is indicated in acute or sub-acute hydrocephalus (secondary to intra ventricular or massive meningeal haemorrhages, severe meningitis, etc.)

## 2. Contra-indications

This device is not intended for any use other than those indicated.

The use of an external drainage and CSF monitoring system with lumbar catheter is not recommended :

- for patients with non-communicating hydrocephalus;

- for patients where lumbar puncture is contraindicated;
- in the presence of large intracranial mass lesions, tumors, hematomas, or cysts;
- in the presence of infections in the surrounding area which includes the skin, subcutaneous tissue, bone and the epidural space;
- for patients which have demonstrated blockage of cerebrospinal fluid to the subarachnoid space due to trauma, hematoma, fracture or tumor.

The use of a lumbar catheter under these conditions for external drainage is at the discretion of the physician.

The use of an external drainage and CSF monitoring system with ventricular catheter is contraindicated if infection of the scalp is present

Use of this external drainage system is contraindicated in patients receiving anticoagulant therapy or in patients receiving anticoagulants or who are known to have a bleeding diathesis.

A patient undergoing external drainage and monitoring must be kept under continuous, close supervision especially to prevent complications (for more information see complications section).

The use of a ventricular or lumbar drainage catheter is contraindicated where trained personnel are not available to supervise monitoring and drainage on a 24-hour-a-day basis.

FIGURE 1 - DE-410 EXTERNAL DRAINAGE AND CSF MONITORING SYSTEM.

- (1) SUPPORT PLATE. (2) LUER CONNECTION + PROTECTIVE CAP. (3) LUER CONNECTION. (4) PATIENT LINE STOPCOCK WITH ACCESS SITE. (5) PATIENT LINE. (6) ANTI-REFLUX VALVE. (7) SYSTEM STOPCOCK WITH ICP ACCESS SITE AND PROTECTIVE CAP. (8) VENTS WITH HYDROPHOBIC FILTERS. (9) 100 ML GRADUATED DRIP CHAMBER. (10) SUSPENSION CORD WITH LOCKING DEVICE. (11) DRIP CHAMBER STOPCOCK WITH ACCESS SITE. (12) CLAMP. (13) DRAINAGE BAG.



FIGURE 2 - EXTERNAL DRAINAGE AND CSF MONITORING SYSTEM WITH STOPCOCKS in the drain only position.

### 3. Device description

Sophysa external drainage and CSF monitoring system is available as single item or included in complete kits for external CSF drainage. Sophysa external drainage and CSF monitoring system is intended to be used with Sophysa external drainage catheters.

The external drainage and monitoring system is a device which provides aseptic shunting of cerebrospinal fluid (CSF) from the ventricular cavities or lumbar sub-arachnoid spaces into a collecting bag. It also allows intracranial pressure (ICP) to be monitored, samples to be taken and medications to be injected.

The system incorporates (Figure 1):

- a support plate graduated in mmH<sub>2</sub>O and mmHg to adjust intracranial pressure (ICP),
- a 100 ml graduated drip chamber (± 2%), equipped with a vent and hydrophobic filter,
- 3 3-way stopcocks each with an access site and a Luer-Lock locking device,
- 2 access sites for injection/sample taking,
- one single-use drainage bag (650 ml), which is graduated and is equipped with protective cap and hydrophobic filter,
- one 160 cm long drainage tube individualised by a green line incorporated into the thickness of the wall with Luer connections and anti-reflux valve.

#### 4. Sterilization - Decontamination

Products are packed individually in double peel-off, sterile, pyrogen-free packaging, sterilized with ethylene oxide.

##### WARNINGS :

Do not use the products if sterile packaging is open or damaged, or if the expiration date has passed.

This product is intended for SINGLE USE ONLY. It is intended to be used once only for a single patient. Do not re-sterilize or re-use after unpacking.

Resterilization can damage product, potentially leading to patient injury. Re-use of this device may change its mechanical or biological features and may cause device failure, allergic reactions or bacterial infections.

##### NOTE:

Sophysa cannot be held responsible for the performance of any product that has been re-sterilized, nor for any complications which might result from this.

#### 5. Patient education

It is physician's responsibility to educate the patient and /or their representative(s) regarding external drainage and monitoring. This should include a description of associated complications, and explanation of potential alternative products and treatments.

#### 6. Instructions

##### 6.1. EXTERNAL DRAINAGE

###### 6.1.1. Drainage catheter placement

Implanting a system for external drainage is a surgical procedure.

Several techniques can be used to introduce the catheter into the ventricles or the sub-arachnoid spaces. The surgeon should choose the one most suited to the needs of the patient and to his experience.

##### WARNINGS :

Prior drainage catheter placement, read all instructions, precautions and warnings mentioned in drainage catheter instructions for use.

Ensure that the drainage catheter is equipped with a Luer connector.

##### 6.1.2. Priming the system


##### WARNINGS :

To avoid to wet the drip chamber vents during priming, suspend the support plate to the infusion stand using the suspension cord.

The tubing has a dead volume of approximately 15 ml, making it necessary to prime the system prior to use. The procedure should take into account the current guidelines on asepsis used in the establishment.


1. Fill a syringe of suitable volume with a saline solution and purge the air from the syringe.
2. Take off the protective cap from the proximal end of the patient line.
3. Take off the access site stopper of the patient line stopcock and connect this to the syringe.
4. Position the stopcock in the proximal priming position (Figure 3.1). Inject the solution until it reaches the proximal end of the patient line.

FIGURE 3.1 - PATIENT LINE STOPCOCK IN THE PROXIMAL PRIMING, SAMPLING OR INJECTION POSITION.

5. Confirm that the system and drip chamber stopcocks are in the correct drain only position (Figure 2).
6. Place the patient line stopcock in the distal priming position (Figure 3.2).

FIGURE 3.2 - PATIENT LINE STOPCOCK IN THE DISTAL PRIMING POSITION.

7. Inject the saline solution and purge the air from the system. Confirm that there is no leak and that the liquid advances correctly to the drip chamber and drainage bag.
8. Close the clamp on the tubing.
9. Once the priming is completed place the distal stopcock in the drain only position, remove the syringe and replace the stopcock access site stopper.

### 6.1.3. Connecting the drainage system to the catheter

After the catheter has been properly placed and the connection fitting fixed into the distal tip of the catheter :

- Ensure the drainage catheter is clamped at the exit site to allow only a minimal amount of CSF to escape during connection to system,
- Attach the catheter connection fitting to the male Luer connector of the patient line,

In order to ensure the system is correctly set up :

1. Check that there are no bubbles in the tubing,
2. Ensure stopcocks are positioned in the drain only position (Figure 2),
3. Remove the clamp from the tubing, then from the catheter,
4. Ensure flow of CSF from the drip chamber of the shunt system,
5. Ensure there is no leak on the tubing and at the connection between the drainage catheter and the drainage system.

Once the system has been correctly set up, if the patient has to be moved position the 3-way stopcock to off (Figure 3.2) and close catheter and tubing clamps.

When the patient reaches the desired location, the drainage system can be positioned.

### 6.1.4. Positioning the system

#### WARNINGS :

Before adjusting the height of the drainage bag relative to the ventricles, allow the CSF to flow along the length of the tubing so as to remove any bubbles of air.


Intracranial pressure in mmH<sub>2</sub>O is controlled by the height (h on figure 4) of the drip point on the drainage bag relative to the ventricles, allow the CSF to flow along the length of the tubing so as to remove any bubbles of air.

Accordingly, to obtain a desired level of intracranial pressure (h), first adjust the height of the drainage system with the cord locking device: point « 0 » on the plate scale (arrow a on figure 4) must be located at the level of the patient's auditory meatus

Then position the drip chamber arrow by sliding the chamber in the groove of the support plate until the appropriate pressure setting in mmH<sub>2</sub>O or mmHg is reached. Secure the chamber position by screwing the adjustment screw at the back of the support plate.

Confirm that the clamp between the drip chamber and the drainage bag is in the open position.

FIGURE 4 - POSITIONING THE SYSTEM.

**WARNINGS :**

Correct positioning of the drainage system is essential. The intracranial pressure is controlled only by the height of the drip chamber relative to the frontal horns of the lateral ventricles.

**PRECAUTIONS :**

If a shunt equipped with a pressure regulation valve is temporarily externalised, the position of the drip chamber must take account of the valve operating pressure; if the drop counter chamber is raised any level above the system zero, this will result in a counter-pressure which will add to the resistance of the valve (for example, with a valve set to 90 mmH<sub>2</sub>O, in order to maintain and ICP of 110 mmH<sub>2</sub>O, the drip chamber must be positioned against the graduation mark 20 on the mmH<sub>2</sub>O scale).

**6.1.5. Supervision during external drainage**

Once the external drainage system is set up and the drainage catheter connected, following warnings must be taken into account.

**WARNINGS :**

A patient undergoing an external drainage an/or ICP monitoring must be kept under closed supervision with a qualified staff.

Make regular checks to ensure no kinking or compression of the tubing exists and that the clamp and the 3-way stopcock are open. Any obstruction of the tubing may stop drainage and have related pathological consequences.

Ensure there is no leak on the tubing and at the connection between the drainage catheter and the drainage system.

Regularly check that neither the patient, nor the drip chamber are lowered or raised accidentally. Any change in the height of the drip chamber relative to the patient may lead to either overdrainage or under-drainage and their related pathological consequences.

Any change in the height of the drip chamber should be made only on medical instructions, and only by qualified staff.

Do not allow alcohol or other solvents to come into contact with the tubing or collecting bag. Certain components will not tolerate exposure to such solvents. Malfunction or erasing of the markings of collecting bag could result if such solvents contact the outer surface of the device. Use only sterile saline solution to clean the device

In order to ensure against ventricular collapse and the possible consequence of tentorial herniation, always perform a drainage maneuver against a positive pressure head on the order of 20 cmH<sub>2</sub>O or 15 mmHg.

**PRECAUTIONS :**

Do not wet the drip chamber or the drainage bag vents with the drained fluid, in order to avoid influencing the quality of drainage and regulation of ICP.

Ensure that the system is maintained vertical during patient transport in order to avoid wetting the drip chamber and drainage bag.

**NOTES :**

The volume graduations on the drainage bag give approximate values of the CSF drain. To measure the volume precisely, regularly collect the volume accumulated in the drip chamber after a pre-set time.

Flow rate may be calculated by dividing the volume collected in the chamber by the length of the time period chosen.

**6.1.6. To change the bag**

1. Close the drip chamber stopcock and then the clamp located on the tubing between the drip chamber and the drainage bag.
2. Disconnect the Luer-lock, close and remove the full bag.
3. Connect and hook on a new bag immediately.
4. Open the clamp and then the drip chamber stopcock.

**WARNINGS :**

If one of the clamps is closed, no drainage takes place. It is essential to check that both clamps are opened following a change of drainage bag.

**6.2. INJECTIONS AND SAMPLING****WARNINGS :**

Ensure that the products to be injected are compatible with silicone and PVC prior to perform any injection through the port included with the tubing.

The current asepsis procedures of the establishment must be followed for injections and sampling.

Two access sites allow samples to be taken:

- from the patient line stopcock
- from the stopcock at the output of the drip chamber.

1. Place the stopcock in the sampling position (Figure 3.1).
2. Carefully disinfect the silicone puncture site.
3. Puncture and remove the sample.
4. Return the stopcock to the drain only position (Figure 2).

Injections are administered only into the patient line stopcock proximal to the anti-reflux valve.

- Place the stopcock in the injection position (Figure 3.1).
- Carefully disinfect the silicone injection site.
- Perform the injection
- Return the stopcock to the drain only position (Figure 2).

**NOTE :**

Any injection should take into account the dead volume between the injection port and the proximal end of the catheter (approximately 4 ml).


**6.3. ICP MONITORING WITH AN EXTERNAL TRANSDUCER**

ICP may be monitored precisely by placing an ICP measurement device filled with liquid (external sensor):

- either onto the free access site of the patient line stopcock,
- or onto the free access site of the stopcock on the support plate against the level « 0 » (recommended).

The ICP may be measured by turning the stopcock selected in order to make the sensor communicate with the CSF and close drainage (Figure 5).

FIGURE 5 - STOPCOCK IN THE ICP MONITORING POSITION.

**WARNINGS :**

This stopcock position must not be used if drainage is the primary indication of the system. Do not forget to re-establish drainage as soon as the pressure has been measured.


Once monitoring has been completed and the sensor removed, place the stopcock in the drainage position. The protective cap must then be replaced on the free access site of the stopcock.

ICP may be measured whilst allowing drainage at the same time. Turn the stopcock selected in order that the central branch points in the direction of the sensor (Figure 6).

**WARNINGS :**

That device configuration for ICP monitoring can generate artifacts. ICP monitoring is optimal when drainage is interrupted.

FIGURE 6 - STOPCOCK IN SIMULTANEOUS ICP MONITORING AND DRAINAGE POSITION.

It is essential to refer to the manufacturer's instructions for the ICP measurement device in order that this is used and functions correctly.

## 7. Complications

There are three main complications associated with CSF external drainage system.

### 7.1. INFECTIONS

One major complication associated with ICP monitoring with a ventricular or lumbar catheter is the risk of infections, particularly meningitis and ventriculitis. The incidence of these infections can be reduced by care in inserting the catheter and stabilizing it:

- by passing it through a subgaleal (ventricular catheter) or subcutaneous (lumbar catheter) tunnel before it emerges;
- by fixing it to the scalp (ventricular catheter) or cutaneous tissue (lumbar catheter) with the fixation tab.

The incidence of these infections can also be reduced by reducing the number of times a drainage system is opened.

Wound infections may occur but usually subside when the catheter is removed.

Limiting the duration of monitoring from a single site to less than five days will reduce the infection rate. If monitoring must continue past five days, consider inserting a new catheter at a fresh site and changing the entire system.

In case of infection, the drainage system should be removed and specific treatment by general or intrathecal routes should be introduced.

### 7.2. OVERDRAINAGE

A second major complication associated with ventricular or lumbar drainage of CSF is overdrainage, which can lead to intracranial hemorrhage and permanent neurological deficit. Overdrainage can occur due to improper system

testing or setup (resulting in system leakage or inappropriate system pressure head heights).

### 7.3. OBSTRUCTION

A third major complication associated with ventricular or lumbar drainage of CSF is obstruction. Poor recording of ICP will result if the catheter, patient line or other components of the monitoring system become clogged with blood clots, brain tissue fragments, or fibrinous debris.

Frequent punctures of the brain to insert the ventricular catheter can predispose to intracerebral hemorrhage and edema causing a further rise in ICP.

In patients with small ventricles, the ventricle walls may collapse around the tip of the catheter resulting in obstruction and predisposing to tentorial herniation. It is therefore extremely important to avoid excessive release of CSF before the catheter is attached to the patient line.

## 8. Guarantee

Sophysa guarantees that this medical device is free from defects in material or manufacturing. Apart from this guarantee, Sophysa does not grant any other guarantee, express or implicit, including marketing or adaptation for a specific use. Sophysa cannot be held responsible for any incident, complication, damage or prejudice resulting directly or indirectly from the use of this device. Sophysa does not authorize anyone to take responsibility on its behalf for its products.

The performance of external drainage and CSF monitoring system is only guaranteed with the range of devices and accessories designed, tested and manufactured by Sophysa. However, it is possible to use other brands of external drainage catheters provided that they are suitable with Sophysa external drainage and CSF monitoring system.

## 9. Processing of the products after use

### 9.1. DESTRUCTION AFTER USE

An unpacked or used external drainage and CSF monitoring system must be destroyed in accordance with the procedures in force in the medical establishment.

### 9.2. RETURN OF PRODUCTS

If the used external drainage and CSF monitoring system needs to be returned to Sophysa for analysis, indicate if necessary whether cleaning has been performed.

In order to properly assess the returned product, it must be accompanied by an explanatory Return to Manufacturer Authorization form.

## 10. Symbols

Manufacturer

**REF**

Catalog number

**SN**

Serial number

**LOT**

Batch code

Use by

Caution, consult accompanying documents

**STERILE EO**

Sterilized using ethylene oxide

Do not resterilize

Do not reuse

CE conformity marking

## 11. References

Catheters for external drainage of CSF

### DE-106

External ventricular catheter

### DE-105

External ventricular catheter (large lumen)

### DE-108

External lumbar catheter

### DE-109

External lumbar catheter with introducer guide

External drainage systems

### DE-410

External drainage and CSF monitoring system

### DE-420

Replacement drainage bag (x5)

Complete kit for external drainage

### DE-406

Complete kit for external ventricular drainage and CSF monitoring (DE-410+DE-106)

### DE-408

Complete kit for external lumbar drainage and CSF monitoring (DE-410+DE-108)

Product Description Subject To Change

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