Instructions for use: P. 3

PRESSIO® MONITORING KIT, VENTRICULAR TUNNELING

- Model PSO-VT -

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





Instructions for use

1. Indications	5
2. Contra-indications	5
3. Description and Operating Principle for Pressio® ventricular tunneling kits	5
4. Configurations for Pressio® ventricular tunneling kits	6
5. Units of measurement.....	6
6. Sterilization of Pressio® kits.....	6
7. Instructions.....	6
Preparation of a Pressio® catheter BEFORE IMPLANTATION: zeroing the pressure sensor	6
Implantation Technique	7
8. Precautions for monitoring, care and nursing of the patient	9
9. Complications / Side effects	9
Infection.....	9
Cerebral hemorrhage	9
Catheter blocked	10
Overdrainage.....	10
10. Guarantee.....	10
11. Processing of the products after use	10
Destruction after use	10
Return of products	10
12. Symbols	10
Technical specifications	11
References	12

Figure 1: Pressio® ICP monitoring kit, ventricular tunneling (Model PSO-VT).



Figure 2: Pressio® ICP & ICT monitoring kit, ventricular tunneling (Model PSO-VTT)

READ THE INSTRUCTIONS FOR USE CAREFULLY BEFORE IMPLANTING THE DEVICE.

1. Indications

The use of a Pressio® monitoring kit is indicated in patients requiring continuous monitoring of intracranial pressure.

The use of a ventricular kit (Model PSO-VT) makes simultaneous drainage of the cerebro-spinal fluid possible.

In addition, the PSO-VTT model enables the intracranial temperature to be read.

NOTE:

THE TEMPERATURE AND PRESSURE VALUES INDICATED ON PRESSIO® MONITORING SYSTEMS DO NOT PREJUDGE THE HEALTH STATUS OF THE PATIENT.

2. Contra-indications

Contra-indications to the use of a Pressio® kit are the following:

- established or suspected infections in the tissues in direct contact with monitoring system components (meningitis, ventriculitis, septicemia or bacteremia) or any infection present in any part whatsoever of the body
- patients on anticoagulant therapy or presenting with bleeding diathesis.

WARNING:

DO NOT USE THE PRESSIO® MONITORING KIT IF THERE ARE NO TRAINED PERSONNEL AVAILABLE TO PROVIDE CONTINUOUS SURVEILLANCE.

3. Description and Operating Principle for Pressio® ventricular tunneling kits

The Pressio® monitoring catheters are single use implantable devices.

The Pressio® monitoring catheters are of the "BF"  type.

Only use with a Pressio® monitoring system.

PRECAUTION:

ONLY CONNECT PRESSIO® CATHETERS TO CALIBRATED PRESSIO® EQUIPMENT.

Pressio® ICP monitoring kit, ventricular tunneling PSO-VT (Figure 1)

The PSO-VT kit makes continuous monitoring of intracranial pressure and simultaneous drainage of the cerebro-spinal fluid possible.

The 3.1 mm diameter catheter designed to be placed in a cerebral ventricle has three distinct lumina:

- The first contains a 1m long catheter with a diameter of 2.1Fr equipped with a 1.2 mm diameter capsule at its end, enclosing a piezoresistive pressure sensor of the silicon strain gauge type.
- The second contains the pre-inserted introduction stylet.
- The third is for the drainage of the cerebro-spinal fluid.

Pressio® ICP & ICT monitoring kit, ventricular tunneling PSO-VTT (Figure 2)

The PSO-VTT kit makes continuous monitoring of the intracranial pressure and temperature and simultaneous drainage of the cerebro-spinal fluid possible.

The 3.1 mm diameter catheter designed to be placed in a cerebral ventricle has three distinct lumina:

- The first contains a 1m long catheter with a diameter of 2.1Fr equipped with a 1.2 mm diameter capsule at its end, enclosing a piezoresistive pressure sensor of the silicon strain gauge type and a thermistor temperature sensor.
- The second contains the pre-inserted introduction stylet.
- The third is for the drainage of the cerebro-spinal fluid.

The two kits can be distinguished visually by their respective markings on the catheter connector.

In both cases the catheter has depth markers at 5, 10 and 15 cm starting at the proximal end, and centimeter markers between the 5 and 10 cm markers.

The end of the catheter opposite the capsule has a connector which makes it possible to connect the catheter to a Pressio® monitoring system (PSO-3000 or PSO-IN00) via a catheter extension cable (PSO-EC20).

The catheter extension cable transmits the measurement signals in the form of analogue signals from the sensor(s) to the Pressio® monitoring system.

The catheter also exchanges a digital signal containing data, such as sensor calibration, value and date on which the pressure sensor was zeroed, with the Pressio® monitoring system.

This data is memorized in the catheter connector and is thus independent of the Pressio® monitoring system used.

PRECAUTION:

DO NOT ALLOW THE CONNECTOR TO COME INTO CONTACT WITH LIQUID.

The external catheter is made of silicon. At its distal end there is a Luer-Lock type connector for connection to an external drainage system.

The capsule enclosing the sensor(s) is made of titanium and the catheter sheath is made of polyamide.

The entire length of the catheter is radiopaque.


The Pressio® ventricular tunneling monitoring kits are latex-free and phthalate-free products.

The drill bit [1] allows intracranial access.

It has a diameter of 3.5mm.

An adjustable stop [2] makes it possible not to exceed the chosen drilling depth.


The Allen wrench [3] enables the position of this stop on the drill bit to be adjusted.

The tunneling trocar [4] makes it possible to create a passage for the catheter under the scalp. Prior insertion of the tunneling catheter into the trocar sheath makes it possible to pass it under the scalp [5].

The fixation tab [6] enables the catheter to be attached to the scalp. It is made of silicone.

4. Configurations for Pressio® ventricular tunneling kits

The Pressio® ventricular tunneling kits come in two different models: PSO-VT and PSO-VTT (cf. §3: "Description and Operating principle for Pressio® ventricular tunneling kits").

In both cases the specific accessories for the introduction of the catheter into the ventricle and its fixation to the scalp are included in the kit, ready-to-use.

Sophysa also offers a sterile disposable hand drill (PSO-DR), enabling the skull to be perforated in combination with the drill bit supplied in each kit.

Pressio® catheters are used with one of the Pressio® monitoring systems and the appropriate accessories (cf. PSO-3000, PSO-IN00, and PSO-MT00 instructions for use).

Specific case for the PSO-VTT

To make the best use of an intracranial temperature and pressure catheter, a Pressio® monitor, Model PSO-3000 Software Version V2, must be used.

NOTES:

VERSION V2 OF THE SOFTWARE CAN BE IDENTIFIED DURING THE SELF-TEST OF THE MONITOR: ON STARTING THE SOPHYSA LOGO AND SOFTWARE VERSION ARE DISPLAYED.

IF YOU DO NOT HAVE A PRESSIO® SOFTWARE VERSION V2 MONITOR, YOU CAN USE THE CATHETER ON ONE OF THE OTHER PRESSIO® MONITORING SYSTEMS, NAMELY THE PRESSIO® PSO-3000 MONITOR SOFTWARE V1 OR THE PRESSIO® PSO-IN00 INTERFACE.

IN THIS CASE ONLY THE PRESSURE WILL BE DISPLAYED. THE TEMPERATURE WILL NOT BE DISPLAYED BUT THE PERFORMANCE OF THE PRESSURE SENSOR WILL BE IDENTICAL.

5. Units of measurement

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$$D(\text{mm}) = Fr / 3$$

$$Fr = D(\text{mm}) \times 3$$

Intracranial pressure

The pressures mentioned are in mmHg.

1 mmHg corresponds to 13.60 mmH₂O and to 133 Pa.

Intracranial temperature

The intracranial temperature can be displayed in degrees Celsius (°C) or degrees Fahrenheit (°F).

$$^{\circ}\text{F} = ^{\circ}\text{C} \times (9/5) + 32$$

$$^{\circ}\text{C} = (^{\circ}\text{F} - 32) \times (5/9)$$

Equivalence in Kelvin (K):

$$\text{K} = ^{\circ}\text{C} + 273.15$$

$$\text{K} = ^{\circ}\text{F} \times (5/9) + 255.37$$

Consult the Pressio® monitor instructions for use to choose the temperature unit to be displayed.

6. Sterilization of Pressio® kits

The Pressio® monitoring kits are packed individually in double peel-off, sterile, pyrogen-free packaging. They are sterilized with ethylene oxide.

WARNINGS:

DO NOT USE THE PRESSIO® MONITORING KITS IF THE STERILE PACKAGING IS OPEN OR DAMAGED, OR IF THE EXPIRY DATE HAS PASSED.

PRESSIO® MONITORING KITS ARE SINGLE USE DEVICES. DO NOT RE-STERILIZE OR RE-USE AFTER UNPACKING AND/OR EXPLANTATION.

THEIR PERFORMANCE (DRIFT, ELECTRICAL SAFETY, ETC.) COULD BE AFFECTED DURING DECONTAMINATION AND THE RISK OF INFECTION WOULD BE SIGNIFICANT.

NOTE:

SOPHYSA CANNOT BE HELD RESPONSIBLE FOR THE PERFORMANCE OF ANY PRODUCT THAT HAS BEEN RE-STERILIZED AND/OR RE-USED, NOR FOR ANY COMPLICATIONS WHICH MIGHT RESULT FROM THIS.

7. Instructions

Before any implantation, it is recommended that the instructions for use of the Pressio® monitoring system used be consulted for installing and setting it up (Pressio® monitor or interface).

Preparation of a Pressio® catheter BEFORE IMPLANTATION: zeroing the pressure sensor

Before implanting the Pressio®, catheter it is necessary to zero the pressure sensor (calibration in relation to atmospheric pressure).

The intracranial temperature sensor is factory calibrated. Therefore the temperature setting does not need prior zeroing.

Ensure that the asepsis of the Pressio® monitoring kit is maintained throughout this step.

1) Powering up the monitoring system

Connect the catheter extension cable (PSO-EC20) to the Pressio® monitoring system.

Before opening the packaging which guarantees the sterility of the Pressio® monitoring kit, turn on the Pressio® monitoring system (PSO-3000 or PSO-IN00). Refer to the instructions for use of the Pressio® Monitor or Interface.

2) Connecting the catheter to the monitoring system

Make the connection between the catheter and the extension cable (PSO-EC20) respecting the guide pin.

Align the blue arrow on the catheter connector and the blue arrow on the cable connector: the catheter connector will be guided safely into the cable connector.

Push up to the stop to ensure that the connection is complete.

The message "ZERO SENSOR" is then displayed on the Pressio® monitor (PSO-3000). The appearance of error codes, "E001", "E002", "E005" or the message, "CONNECT SENSOR" on the Pressio® monitor after connecting the catheter to the extension cable may be a sign of an incorrect connection.

In this case, make sure that the catheter connector is pushed right up to the stop in the catheter extension cable, and that the latter is correctly connected to the Pressio® monitor.

3) Preparation of the sterile water cup

Prepare a sterile saline solution (or even sterile water) and a cup in a sterile environment. Pour the sterile saline solution (or sterile water) into the cup.

PRECAUTION:

PRESSURE SENSOR ZEROING MUST BE PERFORMED IN SHALLOW LIQUID (~5MM). USING A DEEP RECEPTACLE COULD LEAD TO HYDROSTATIC PRESSURE GREATER THAN ATMOSPHERIC ZERO, WHICH WOULD THEN LEAD TO AN ERRONEOUS REFERENCE ZERO.

4) Zeroing the catheter

Immerse the sensor (metal part) in the liquid completely.

Be careful not to touch the cup with the sensor.

When the sensor is ready, press the zeroing key .

NOTE:

DURING ZEROING KEEP THE SENSOR IMMERSED AND AVOID MOVING IT IN THE CUP.

PRECAUTIONS:

AVOID MAKING ANY CONTACT WITH THE END OF THE CATHETER DURING SENSOR ZEROING. CALIBRATION IN RELATION TO ATMOSPHERIC PRESSURE COULD BE FALSIFIED.

ENSURE THE END OF THE CATHETER IS IMMERSED IN THE STERILE WATER OR A STERILE SALINE SOLUTION DURING ZEROING. SENSOR ZEROING IN CONDITIONS OTHER THAN THOSE RECOMMENDED MAY CAUSE INACCURATE PRESSURE READINGS.

Movement of a task progress bar makes it possible to follow the correct performance of the zeroing.

The procedure takes about 4 seconds.

5) Verification that the zeroing procedure took place correctly

When the progress bar reaches the end, the message "ZERO SENSOR OK – INSERT SENSOR" shows that the procedure has taken place correctly.

Check that the Pressio® monitor screen displays the value "0 mmHg" (or the patient monitor screen if the Pressio® interface is used).

If during zeroing, the sensor suffers parasite pressure variations or is damaged, the screen will display the message "ZERO SENSOR FAILURE – REPLACE SENSOR". In this case, disconnect and then reconnect the catheter, and try the zeroing sequence again.

If this fails again, change the catheter.

Also refer to the instructions for use for the Pressio® monitoring system used (PSO-3000 or PSO-IN00, 5 "Anomaly in operation: symptoms / suggested measures").

NOTE:

SENSOR ZEROING CAN ONLY BE DONE ONCE. IF A NEW ATTEMPT TO ZERO THE SAME SENSOR IS MADE, THE MESSAGE, "SENSOR ALREADY ZEROED" IS DISPLAYED.

Storing data in the catheter connector memory (cf. §3) makes it possible to use Pressio® monitors or interfaces (PSO-3000 models software V2, PSO-3000 software V1 and PSO-IN00) other than the model on which the zeroing was performed.

Specific case for the PSO-VTT

NOTE:

WHEN THE TEMPERATURE READ BY THE SENSOR IS LOWER THAN 20°C OR HIGHER THAN 45°C, THE PRESSIO® MONITOR DISPLAYS: "-- --". THIS MEANS THAT THE MEASUREMENT IS IN THE SENSOR MEASUREMENT RANGE BUT OUTSIDE THE DISPLAY RANGES FOR THE PRESSIO® SYSTEM.

Implantation Technique

Implantation of the Pressio® catheter must take account of current aseptic neurosurgical practices.

Use the insertion accessories supplied by Sophysa in the Pressio® monitoring kit.

The surgeon will choose the technique depending on his own experience and the clinical status of the patient. The final implantation of the device should satisfy the conditions for optimal positioning of the sensor in the ventricle.

PRECAUTION:

DO NOT PERFORM THE IMPLANTATION OF A MONITORING CATHETER WITHOUT HAVING A REPLACEMENT KIT AVAILABLE IN CASE IT IS REQUIRED.

Choice of implantation area

Indication for implantation area: the standard right and left prefrontal areas are the main implantation areas. This region allows the patient to turn his/her head while remaining in the decubitus position without interfering with the intracranial pressure monitoring function. Additionally, in most cases, the incision is made behind the hairline, which is acceptable from an esthetic point of view.

It is recommended that the catheter be tunneled under the scalp to improve its fixation and reduce the risks of infection.

The site where the catheter emerges is generally found 5cm from the burr hole in a posterior position.

Intracranial access

Once the implantation site has been chosen, the area is shaved and prepared aseptically. A local anesthetic is applied in the incision area. This is generally 2 to 3 centimeters in front of the coronal suture on the mid-pupillary line. An incision of about one centimeter is made down to the bone.

Make sure the bony plate is well exposed and perform hemostasis on the wound edges.

The adjustable stop on the drill bit supplied in the kit may be positioned as required by loosening the locking screw with the Allen wrench.

Position the adjustable stop depending upon the drilling depth chosen and retighten the locking screw to maintain this position.

WARNING:

INCOMPLETE TIGHTENING OF THE LOCKING SCREW WILL PREVENT THE ADJUSTABLE STOP FROM PLAYING ITS ROLE, WITH THE RISK OF DRILLING TOO DEEPLY.

The drill bit is then fixed to a helical drill and the perforation is performed through the internal and external skull plates. The surgeon must ensure that any possibility of a parenchymal lesion when he crosses the internal plate is avoided.

After having crossed the internal plate, the drill bit is withdrawn and the hole is irrigated with sterile normal saline.

Incise the dura mater.

Technique recommended for tunneling

- Insert the proximal end of the catheter in the tunneling trocar sheath.
- Make a small incision at the site chosen for the emergence of the catheter (Figure 3.1).
- Starting at the site chosen for the emergence of the catheter, insert the tunneling trocar between the scalp and the skull towards the burr hole (Figure 3.2).
- Pull the catheter out of the tunnel making sure that a length of at least 30 cm comes out of the implantation site, and remove the trocar with the sheath (Figure 3.3).

PRECAUTION:

THE END OF THE TUNNELING TROCAR IS SHARP. INTRODUCE THE CATHETER VERY CAREFULLY.

Catheter implantation

WARNING:

LIMIT THE REPETITION OF INTRACEREBRAL CATHETER IMPLANTATIONS: FREQUENT INCISIONS INTO THE BRAIN TO ENABLE THE INSERTION OF THE CATHETER COULD PREDISPOSE IT TO EDEMA AND INTRACEREBRAL HEMORRHAGE, RESULTING IN A RISE IN INTRACRANIAL PRESSURE.

Ensure that the pressure sensor has been zeroed beforehand (cf. § "Preparation of a Pressio® catheter BEFORE IMPLANTATION: zeroing the pressure sensor").

WARNING:

DO NOT PERFORM THE "SENSOR ZERO" ONCE THE CATHETER IS IMPLANTED ; IT IS ESSENTIAL TO PERFORM THE "SENSOR ZERO" BEFORE IT IS IMPLANTED TO CALIBRATE IT WITH ATMOSPHERIC PRESSURE.

- Holding it by the stylet pre-inserted in the specifically designed lumen, implant the catheter in the direction of the ventricle in accordance with standard techniques.
- Remove the protective cap from the Luer-Lock connector, at the end of the catheter drainage tubing. Check correct positioning in the ventricle by observing the return of the cerebro-spinal fluid at the drainage connector opening.
- If a series of air bubbles and liquid segments appear in the catheter, purge it by letting possible air bubbles move down to the end of the drainage tubing. If the ventricles of the patient are swollen, it could be wise to move the catheter forward by several millimeters beyond the point where the first liquid sample was taken. In this way the end of the catheter will remain in the ventricle during decompression.
- If the ventricular approach fails it is, however, possible to monitor the ICP by leaving the catheter in place in the parenchyma it has traversed. The ICP values measured will be the intra-parenchymal pressure values where the sensor is located. In this case, close the Luer-Lock to limit the risk of infection.
- While holding the catheter in place at the implantation site, remove the stylet by gently bending the catheter between the two arrow-shaped markers (Fig. 3-5).

Connection to an external drainage system

- The catheter can be connected to various external CSF drainage systems. Connect the catheter to the external drainage tubing using a Luer-Lock connector.
- If the catheter is not connected to an external drainage system, close the Luer-Lock to limit the risk of infection.

PRECAUTIONS:

REFER TO THE EXTERNAL DRAINAGE SYSTEM INSTRUCTIONS FOR ITS CONDITIONS FOR USE.

Catheter fixation

- Hold the catheter in place at the implantation site and pull very gently on the end located at the side of the connector until it forms a right angle and rests flat against the skull.
- Close the incision above the burr hole in compliance with standard hospital procedures.
- Use the fixation tab supplied to fix the catheter on the scalp at the site where it emerges (Figure 3.6).
- To keep the catheter in place and reduce the tension, roll the catheter and attach the loop thus formed. Make sure no traction is exercised on the fixation tab during these stages.

WARNING:

MAKE SURE THE LOOP FORMED BY THE CATHETER IS NOT TOO TIGHT (RADIUS OF MORE THAN 15 MM) OR THE DRAINAGE OF THE CSF MAY BE SLOWED.

At this stage the Pressio® monitor displays the mean intracranial pressure in mmHg.

Specific case for the PSO-VTT

At this stage, the Pressio® monitor displays the mean intracranial pressure in mmHg. The temperature is also displayed and will stabilize at its precise value in a maximum of 150 seconds.

NOTE:

A TEMPERATURE MEASUREMENT CAN ONLY BE ACCURATE IF THE SENSOR IS IMPLANTED. THE SENSOR IS NOT SUITABLE FOR MEASURING TEMPERATURE IN THE AIR.

In both cases (PSO-VT and PSO-VTT), instability in the pressure value, the appearance of error codes "E001", "E002", "E005" or the message "CONNECT SENSOR" on the Pressio® monitor after connecting the catheter to the extension cable may be a sign of an incorrect connection.

In this case, make sure that the catheter connector is pushed right up to the stop in the catheter extension cable, and that the latter is correctly connected to the Pressio® monitor.

The appearance of error codes "- -", "999", "E001", "E002", "E005" or the message "CONNECT SENSOR" on the Pressio® monitor after implantation of the catheter, even though the catheter is correctly connected, may be a sign that the sensor on the end of the catheter is damaged.

In this case, try to reposition the catheter, and if it fails again, try with a new catheter.

The maximum recommended duration for the catheter to be implanted is 5 days. After the first 24 hours the Pressio® monitor displays the duration of implantation for the catheter in the form of a message "Implantation: X days". This is then shown throughout the implantation period.

This message starts to flash on the Pressio® monitor from the 6th day of implantation meaning that the catheter has exceeded the recommended implantation duration. If the Pressio® interface is used, the diode showing the implantation duration for the catheter will flash as soon as the catheter exceeds the recommended implantation duration.

Once the catheter is implanted it can be disconnected and then reconnected to the same or another unit from the Pressio® range (PSO-3000 software V2, PSO-3000 software V1, PSO-IN00) without the necessity of recalibration to atmospheric pressure. In fact, the monitoring catheter keeps all the zeroing data in the memory.

Recommended stages for explantation of the catheter

- Stop the Pressio® monitor.
- Disconnect the catheter from the extension cable.
- Proceed with the explantation of the catheter.
- Check the integrity of the explanted catheter.

8. Precautions for monitoring, care and nursing of the patient

WARNING:

DO NOT USE THE PRESSIO® MONITORING KIT IF THERE ARE NO TRAINED PERSONNEL AVAILABLE TO PROVIDE CONTINUOUS SURVEILLANCE.

PRECAUTIONS:

THE PATIENT SHOULD ONLY BE MONITORED BY AN ANESTHETIST, A NEUROSURGEON OR OTHER QUALIFIED PERSON.

IT IS RECOMMENDED THAT THE CLIP BE USED ON THE CATHETER EXTENSION CABLE TO PREVENT ANY DISCONNECTION OF THE CATHETER.

During monitoring, look after the catheter implantation site in compliance with standard hospital procedures.

PRECAUTIONS:

HANDLE THE PATIENT WITH CARE TO PREVENT ANY DISCONNECTION OF THE CABLE OR ANY MOVEMENT OF THE IMPLANTED CATHETER. CHECK THE CONNECTION OF THE CATHETER TO THE EXTENSION CABLE AND TO THE PRESSIO® MONITORING SYSTEM AFTER HANDLING THE PATIENT.

THE PRESSIO® CATHETERS ARE NOT SENSITIVE TO THE EFFECTS THAT COULD OCCUR DURING AN MRI EXAMINATION BUT THEY COULD BE A SOURCE OF ARTIFACTS.

DO NOT USE A PRESSIO® MONITORING SYSTEM AND THE IMPLANTED PRESSIO® CATHETER AT THE SAME TIME AS A HIGH FREQUENCY ELECTRO-SURGICAL INSTRUMENT OR A DEFIBRILLATOR. THE CATHETER AND/OR MONITORING SYSTEM COULD BE DAMAGED OR THEIR OPERATION COULD BE DISTURBED.

9. Complications / Side effects

Complications which may result from the implantation of a Pressio® monitoring system include the inherent risks in any surgical intervention and the insertion of a foreign body.

These complications require the rapid intervention of a doctor.

Infection

The major complication associated with this type of monitoring is infection.

The infection risks can be reduced by adhering to asepsis techniques for the handling and implantation of the catheter and also adhering to the maximum duration for catheter implantation (5 days). If monitoring must be continued after this timescale, it is recommended that a new system is put in place on another site.

Adherence to the tunneling technique recommended in these instructions for use will also make it possible to reduce the risks of infection.

If there is infection, removal of the system is indicated in conjunction with the start of a specific treatment by a general or intrathecal route.

Cerebral hemorrhage

Cerebral hemorrhage may also be observed during monitoring of this type.

The probability of this complication may be reduced by limiting the number of cerebral incisions during the introduction procedure and by ensuring that this procedure is only performed by trained, competent professionals.

Catheter blocked

The implantation of the end of the catheter within reach of the choroid plexus may cause blockages in the inlet openings for the CSF. This phenomenon can also occur when the CSF has a high protein content, contains blood, or if there is intra-ventricular debris, blood clots or tissues.

This obstruction may have an impact on CSF drainage but will not have any effect on the reliability on the measurement of intracranial temperature or pressure.

Overdrainage

Overdrainage is a complication specific to the use of an intraventricular ICP monitoring catheter.

Overdrainage can result in a collapse of the ventricles and the appearance of a subdural hematoma.

10. Guarantee

The performances of the Pressio® monitoring kit are only guaranteed with the range of Pressio® monitoring systems and the accessories designed, tested and manufactured by Sophysa.

Sophysa guarantees that this medical device is free from defects in material or manufacturing. Apart from this guarantee, Sophysa does not grant any other guarantee, express or implicit, including marketing or adaptation for a specific use. Sophysa cannot be held responsible for any incident, complication, damage or prejudice resulting directly or indirectly from the use of this device. Sophysa does not authorize anyone to take responsibility on its behalf for its products.

11. Processing of the products after use

Destruction after use

An unpacked, used or explanted Pressio® monitoring kit must be destroyed in accordance with the procedures in force in the medical establishment.

Return of products

If an explanted Pressio® monitoring kit needs to be returned to Sophysa for analysis, indicating if necessary whether cleaning has been performed.

In order to properly assess the returned product, it must be accompanied by an explanatory Return to Manufacturer Authorization form.

Note:

FOR MORE EFFECTIVE ANALYSIS, IT IS RECOMMENDED THAT THE DEVICE IS NOT CLEANED.

12. Symbols

	Catalog reference
	Manufacturer
	Sterilization Method using Ethylene Oxide
	Do not re-use
	Do not re-sterilize
	Use until
	Batch code
	Serial number
	CE Conformity Marking
	Refer to the instructions for use
	BFTYPE EQUIPMENT: Procuring an appropriate degree of protection against electric shock, having a Type F insulated applied section (floating).
	Temperature conditions for storage and transport
	Store in a dry place
	Fragile, handle with care

Technical specifications

Diameter of the catheter with 3 lumina	9.3 Fr (3.1 mm)
Catheter diameter	2.1 Fr (0.7 mm)
Drainage catheter length	700 mm
Catheter length	1000 mm
Marking of the catheter at 5, 10 and 15 cm on the proximal end, and every centimeter between the 5 and 10 cm markers.	
Pressure	
Type of sensor	Piezoresistive
Reference pressure	Atmospheric pressure
Display range (complete system)	-40 to +100 mmHg
Accuracy of the pressure measurement (Linearity and hysteresis) (complete system)	± 2% reading in the range 0 to +100 mmHg
Bandwidth	> 100 Hz
Functional range of overpressure without damage	-700 mmHg to + 1250 mmHg
Temperature coefficient	0.1 mmHg/°C max.
Input resistance	667 Ω
Output resistance	810 Ω
Excitation voltage	1 to 8 V AC or DC
Drift from zero	Less than 1 mmHg during the first 24 hours at 37°C Less than 2 mmHg in the first week at 37°C
Temperature	
Type of sensor	Thermistor
Display range (complete system)	+20°C to +45°C +68°F to +113°F
Maximum tolerated error (Accuracy) in the reference conditions (complete system)	± 0.2°C max from 25°C to 45°C ± 0.4°C max from 20°C to 25°C ± 0.4°F max from 77°F to 113°F ± 0.7°F max from 68°F to 77°F
Resolution (complete system)	0.1°C / 0.1°F
Leakage current	<10 μA to 120 V _{AC}
Environmental conditions during use (complete system)	
Ambient temperature	10°C to 40°C (50°F to 104°F)
Atmospheric pressure	500 hPa to 1060 hPa
Relative humidity	30% to 75%
Conditions for storage / transport	
Temperature	0°C to 50°C (32°F to 122°F)
Atmospheric pressure	500 hPa to 1060 hPa
Relative humidity	20% to 95%
Weight	16g

References

Name/Description	Reference
PRESSIO® MONITORING KITS	
<i>Only use with a Pressio® monitoring system.</i>	
Pressio® ICP monitoring kit, parenchymal with bolt	PSO-PB
Pressio® ICP & ICT monitoring kit, parenchymal with bolt	PSO-PBT
Pressio® ICP monitoring kit, parenchymal tunneling	PSO-PT
Pressio® ICP & ICT monitoring kit, parenchymal tunneling	PSO-PTT
Pressio® ICP monitoring kit, ventricular tunneling with external CSF drainage function	PSO-VT
Pressio® ICP & ICT monitoring kit, ventricular tunneling with external CSF drainage function	PSO-VTT
PRESSIO® MONITORING SYSTEMS	
ICP monitor Mains power supply cable and Catheter extension cable included	PSO-3000
ICP interface control unit	PSO-IN00
PRESSIO® ACCESSORIES	
Catheter extension cable	PSO-EC20
Patient monitor cable <i>Only use with a Pressio® monitoring system.</i>	PSO-MCxx
<i>PHILIPS (AGILENT) 12 PINS</i>	-MC01
<i>SIEMENS (SIRECUST) 10 PINS</i>	-MC02
<i>SPACELABS & MINDRAY 6 PINS</i>	-MC03
<i>GE DATEX-Ohmeda 10 PINS</i>	-MC04
<i>GE Solar (MARQUETTE) 11 PINS</i>	-MC05
<i>HELLIGE 10 PINS</i>	-MC06
<i>SIEMENS 7 PINS</i>	-MC07
<i>NIHON KOHDEN 5 PINS</i>	-MC08
<i>DATASCOPE 6 PINS</i>	-MC10
Intracranial temperature module	PSO-MT00
Enables the temperature value to be displayed on the patient monitor <i>Use only with a Pressio® PSO-3000 monitoring system.</i> <i>To be used with appropriate patient monitor cables (pressure and temperature).</i>	
Patient monitor cable / temperature <i>Use only with the Pressio® intracranial temperature module.</i>	PSO-MCT-y
<i>PHILIPS (AGILENT) - 2 PINS</i>	-MCT-A
<i>SIEMENS - 7 PINS</i>	-MCT-B
<i>SPACELABS - 10 PINS</i>	-MCT-C
<i>GE Solar (MARQUETTE), GE DATEX-Ohmeda - 11 PINS</i>	-MCT-E
<i>HELLIGE, DATEX-Ohmeda, NIHON KOHDEN, MINDRAY & DATASCOPE - JACK 6.35 mm</i>	-MCT-F
Disposable hand drill	PSO-DR
Serial transmitter	PSO-TX00

Technical specifications and List of product references may be modified without notice.
Availability may vary according to country.

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Sophysa USA Inc.: