
CSF SHUNT ACCESSORIES



 **SOPHYSA**
The way to the future



Sophysa has developed a range of accessories dedicated to CSF shunt implantation.

This line offers simple and adapted solutions to make operative gesture easier.

<p>CO-1010 (x5)</p> <p>Titanium symmetrical connector</p> <p>OD: 1.9 mm ID: 1.2 mm</p>	<p>CO-2010 (x5)</p> <p>Titanium asymmetrical connector</p> <p>OD: 1.9 mm OD: 1.2 mm ID: 1.2 mm ID: 1 mm</p>	<p>CO-3010</p> <p>Polypropylene 3-way connector</p> <p>OD: 2 mm ID: 1 mm</p>
<p>CO-LPS50 (x5)</p> <p>Polypropylene Luer connector</p>	<p>TU-9321</p> <p>Tuohy needle - 14 G</p>	
<p>TUN-45 L=45 cm TUN-65 L=65 cm</p> <p>Stainless steel sterilisable catheter passer</p>		
<p>TUG-45 (x25) L=45 cm TUG-65 (x25) L=65 cm</p> <p>Polyethylene disposable sheath</p>		



Distribution

Head Office :

22, rue Jean Rostand
Parc Club Orsay Université
F- 91893 Orsay Cedex
Tel. : +33 (0) 1 69 35 35 00
Fax : +33 (0) 1 69 35 36 90
www.sophysa.com

**Verbindungsbüro
Deutschland :**

Tel. : +49 (0) 681 99 63 385
Fax : +49 (0) 681 99 63 111