

Patient Identification Card

Patient implanted with a Sophy® Mini Monopressure valve

Keep this card with you at all times

Patient's name:

Emergency contact (Name + Tel):
.....

Hospital: City:

Doctor: Tel.:



Date of implantation:

Type of shunt:

V-P

V-A

L-P

Valve implantation:

Skull

Thoracic region

Lumbar region

Anti-siphon:

Yes

No

Stick the traceability label here
with the valve diagram

If you experience symptoms such as headache, vomiting, visual or hearing disorders, sleepiness, walking difficulties, failing mental function or loss of coordination; prolonged, unexplained fever, redness, edema or a skin anomaly along the course of the shunt, contact your surgeon.

Non-clinical testing has demonstrated that the Sophy® Mini Monopressure valve is MR Conditional with the following conditions:

- Static magnetic field of 1.5 or 3.0 tesla
- Maximum spatial gradient magnetic field of 1900 gauss/cm (19 T/m)
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 4W/kg for 15 minutes (First Level Controlled Operating Mode)

Sophysa SA | www.sophysa.com

5, rue Guy Moquet - 91400 Orsay - France
+33 (0)1 69 35 35 00

