



POSITIONNEMENT DU PATIENT POUR LA PRISE DE RADIOGRAPHIE ET POINTS DE REPÈRES DE LA VALVE




Sophy® est une marque enregistrée par Sophysa

©Sophysa 2020

DFSh-SMB-02(FR) Rev000 - 05/2020



LECTURE DU RÉGLAGE DE LA VALVE



	SM8-140	SM8	SM8-300	SM8-400
P1	10 mm H ₂ O	30 mm H ₂ O	50 mm H ₂ O	80 mm H ₂ O
P2	25 mm H ₂ O	50 mm H ₂ O	75 mm H ₂ O	120 mm H ₂ O
P3	40 mm H ₂ O	70 mm H ₂ O	100 mm H ₂ O	150 mm H ₂ O
P4	60 mm H ₂ O	90 mm H ₂ O	125 mm H ₂ O	190 mm H ₂ O
P5	80 mm H ₂ O	110 mm H ₂ O	150 mm H ₂ O	230 mm H ₂ O
P6	100 mm H ₂ O	140 mm H ₂ O	180 mm H ₂ O	270 mm H ₂ O
P7	120 mm H ₂ O	170 mm H ₂ O	220 mm H ₂ O	330 mm H ₂ O
P8	140 mm H ₂ O	200 mm H ₂ O	300 mm H ₂ O	400 mm H ₂ O