


# PATIENT POSITIONING FOR THE RADIOGRAPHIC EXAM AND LANDMARKS OF THE VALVE



Polaris® is a registered trademark by Sophysa


©Sophysa 2021

DFSh-SPV-02(EN) Rev000 - 01/2021



## READING THE OPERATING PRESSURE

POLARIS®



P1

P2

P3

P4

P5

	SPV-140	SPV	SPV-300	SPV-400
P1	10 mm H <sub>2</sub> O	<b>30</b> mm H <sub>2</sub> O	50 mm H <sub>2</sub> O	80 mm H <sub>2</sub> O
P2	40 mm H <sub>2</sub> O	<b>70</b> mm H <sub>2</sub> O	100 mm H <sub>2</sub> O	150 mm H <sub>2</sub> O
P3	80 mm H <sub>2</sub> O	<b>110</b> mm H <sub>2</sub> O	150 mm H <sub>2</sub> O	230 mm H <sub>2</sub> O
P4	100 mm H <sub>2</sub> O	<b>150</b> mm H <sub>2</sub> O	220 mm H <sub>2</sub> O	330 mm H <sub>2</sub> O
P5	140 mm H <sub>2</sub> O	<b>200</b> mm H <sub>2</sub> O	300 mm H <sub>2</sub> O	400 mm H <sub>2</sub> O